BOARD OF TESTS FOR ALCOHOL AND DRUG INFLUENCE

INSTRUCTIONS FOR LAW ENFORCMENT OFFICER FOR USE OF BLOOD ALCOHOL/DRUG SPECIMEN COLLECTION KIT

- STEP 1 If kit was not sealed when first opened, discard entire kit and obtain another kit. If kit has expired, kit should be discarded and another kit obtained (hospital can use appropriate substitute items, preferably gray stopper vials).
- The officer should first fill out information requested on the **BLOOD TEST OFFICER'S AFFIDAVIT** form. It is provided with this kit and must be submitted with the kit.
- **STEP 3** The officer must witness blood withdrawal.
- Open the styrofoam blood tube holder and allow a qualified blood collector to collect blood specimens from the arrestee.

INSTRUCTIONS FOR QUALIFIED BLOOD COLLECTOR

- STEP 5 Using the items provided in this kit or appropriate substituted items, withdraw blood specimens allowing all three (3) blood tubes to fill. To assure proper blood mixing with anti-coagulant/preservative powder, invert the tubes after collection. After collection, discard used needle, needle holder, and alcohol free prep pads.
- STEP 6 The blood collector must sign the BLOOD TEST OFFICER'S AFFIDAVIT form.
- The officer is to complete the information requested on the three (3) blood specimen tubes seals and one (1) internal kit seal. All specimen seals must be completed.
- Using the three (3) blood specimen seals, wrap one (1) seal lengthwise over the top of each of the three (3) tubes.
- Return the three (3) sealed tubes to the styrofoam holder. Ensure that the tubes are properly positioned to prevent breakage.
- STEP 10 Close the styrofoam holder and seal with the internal kit seal. Overlap the seal on the top of the styrofoam holder so that the information completed on the seal is positioned on top of the styrofoam holder.
- **STEP 11** The officer is to complete the information requested on the label on the top of the styrofoam holder.
- Place the sealed styrofoam holder in the re-sealable plastic bag, and close the end of the bag. <u>Do not remove the liquid absorbing cotton pad from the re-sealable plastic bag.</u>
- **STEP 13** Return sealed styrofoam holder to mailing container. <u>Do not seal the kit mailing box at this time.</u>
- Affix enclosed biohazard label where indicated on top of the kit mailing box, then return the two remaining kit seals to mailer box; do not use these at this time.
- STEP 15 Complete the remainder of the **BLOOD TEST OFFICER'S AFFIDAVIT** form and follow the directions at the bottom for distribution of form. Place the original in the mailing box. Close the mailing box and affix the two remaining seals on the mailing box and mail or deliver to the OSBI laboratory.

				BL	OOD	TES	ST O	FFIC	CER'S	S AFF	ID/	AVIT	•				Section 1
UNDER 21 YES NO CO#			CO#	CITY# CITATION NUMBER				BER	INCIDENT NUMBER (IF APPLICABLE)								
ATE OF INCIDENT/ARREST	TIME OF IN	CIDENT/ARRE	EST (MILITAI	RY)	INCI	DENT LOC	CATION				CIT	Υ			COUN	NTY	
NAME OF SUBJECT (LA	ST, FIRST, MI	DDLE)						DOB (N	1O/DAY/YR)	HE	GHT	WE	IGHT	RACE		SEX
ADDRESS OF SUBJECT								CITY					STA	T.	ZIP CO	205	
ADDRESS OF SUBJECT								CITY					SIA	IIE.	ZIP CC	JDE	
DL NUMBER	EXP. (MO/	YR)	STATE	CLASS	RESTRICTIO	NS EN	NDORSEMEN'		MV/CDL	HAZ MAT PLACARDEI Y	D	EHICLE MA	KE	VEHICLE M	ODEL	TAG #	
For Implied Consent control of a motor ve	hicle upon	a public roa	d, highwa	y, street, t	urnpike, pub	lic parkir	ng lot, or oth	her publi	place or	upon a private	e road, s	street, alle				-	
single or multi-family Describe driving beha																	
Describe arriving being	ivior or circo	amstances (, 11 c, com	,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,													
Describe person's co	ndition (odo	r, actions, S	FST, etc.):														
THE PERSON V	VAS INFO	DRMED	OF THE	IMPLIE	D CONS	ENT TI	EST REQ	UEST?		,	YES		N	10			
OTHER WITNESSES	i :																Section 2
A. Name:				Title:				Addres	s:					Phone:			
B. Name:				Title:				Addres	s:					Phone:			
ight Arm Blood Draveft Arm Blood Drav	aw	Notice:	The arm		ion of the	blood di	raw is not			itate and is p	orovide			est of law er		nent fo	
COLLECTION DATE	TIME (MI	LITARY)	BLO	OD KIT NU	MBER (6 dig	its)		#VIA	LS COLLEC	TED PLAC	CE OF CO	OLLECTION	(HOS	PITAL NAME)		Section 4
SIGNATURE OF PERSO	ON WITHOR	AWING BLO	NOD.					DDIN	T NAME A	ND TITLE							
SIGNATURE OF PERSO	JIN WITHOR	AVVING BLO	ЮВ					PNIN	I NAIVIE A	ND TITLE							
COUNTY OF OFFENSE								DELI	/ERED OR	MAILED BY							
WITNESSED BY OFFIC	ER							PRIN	T NAME A	ND TITLE							BADGE #
NOTE: DO NOT S	SEIZE THE	PERSON	I'S DRIV	ER'S LIC	ENSE BAS	ED UP	ON CONS	ENT TO	THE ST	ATE'S BLO	OD TE	ST.					Section 5
1	iala Tial a	12.0.5.6		42C (() -						- I£ (ما داده			- f :-			Section 6
In accordance w	nin niie	12 0.5. 5	ection 4	426. 15	tate unde	er pe na	aity or pe	erjury	maer tn	ie iaws oi t	JKIan	oma ma	it trie	e roregoin	ig is tr	ue an	a correct.
Date	Plac	e (locatio	n when	signed)	Sig	nature	of Law E	nforce	ment Of	ficer*	Pr	int Nam	e				
		• (3.624,			00										
Agency (list troo	p if appli	cable)			Ad	dress				Z	Zip		_	Agency/T	roop I	Phone	Number
*Upon submiss					-				-	_				his/her si	ignatı	ıre:	
I am aware of a								uce to	custor	ners (USB	i CSD	QIVIA 1	⊥).				
Complete all se					-		-	l distri	bute as	follows:							
		-	-	-	n whose	blood	l was wit	thdrav	n. _Γ								
OF FORM:		the orig															
		one co	-									L/	AB L	JSE ONL	Y	1:	1/2024

BA10K:BTOA1.4 11/24

For I	Implied	Consent	Test	Requests	Only:
-------	---------	---------	------	----------	-------

What Happens Next?

The State of Oklahoma will send you an Order of Revocation and/or Disqualification if action is taken. It will be mailed to the address you have on file with Service Oklahoma. To confirm or change your mailing address on file, contact Service Oklahoma.

The Order of Revocation and/or Disqualification will contain important information about your reinstatement requirements and your right to appeal the revocation and/or disqualification of your driving privileges.

NOTICE OF RIGHT TO REQUEST INDEPENDENT TESTING

State law requires the retention of blood specimens. Your specimen will be retained for **sixty (60)** days from the date of collection. During this period, you may request to have your specimen submitted to an approved laboratory for an independent test at your expense. Contact the OSBI Forensic Science Center at 405-330-6724 to request that your specimen be forwarded to an approved laboratory.