

Dr. Jarrad Wagner, Ph.D. F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

ATTESTATION

I, the undersigned <u>Joshua Smith</u> (Director/Records Custodian), attest under oath pursuant to 12 O.S. §2902 (11)(a) this record was made, at or near the time of the occurrence of the matters set forth by or from information transmitted by a person having knowledge of those matters; was kept in the course of the regularly conducted business activity of the Board; and was made pursuant to the regularly conducted activity. The duplicate record provided is listed below:

1) The attached "Service Overview" is a true and accurate copy of maintenance performed on the Intoxilyzer 8000, serial number 80-003294, in compliance with the agency's Maintenance Policy and Procedures.

Signed

Name of Position: Director/Records Custodian

Date of Attestation: 01/24/2025

Service Overview

3600 N. Martin Luther King Bldg #9 Oklahoma City, OK 73111 405-425-2460

BENCH CHECK DATE: 1/24/2025

SERIAL: 80-003294

DOCUMENTATION WITHIN SERVICE PACKET:

- BOT Certificate of Maintenance, Calibration and Operation
- Bench Check
- BOT Dry Gas Certificate of Analysis
- Ilmo Specialty Gases Certificate of Analysis
- Mock Subject Test
- BOT Certificate of Analysis on Certified Solutions (4)
- Guth Laboratories Certificate of Analysis (4)



Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

CERTIFICATE of MAINTENANCE, CALIBRATION, and OPERATION

This certifies the Intoxilyzer 8000 referenced below, manufactured by CMI, Inc., a subsidiary of MPD, Inc. of Owensboro, Kentucky, is maintained by the Board of Tests for Alcohol and Drug Influence in accordance with 40 O.A.C. 30-1-3. The instrument reported below is hereby deemed to be in compliance with Board rules via linearity accuracy examination and operational/functional procedure evaluation also referred to as a "Bench Check". The Intoxilyzer 8000 appears on the National Highway Traffic Safety Administration's Conforming Products List of Evidential Breath Alcohol Measurement Devices (82 FR 50940 et seq) as an approved device. In accordance with 40 O.A.C. 25-1-2, this instrument is approved and certified for use in the State of Oklahoma as an automated evidential breath-alcohol measurement device by the Board of Tests pursuant to 47 O.S. § 759.

All analyses conducted by this instrument are considered precise, valid, and in accordance with 40 O.A.C. 30-1-3(b) Operating Procedure(s). A printed breath test operator affidavit containing the seal of the Board by this instrument confirms a calibration check successfully verified the instrument's calibration at the time of test. The printed breath test affidavit also serves as proof that safeguards in place to protect the breath samples, including but not limited to Ambient Air Interference, Radio Frequency Interference, Residual Mouth Alcohol, Deficient Samples, etc. were not triggered and the test is a true and accurate scientific measurement of the person's alcohol concentration.

- Maintenance was performed on the Intoxilyzer 8000, Serial No. <u>80-003294</u> on <u>1/24/2025</u> in accordance with 40 O.A.C. 30-1-3 and internal BOT maintenance policy and procedures.
- A dry gas canister appearing on the National Highway Traffic Safety Administration's Conforming Products List of Calibrating Units for Breath Alcohol Testers (77 FR 64588 et seq) and approved in accordance with 40 O.A.C. 25-1-3 and BOT maintenance policy was installed.
 - Manufacturer: ILMO Unit: 208 ppm / 105 L
 - Cylinder Lot No. Canister No.: 31124080A2-031
 - BOT analysis average test result: 0.082
 - Cylinder expiration date: 12/5/2026

In accordance with 40 O.A.C. 25-1-2.1, the approved reference standards and equipment were used:

Simulator Model		Plan 1995 has all daily to be an	O make a principal
Guth 2100	Guth 2100	Guth 2100	Guth 2100
Serial #		44.1.00	
DR3591	DR3694	DR2455	DR3641
Concentration			
.02	.04	.10	.20
LOT#			
24150	23400	24390	24080
Manufactured Date			1
4/16/2024	10/24/2023	10/29/2024	2/13/2024
Expiration Date			
4/16/2026	10/24/2025	10/29/2026	2/13/2026
Solution Commission Date		Washington to	
1/15/2025	1/15/2025	1/22/2025	1/22/2025

Completed by BOT Personnel:

Christopher Pape

Breath-Alcohol Testing Program Administrator

Title

Certified by:

Signature/Name

Joshua Smith
Signature Name

State Director Title

strument type Intoxilyzer		STRUMENT MODEL		INSTRUMENT LOCATION 80-003294				
•		I-8000		80 - 00	ERIAL NUMBER			
			RESULTS					
Diagnostics								
PASS								
Mouth Alcohol T	est							
PASS								
RFI Test								
PASS			Duplicate Bre	ath Sampling R	esults			
Abort Test			Dupiloute Die		100/00/2007			
PASS		Valu	ue Test	1 Test 2	Control Test			
Improper Sample	Test			The Profession	200			
PASS		.0	2 0.01	9 0.019	0.020			
Interferent Detect	Test	.0	4 0.03	8 0.036	0.039			
PASS		.1	0 0.09	7 0.095	0.097			
No Sample Given	Test	.2	0 0.19	5 0.192	0.198			
PASS		. 4	0 0.13		0.150			
Insufficient Sample PASS	e Test			PASS				
Air Blank Contaminat	tion Test							
Range Exceeded PASS	Test							

0000023020

PAPE, CHRISTOPHER



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis - Dry Gas

ALPHA INSTRUMENT
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-002591
01/15/2025 13:08

Test	g/210L	Time
Air Blank	0.000	13:09
Cal Check	0.081	13:09
Air Blank	0.000	13:10
Cal Check	0.082	13:10
Air Blank	0.000	13:10
Cal Check	0.082	13:11
Air Blank	0.000	13:11
Cal Check	0.082	13:12
Air Blank	0.000	13:12
Cal Check	0.082	13:12
Air Blank	0.000	13:13
Cal Check	0.082	13:13
Air Blank	0.000	13:14
Cal Check	0.082	13:14
Air Blank	0.000	13:15
Cal Check	0.081	13:15
Air Blank	0.000	13:15
Cal Check	0.082	13:16
Air Blank	0.000	13:16
Cal Check	0.082	13:17
Air Blank	0.000	13:17
Cal Check Stats		
Average	0.0818	
Std Dev	0.0004	
Rel Std Dev(%)	0.5154	
	mark thanks	

Known Breat	h-Alcohol Equivalent (g/210L) Ratio of Nitrogen and 0.08 g/210L = 208 PPM ÷ 2605*	Ethanol Gas (208 ppm/105 L)
Manufacturer: ILMO	Manufacturer Cert. No.: 16926	Cylinder Size: 105L
Lot No.: 31124080A2	Cylinder No.: 031	Expiration Date: 12/05/2026
prescribed by the State Director of Tests p calibrating unit in accordance with 47 O.S	oursuant to 40 OAC 30-1-3. The above NIST traceable . 759, 40 OAC 25-1-3, and appears on the current confo	000 and was performed in accordance with the procedures reference standard meets the requirements for use as a priming products list of calibrating units for breath alcohol et seq]. * See Dubowski, K.M. and Essary N., Journal of



7 Eastgate Dr. • P.O. Box 790 • Jacksonville, IL 62651-0790 217-245-2183 • Fax: 217-243-7634 • www.ilmoproducts.com

Certificate of Analysis

Certificate ID:

16926

Part #:

BAC105L080T

Cylinder Size:

105L

Lot Number:

31124080A2

Expiration:

12/5/2026

0.080 BAC (For the calibration of instruments used to determine breath alcohol concentration)

Contents:

105 Liters @ 1000 psig 70°F (21°C)

Analytical

Reported

Analytical

Component:

Concentration:

Accuracy

(U, k=2): abs

Method:

Ethanol Nitrogen 208 ppm Balance +/-0.002 BAC(G/210L) NDIR [5.2 ppm]

CMI Inc. 316 East Ninth Street

Distributed by:

Owensboro, KY 42303 Phone 866-835-0690 www.alcoholtest.com

*Reference Standard: NIST Traceable Certified Reference Material - 262.3 µmol/mol Ethanol in Nitrogen - Serial No. ND7017 Lot No. 080722E2

Store in dry area, away from sources of heat, ignition and direct sunlight. Do not allow storage area to exceed 52 °C (125 °F).

Specialty Gas Lab Tech

11-14-2024

Issuance Date

Accreditation #61895

The calibration results within this certificate were obtained at the facility listed above using equipment and standards capable of producing analytical results traceable to NIST, and apply only to the items contained on this certificate. ILMO Products Company makes no warranty or representation as to the suitability of the use of any information provided for any particular purpose. The information use is at the sole discretion and risk of the user. Liability shall be limited to established replacement cost of this material or service. This certificate applies only to the items described and shall not be reproduced other than in full, without written approval from the issuing facility.

OFFICER'S IMPAIRED DRIVING AFFIDAVIT

			AFFI	DAVII						Section
RREST DATE	ARREST			COUNTY #		CITY#	CITATI	ON#		
01/24/2025 RREST LOCATION	104	10		55	CITY			COUN		KA/
BOT	Y 11 11 11 11 11 11 11 11 11 11 11 11 11						les	the second secon	LAHON	
BJECT NAME BREATHTEST,	FRANK		lore			3/15/1	981 STATE	X M	509	175
DORESS 3600 MARTIN RIVER LICENSE #	LUTHER KIN	IG	OKL	AHOMA C	and a Caralter for the	STATE	OK		731:	L1 NDORSEMEN
H083463197				10/31	/202		D	(19.9)		CB-4778943767
HICLE MAKE MOCK		MOCK				MOCK		OK	42 C	N N
the above date, time, and local hway, street, turnpike, public pauence of alcohol and/or other in escribe driving behavior or circuit MOCK SUBJECT	arking lot, or other public place, toxicating substances as prohi mstances): $ ext{TEST}$	or upon a private	I had reasonable gro road, street, alley or I	unds to believe the lane which provides	person was access to o	driving or in actua ne or more single	or multi-family	rol of a motor v dwelling within	ehicle upon a this state whi	public road. le under the
scribe person's condition (odor,		LLED NE	W CYLINDE	R						
			VAS INFORMED OF	A Committee of the second			1007			W 1.16
REATH-ALCOHOL ANALYSIS I	RECORD, REPORT OF THE P		IN SECTION 1 ABOV	E, and LOG OF TE	ST AND MA	AINTENANCE RE	CORD.			Section
I-8000			003294	BOARD	N. T. Carlos T.	ESTS				
TANDARD			TARGET VALU		PERMIT N					
Dry ANUFACTURER			0.080		DATE CYL	∠ U INDER INSTALLI	EO	DATE CYLINI	DER EXPIRES	
ILMO				080A2	01/	24/202	5	12/0	5/202	6
SERVATION START TEST	DATE /24/2025		A 15-minu	ite (or greater						
ESULT TABLE	12412023			vas observed			ing	62	THE ST	A
Test	g/210L	Time	to	the time of te	st admir	istration.	· .	13	**/*	1
Air Blank Subject Test Breath Volume Air Blank Wait Air Blank Subject Test Breath Volume Air Blank	0.00 0.00 1.656 LITERS 0.00 0.00 2.074 LITERS 0.00	11:24 11:25 11:25 11:27 11:27 11:28	complian	ered said I ace with the edure of th	applic	able oper	ating	ONEAT SEA	1007	
Cal Check	0.08	11:29	NAME	D T COO DUE					Test Re	sult
Air Blank	0.00	11:29	BADGE #	RISTOPHER	PERMIT#			1	BrAC g/2	210L
			0000000	0484	000	002302	0		0.0	0
			BOARD O	F TESTS					0.0	-
				appens ne						Section
Upon receipt and review privileges. If you hold a commercial You will continue to have The State of Oklahoma To confirm or change you have order and Notice of the Order and N	I driver license or were d e driving privileges, if oth will send you the Order o our address on file, conta	riving a common nerwise eligible of Revocation a ct Service Okl	vit and any assoc ercial motor vehic , unless or until a and/or Disqualifica ahoma.	ciated test result sle, the State of an Order of Revo ation to the addr	s, if applic Oklahoma ocation an ress you h	may disqualif d/or Disqualific ave on file with	y your privile cation is issu h Service Ol	ege to opera	ite a comm	ercial motor
accordance with Title 1	2 O.S. Section 426 "Le	state under no	nalty of periun	v under the law	s of Okla	homa that the	e foregoing	is true and	correct."	Secti
accordance with ride I			or porjul							
ate:	Place (location when	signed):		Carlos and Carlos and	ignature	of arresting of	flicer:			
AME PAPE, CHRIS	TOPHER			BOARI	OF	TESTS				
ADGE	AGENCY ADDRESS	TON				ZIP CODE			1460	
0000000484	INTER-AGEN	ICY				000	00 40)54252	4460	
THER WITNESSES:								1,172		Secti
NAME			TITLE	ADDRESS			F	HONE		
1			W. W. W.	LDBBBBB				LIONE		

2

8/17

OFFICER'S IMPAIRED DRIVING AFFIDAVIT

					201 20							Section
01/24/2025		ARREST TIM			COUNTY#	eger in	CITY#	CITA	# MOITA			
ARREST LOCATION BOT						CITY				OK	LAHC	MA
BREATHTEST,	FRANK						3/15/		M		509	9 175
ADDRESS 3600 MARTIN	LUTHER	KING		OKLAH	HOMA CI	TY		STA	K		73	
DRIVER LICENSE # H083463197					10/31/		2 OK	D		1	ICTIONS	ENDORSEMENT
MOCK			MOCK				MOCK		70.7	OK.	N CWA	N N

On the above date, time, and location, the above named person was arrested, and I had reasonable grounds to believe the person was driving or in actual physical control of a motor vehicle upon a public road, highway, street, turnpike, public parking lot, or other public place, or upon a private road, street, alley or lane which provides access to one or more single or multi-family dwelling within this state while under the influence of alcohol and/or other intoxicating substances as prohibited by law. (Describe driving behavior or circumstances):

MOCK SUBJECT TEST

Describe person's condition (odor, actions, etc.):

ROUTINE MAINTENANCE, INSTALLED NEW CYLINDER

	the last and the last and the	THE PERSON V	VAS INFORMED OF	THE IMPLIED CO	ONSENT TEST REQUEST	
BREATH-ALCOHOL ANALYSIS	RECORD, REPORT OF THE	PERSON NAMED	IN SECTION 1 ABOV	E, and LOG OF 1	EST AND MAINTENANCE RECORD.	Section 2
I-8000			NUMBER -003294	SPECIALIST BOARD	OF TESTS	
STANDARD Dry			0.080		PERMIT NUMBER 23020	
MANUFÁCTURER ILMO			31124	080A2	01/24/2025	DATE CYLINDER EXPIRES 12/05/2026
0BSERVATION START TEST	DATE L/24/2025				r) deprivation/observation	7
Test	g/210L	Time			prior to and continuing est administration.	OF THE STATE OF
Air Blank Subject Test Breath Volume	0.00 0.00 1.656 LITERS	11:24 11:25	complian	ce with th	Breath-Alcohol Test in e applicable operating	
Air Blank Wait Air Blank	0.00	11:25 11:27 11:27	proce	edure of the		
Subject Test Breath Volume	0.00 0.00 2.074 LITERS	11:28	SIGNATURE OF	OPERATOR		1907
Air Blank Cal Check Air Blank	0.00 0.08 0.00	11:28 11:29 11:29		RISTOPHE		Test Result BrAC g/210L
	-2/2 2.20		BADGE # 0000000 AGENCY	than Donne	0000023020	0.00
			A STATE OF THE STA	F TESTS		0.00

What happens next?

Section 3

- Upon receipt and review of the Officer's Impaired Driving Affidavit and any associated test results, if applicable, the State of Oklahoma may revoke your Oklahoma driving privileges.
- 2. If you hold a commercial driver license or were driving a commercial motor vehicle, the State of Oklahoma may disqualify your privilege to operate a commercial motor vehicle.
- You will continue to have driving privileges, if otherwise eligible, unless or until an Order of Revocation and/or Disqualification is issued by the State of Oklahoma.
 The State of Oklahoma will send you the Order of Revocation and/or Disqualification to the address you have on file with Service Oklahoma.

To confirm or change your address on file, contact Service Oklahoma.

5. The Order and Notice of Revocation will inform you of your requirements to regain driving privileges in Oklahoma.

In accordance with Title 12 O.S. Section 426, "I state und	er penalty of perjury under the laws of Oklahoma that the foregoing is true and correct."	Section 4
Date: Place (location when signed):	Signature of arresting officer:	
PAPE, CHRISTOPHER	BOARD OF TESTS	
BADGE AGENCY ADDRESS TNTER-AGENCY	ZIP CODE PHONE 00000 4054252460	



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis - Reference Solution

ALPHA INSTRUMENT Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-002591 01/15/2025 09:13

Test	g/210L	Time
Air Blank	0.000	09:14
Cal Check	0.021	09:14
Air Blank	0.000	09:15
Cal Check	0.021	09:16
Air Blank	0.000	09:16
Cal Check	0.021	09:17
Air Blank	0.000	09:17
Cal Check	0.021	09:18
Air Blank	0.000	09:19
Cal Check	0.020	09:19
Air Blank	0.000	09:20
Cal Check	0.021	09:21
Air Blank	0.000	09:21
Cal Check	0.021	09:22
Air Blank	0.000	09:22
Cal Check	0.021	09:23
Air Blank	0.000	09:24
Cal Check	0.021	09:24
Air Blank	0.000	09:25
Cal Check	0.020	09:26
Air Blank	0.000	09:26
Cal Check Stats		
Average	0.0208	
Std Dev	0.0004	
Rel Std Dev(%)	2.0270	

Manufacturer: Guth	Manufactured Date: 4/16/2024	Concentration: 0.020
Lot No.: 24150	Expiration Date: 4/16/2026	Equipment Model: Guth 2100
prescribed by the State Director of Tests p to 40 OAC 25-1-2.1. Approved breath-alco	Infrared Spectroscopy utilizing the Intoxilyzer 8000 and vorsuant to 40 OAC 30-1-3. The above NIST traceable reshol measurement equipment and reference solutions. Equits for Breath Alcohol Testers (72 FR 34747), published to Date:	ference solution meets the requirements for use pursuant ulpment used appears on the current or supplemented by the National Highway Traffic Safety Administration.
Name Christopher Pape	Permit No.: 230	20

RCVV 12/3/24 Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 24150 Alcohol Reference Solution for Simulator were analyzed by chromatography on April 17, 2024, using a Perkin Elmer Gas gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.0241% (w/vol) ethyl alcohol. The expiration date for this lot number is April 16, 2026 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.020 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

> Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN03122113 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis - Reference Solution

ALPHA INSTRUMENT
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-002591
05/02/2024 09:26

Test	g/210L	Time
Air Blank	0.000	09:26
Cal Check	0.039	09:27
Air Blank	0.000	09:28
Cal Check	0.039	09:28
Air Blank	0.000	09:29
Cal Check	0.039	09:30
Air Blank	0.000	09:30
Cal Check	0.039	09:31
Air Blank	0.000	09:31
Cal Check	0.039	09:32
Air Blank	0.000	09:33
Cal Check	0.040	09:33
Air Blank	0.000	09:34
Cal Check	0.039	09:35
Air Blank	0.000	09:35
Cal Check	0.039	09:36
Air Blank	0.000	09:37
Cal Check	0.039	09:37
Air Blank	0.000	09:38
Cal Check	0.039	09:38
Air Blank	0.000	09:39
Cal Check Stats	77500	05.35
Average	0.0391	
Std Dev	0.0003	
Rel Std Dev(%)	0.8087	

Manufacturer: Guth	Manufactured Date: 10/24/2023	Concentration: 0.04	
Lot No.: 23400	Expiration Date: 10/24/2025	Equipment Model: Quality 04.00	
to 40 OAC 25-1-2.1. Approved breath-alcol Conforming Products List of Calibrating Uni Signature	Infrared Spectroscopy utilizing the Intoxilyzer 8000 and warsuant to 40 OAC 30-1-3. The above NIST traceable referrol measurement equipment and reference solutions. Equits for Breath Alcohol Testers (72 FR 34747), published by Date:	refrice solution meets the requirements for use pursual ipment used appears on the current or supplemented the National Highway Traffic Safety Administration.	
Name Christopher Pape	Permit No.: 2202	Permit No.: 23020	



RCVD 12/19/23

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 23400 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on October 25, 2023, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.0484% (w/vol) ethyl alcohol. The expiration date for this lot number is October 24, 2025 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.040 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN06032102 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis – Reference Solution

09:30

ALPHA INSTRUMENT Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-002591 01/15/2025

Test · · g/210L Time Air Blank 0.000 09:31 Cal Check 0.100 09:31 Air Blank 0.000 09:32 Cal Check 0.100 09:33 Air Blank 0.000 09:33 Cal Check 0.099 09:34 Air Blank 0.000 09:35 Cal Check 0.099 09:35 Air Blank 0.000 09:36 Cal Check 0.099 09:37 Air Blank 0.000 09:37 Cal Check 0.099 09:38 Air Blank 0.000 09:38 Cal Check 0.100 09:39 Air Blank 0.000 09:40 Cal Check 0.099 09:40 Air Blank 0.000 09:41 Cal Check 0.099 09:42 Air Blank 0.000 09:42 Cal Check 09:43 0.100 Air Blank 0.000 09:43 Cal Check Stats Average 0.0994 Std Dev 0.0005 Rel Std Dev(%) 0.5195

NIST Traceable Breath-A	lcohol Equivalent (g/210L +/- 3%) R	atlo of ethanol and w	ater free of test interfering substances.
Manufacturer: Guth	Manufactured Date: 10	/29/2024	Concentration: 0.100
Lot No.: 24390	Expiration Date: 10/29	9/2026	Equipment Model: Guth 2100
prescribed by the State Director of Tests to 40 OAC 25-1-2.1. Approved breath-ald	pursuant to 40 OAC 30-1-3. The abo cohol measurement equipment and re	ve NIST traceable refe ference solutions. Equ	as performed in accordance with the procedures erence solution meets the requirements for use pursuant ipment used appears on the current or supplemented y the National Highway Traffic Safety Administration.
Name Christopher Pape		Permit No.: 23020	

RCVP 12/3/24

Certified Alcohol Reference Solution for Simulator

cp

Random Samples of Lot Number 24390 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on October 30, 2024, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1205% (w/vol) ethyl alcohol. The expiration date for this lot number is October 29, 2026 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN03072301 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis - Reference Solution

09:45

ALPHA INSTRUMENT Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-002591 01/15/2025

Test	g/210L	Time
Air Blank	0.000	09:45
Cal Check	0.203	09:46
Air Blank	0.000	09:47
Cal Check	0.202	09:47
Air Blank	0.000	09:48
Cal Check	0.202	09:49
Air Blank	0.000	09:49
Cal Check	0.201	09:50
Air Blank	0.000	09:50
Cal Check	0.202	09:51
Air Blank	0.000	09:52
Cal Check	0.201	09:52
Air Blank	0.000	09:53
Cal Check	0.202	09:54
Air Blank	0.000	09:54
Cal Check	0.201	09:55
Air Blank	0.000	09:56
Cal Check	0.202	09:56
Air Blank	0.000	09:57
Cal Check	0.201	09:58
Air Blank	0.000	09:58
Cal Check Stats		
Average	0.2017	
Std Dev	0.0006	
Rel Std Dev(%)	0.3346	

Manufacturer: Guth	Manufactured Date: 2/13/2024	Concentration: 0.200	
Lot No.: 24080	Expiration Date: 2/13/2026	Equipment Model: Guth 2100	
prescribed by the State Director of Test to 40 OAC 25-1-2.1. Approved breath-	I via Infrared Spectroscopy utilizing the Intoxilyzer 8000 and v sts pursuant to 40 OAC 30-1-3. The above NIST traceable re alcohol measurement equipment and reference solutions. Ec g Units for Breath Alcohol Testers (72 FR 34747), published l	ference solution meets the requirements for use pursuant juipment used appears on the current or supplemented	
Signature	Date: 1/15/20	Date: 1/15/2025	



RCVD 12/3/24 CP

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 24080 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on February 14, 2024, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.2419% (w/vol) ethyl alcohol. The expiration date for this lot number is February 13, 2026 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.200 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN03132302 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.

<u>Unit Serial #</u>	Date:
80-003294	1/23/25 ^{17/17}
Old Cylinder #: 33922080A. Unit Old cylinder disposed in inv n unit (After Lab Evaluation) Old c	ventory
Notes from the Field ✓ Routine Swap	ssues in the Field
namber Hold Vacuum?:	
Snap Test ok?	
	Is regulator close to Tank Sensor Value?
Repairs Made:	THE DESCRIPTION OF THE PROPERTY OF THE PROPERT
Ready for Bench Check? Inician Report Completed Completed To complete Completed CMI Workord Tation and Operation	
	Old Cylinder #: Old Cylinder disposed in into unit (After Lab Evaluation) Old color in unit (After Lab Evaluation) Old