

Dr. Jarrad Wagner, Ph.D. F-ABFT Chairman Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

ATTESTATION

I, the undersigned <u>Joshua Smith</u> (Director/Records Custodian), attest under oath pursuant to 12 O.S. §2902 this record was made, at or near the time of the occurrence of the matters set forth by or from information transmitted by a person having knowledge of those matters; was kept in the course of the regularly conducted business activity of the Board; and was made pursuant to the regularly conducted activity. The duplicate record provided is listed below:

 The attached "Service Overview" is a true and accurate copy of maintenance performed on the Intoxilyzer 8000, serial number 80-003358, in compliance with the agency's Maintenance Policy and Procedures.

Signed

Joshua Smith

Name of Position: Director/Records Custodian

JORDAN FAITH THOMPSON NOTARY PUBLIC STATE OF OKLAHOMA COMMISSION NO. 25010889 EXPIRES 09-16-2029 State of Oklahoma

County of Oklahoma

Signed or attested before me on \2/11/25 by

Joshua Smith.

(Signature of notarial officer)

Service Overview



3600 N. Martin Luther King Bldg #1 Oklahoma City, OK 73111 405-425-2460

BENCH CHECK DATE: 12/4/2025

SERIAL: 80-003358

DOCUMENTATION WITHIN SERVICE PACKET:

- BOT Certificate of Maintenance, Calibration and Operation
- Bench Check
- BOT Dry Gas Certificate of Analysis
- Ilmo Specialty Gases Certificate of Analysis
- Mock Subject Test
- BOT Certificate of Analysis on Certified Solutions (4)
- Guth Laboratories Certificate of Analysis (4)



Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

CERTIFICATE of MAINTENANCE, CALIBRATION, and OPERATION

This certifies the Intoxilyzer 8000 referenced below, manufactured by CMI, Inc., a subsidiary of MPD, Inc. of Owensboro, Kentucky, is maintained by the Board of Tests for Alcohol and Drug Influence in accordance with 40 O.A.C. 30-1-3. The instrument reported below is hereby deemed to be in compliance with Board rules via linearity accuracy examination and operational/functional procedure evaluation also referred to as a "Bench Check". The Intoxilyzer 8000 appears on the National Highway Traffic Safety Administration's Conforming Products List of Evidential Breath Alcohol Measurement Devices (82 FR 50940 et seq) as an approved device. In accordance with 40 O.A.C. 25-1-2, this instrument is approved and certified for use in the State of Oklahoma as an automated evidential breath-alcohol measurement device by the Board of Tests pursuant to 47 O.S. § 759.

All analyses conducted by this instrument are considered precise, valid, and in accordance with 40 O.A.C. 30-1-3(b) Operating Procedure(s). A printed breath test operator affidavit containing the seal of the Board by this instrument confirms a calibration check successfully verified the instrument's calibration at the time of test. The printed breath test affidavit also serves as proof that safeguards in place to protect the breath samples, including but not limited to Ambient Air Interference, Radio Frequency Interference, Residual Mouth Alcohol, Deficient Samples, etc. were not triggered and the test is a true and accurate scientific measurement of the person's alcohol concentration.

- 1. Maintenance was performed on the Intoxilyzer 8000, Serial No. <u>80-003358</u> on <u>12/4/2025</u> in accordance with 40 O.A.C. 30-1-3 and internal BOT maintenance policy and procedures.
- 2. A dry gas canister appearing on the National Highway Traffic Safety Administration's Conforming Products List of Calibrating Units for Breath Alcohol Testers (77 FR 64588 *et seq*) and approved in accordance with 40 O.A.C. 25-1-3 and BOT maintenance policy was installed.
 - Manufacturer: ILMO Unit: 208 ppm / 105 L
 - Cylinder Lot No. Canister No.: 25925080A1-060
 - BOT analysis average test result: 0.079
 - Cylinder expiration date: 11/5/2027

3. In accordance with 40 O.A.C. 25-1-2.1, the approved reference standards and equipment were used:

Simulator Model			
Guth 12V500	Guth 12V500	Guth 12V500	Guth 12V500
Serial #			7th
MP9201	MP9202	MP9203	MP9204
Concentration			
.02	.04	.10	.20
LOT#		10/11/	NO NO NO
25010	25090	25100	25020
Manufactured Date			M
1/7/2025	3/11/2025	3/18/2025	1/14/2025
Expiration Date			
1/7/2027	3/11/2027	3/18/2027	1/14/2027
Solution Commission Date			
12/4/2025	12/4/2025	12/4/2025	12/4/2025

Completed by BOT Personnel:

Christopher Pape

Breath-Alcohol Testing Program Administrator

Title

Signature/Na

Signature

Joshua Smith Name State Director Title

12/04/2025	12:24:1	INSTRUMENT LOCATION 80-003358
NSTRUMENT TYPE Intoxilyzer	INSTRUMENT MODEL I - 8000	INSTRUMENT SERIAL NUMBER 80-003358
	TES	ST RESULTS

Diagnostics

PASS

Mouth Alcohol Test

PASS

RFI Test

PASS

Abort Test

PASS

Improper Sample Test

PASS

Interferent Detect Test

PASS

No Sample Given Test

PASS

Insufficient Sample Test

PASS

Air Blank Contamination Test

PASS

Range Exceeded Test

PASS

Duplicate Breath Sampling Results

Value	Test 1	Test 2	Control Test
.02	0.019	0.019	0.020
.04	0.039	0.038	0.039
.10	0.097	0.096	0.099
.20	0.193	0.192	0.197

PASS

CERTIFYING TECHNICIAN

I swear under penalty of perjury that in regards to the above listed instrument, I have complied with the Intoxliyzer Bench Checking Procedure approved by the State of Oklahoma.

DATE

| 12/4/2025
| PAPE, CHRISTOPHER | 0000023020



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis - Dry Gas

ALPHA INSTRUMENT 2
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-004923
11/10/2025 15:55

Test	g/210L	Time
Air Blank	0.000	15:55
Cal Check	0.079	15:56
Air Blank	0.000	15:56
Cal Check	0.079	15:57
Air Blank	0.000	15:57
Cal Check	0.079	15:57
Air Blank	0.000	15:58
Cal Check	0.079	15:58
Air Blank	0.000	15:59
Cal Check	0.079	15:59
Air Blank	0.000	16:00
Cal Check	0.079	16:00
Air Blank	0.000	16:01
Cal Check	0.079	16:01
Air Blank	0.000	16:01
Cal Check	0.079	16:02
Air Blank	0.000	16:02
Cal Check	0.079	16:02
Air Blank	0.000	16:03
Cal Check	0.079	16:03
Air Blank	0.000	16:04
Cal Check Stats		
Average	0.0790	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

	th-Alcohol Equivalent (g/210L) Ratio of Nitrogen and 0.08 g/210L = 208 PPM ÷ 2605*	Edition and (200 blum 100 E)
Manufacturer: ILMO	Manufacturer Cert. No.: 17552	Cylinder Size: 105L
Lot No.: 25925080A1	Cylinder No.: 060	Expiration Date: 11/05/2027
prescribed by the State Director of Tests p	nursuant to 40 OAC 30-1-3. The above NIST traceable	reference standard meets the requirements for use as a
calibrating upit in accordance with 47 O.S. esters published by the U.S. Department Anal. Tox col. (1996) 20(6): pp. 484-491.	. 759, 40 OAC 25-1-3, and appears on the current confort Transportation in the Federal Register [77 FR 64588	prming products list of calibrating units for breath alcohol et seq]. * See Dubowski, K.M. and Essary N., <i>Journal of</i>
calibrating upit in accordance with 47 O.S.	. 759, 40 OAC 25-1-3, and appears on the current confort Transportation in the Federal Register [77 FR 64588	orming products list of calibrating units for breath alcohol et seq]. * See Dubowski, K.M. and Essary N., <i>Journal of</i>



7 Eastgate Dr. • P.O. Box 790 • Jacksonville, IL 62651-0790 217-245-2183 • Fax: 217-243-7634 • www.ilmoproducts.com

Certificate of Analysis

Certificate ID:

17552

Part #:

BAC105L080T

Cylinder Size:

105L

Lot Number:

25925080A1

Expiration:

11/5/2027

0.080 BAC (For the calibration of instruments used to determine breath alcohol concentration)

Contents:

105 Liters @ 1000 psig 70°F (21°C)

Reported

Analytical

Analytical

Component:

Accuracy

Concentration:

(U, k=2): abs

Method:

Ethano1 Nitrogen 208 ppm Balance +/-0.002 BAC(G/210L) NDIR

[5.2 ppm]

Distributed by:

CMI Inc.

316 East Ninth Street Owensboro, KY 42303 Phone 866-835-0690 www.alcoholtest.com

*Reference Standard:

NIST Traceable Certified Reference Material - 260.6 µmol/mol Ethanol in Nitrogen - Serial No. GN0015020 Lot No. 00424E1 Store in dry area, away from sources of heat, ignition and direct sunlight. Do not allow storage area to exceed 52 °C (125 °F).

Spegialty Gas Lab Tech



The calibration results within this certificate were obtained at the facility listed above using equipment and standards capable of producing analytical results traceable to NIST, and apply only to the items contained on this certificate. ILMO Products Company makes no warranty or representation as to the suitability of the use of any information provided for any particular purpose. The information use is at the sole discretion and risk of the user. Liability shall be limited to established replacement cost of this material or service. This certificate applies only to the items described and shall not be reproduced other than in full, without written approval from the issuing facility.

OFFICER'S IMPAIRED DRIVING

			AFFII	DAVII						Section
RREST DATE 12/04/2025	1	ARREST TIME		COUNTY #		CITY#	CITATIC	N #		
REST LOCATION					CITY			COUN	TY LAHC	MA
BJECT NAME	ON LEVER	NE				OF BIRTH 8/09/19	61 N	K	HEIGHT 510	WEIGHT
DRESS 600 MARTIN			OKLA	HOMA C			STATE		731	
IVER LICENSE #	Dolling			EXPIRATION D	ATE	STATE OK	CLASS A	REST	RICTIONS	ENDORSEMEN
HICLE MAKE MOCK		MOCK		100/00		TAG# MOCK		STATE	N CMV	PLACARD N
the above date, time, and loca iway, street, turnpike, public prence of alcohol and/or other in scribe driving behavior or circu IOCK SUBJECT cribe person's condition (odor, POST REPAIR Markets)	urking lot, or other publicoxicating substances mstances): PEST actions, etc.):	ic place, or upon a private as prohibited by law.	e road, street, alley or lai	ne which provides	person was c	Ifving or in actual p e or more single or	hysical contro multi-family d	ol of a motor v	vehicle upon	i a public road, while under the
EATH-ALCOHOL ANALYSIS I	RECORD, REPORT O	330000	WAS INFORMED OF THE				ORD.			Section
UIPMENT - 8000		SERIAL		SPECIALIST BOARD	OF T		54(45)4			
andard Ory		100	TARGET VALUE		PERMIT NU 2302	MBER				
NUFACTURER			LOT#		DATE CYLI	NDER INSTALLED	Ti C	DATE CYLINI		
LMO SERVATION START TEST	DATE		259250		-	14/2025 tion/observat		11/0	5/20:	27
Air Blank Subject Test Breath Volume Air Blank Wait Air Blank Subject Test Breath Volume Air Blank Cal Check Air Blank	0.00 0.00 2.140 LITER 0.00 0.00 2.179 LITER 0.00 0.08 0.00	14:16 14:18 14:18 14:18		OPERATOR ISTOPHER 484	e applica e Board	able operat	200	THE THE STATE OF T	Test R BrAC g	esult /210L
Upon receipt and review privileges. If you hold a commercial You will continue to have The State of Oklahoma To confirm or change you The Order and Notice of	driver license or versions or versions or versions or driving privileges will send you the Court address on file,	vere driving a comm i, if otherwise eligible order of Revocation i contact Service Okl	ivit and any associa ercial motor vehicle s, unless or until an and/or Disqualificati ahoma.	o, the State of o Order of Revo ion to the addr	s, if applica Oklahoma i ocation and ess you ha	may disqualify y /or Disqualificat ve on file with S	our privileç ion is issue Service Okl	je to opera	ite a comr	mercial motor
accordance with Title 1	2 O.S. Section 42	26, "I state under po	enalty of perjury	under the law	s of Oklah	oma that the fo	oregoing is	true and	correct."	Section
ate:	Place (location	when signed):		S	ignature of	arresting offic	er:			
PAPE, CHRIS	TOPHER AGENCY ADDRESS			BOARI	OF '	TESTS ZIP CODE	PHONE			
0000000484	INTER-A	GENCY				0000	0 40	54252	2460	****
HER WITNESSES:			TITLE	ADDRESS			lou	ONE		Section
Lavaterie.			111166	MUNICOS			len	No. of April		

TITLE

ADDRESS

PHONE

1

2

NAME

OFFICER'S IMPAIRED DRIVING AFFIDAVIT

										Section
12/04/2025	ARREST TIME 1315		COUNTY#	C	ITY#	CIT	ATION#	·		
RREST LOCATION BOT			CI	ITY				OKLA	НОМ	A
NORTH, GAYLON LEVE	RNE				OF BIRTH / 09/1	1961	SEX M	HEI 5	знт 10	WEIGHT 165
3600 MARTIN LUTHER	KING A	OKLAHO	MA CIT	Ϋ́	olt e	2000	OK.	2010.0	313	3
T083454871			6/30/2	019	OK	CLASS		RESTRICTIO		IDORSEMENT M
PHICLE MAKE MOCK	MOCK	-	- N	2.0	ag# MOCK		(75.00)	OK N	OV. 19.32	acard N

On the above date, time, and location, the above named person was arrested, and I had reasonable grounds to believe the person was driving or in actual physical control of a motor vehicle upon a public road, highway, street, turnpike, public parking lot, or other public place, or upon a private road, street, alley or lane which provides access to one or more single or multi-family dwelling within this state while under the influence of alcohol and/or other intoxicating substances as prohibited by law.

(Describe driving behavior or circumstances): MOCK SUBJECT TEST

Describe person's condition (odor, actions, etc.):

POST REPAIR MAINTENANCE, INSTALLED NEW CYLINDER

THE PERSON WAS INFORMED OF THE IMPLIED CONSENT TEST REQUEST

BREATH-ALCOHOL ANALYSIS	RECORD, REPORT OF THE P	ERSON NAMED	IN SECTION 1 ABOV	/E, and LOG OF T	EST AN	MAINTENANCE RECORD.	Section
EQUIPMENT			NUMBER	SPECIALIST	(00.29935) (5)	CONTRACTOR MANAGEMENT OF THE PARTY OF THE PA	
I-8000		80-	-003358	BOARD	OF	TESTS	
STANDARD			TARGET VALU	JÉ .	PERM	IT NUMBER	
Dry			0.080		23	020	
MANUFACTURER			LOT#	NOSS WILLIAMS WAY IN	DATE	CYLINDER INSTALLED	DATE CYLINDER EXPIRES
ILMO			25925	080A1	12	/04/2025	11/05/2027
OBSERVATION START TEST	DATE		A 15-min	ite (or greate	r) den	rivation/observation	T
1345 12	2/04/2025						
RESULT TABLE	***					to and continuing	THE STAR
Test	g/210L	Time	to	the time of to	est adı	ninistration.	0° * * ^ * * 0
			=				~ * * * * * * * * * * * * * * * * * * *
Air Blank	0.00	14:15	I administ	tered said I	3reat	h-Alcohol Test in	***
Subject Test	0.00	14:16	2.5			licable operating	10 × 00 × 10 × 10
Breath Volume	2.140 LITERS	(and the state of t	日本本文化 直動社 4代本本 15
Air Blank	0.00	14:16	proc	edure of the	e Bo	ard of Tests.	
Wait	0.2020.0003	14:18		()		/ >	
Air Blank	0.00	14:18		16	-	Re-	
Subject Test	0.00	14:18	4				1907
Breath Volume	2.179 LITERS		SIGNATURE C	F OPERATOR			1907
Air Blank	0.00	14:19		X-20			
Cal Check	0.08	14:19	NAME				Test Result
Air Blank	0.00	14:20	PAPE, CH	RISTOPHER			BrAC g/210L
			BADGE#		PERM		DIAC 9/210L
			0000000	0484	00	000023020	0 00
			AGENCY	NESS RECEIVANT OF THE PROPERTY OF THE			0.00
			BOARD O	F TESTS			

What happens next?

Section :

- 1. Upon receipt and review of the Officer's Impaired Driving Affidavit and any associated test results, if applicable, the State of Oklahoma may revoke your Oklahoma driving privileges.
- 2. If you hold a commercial driver license or were driving a commercial motor vehicle, the State of Oklahoma may disqualify your privilege to operate a commercial motor vehicle.
- You will continue to have driving privileges, if otherwise eligible, unless or until an Order of Revocation and/or Disqualification is issued by the State of Oklahoma.
 The State of Oklahoma will send you the Order of Revocation and/or Disqualification to the address you have on file with Service Oklahoma.

To confirm or change your address on file, contact Service Oklahoma.

5. The Order and Notice of Revocation will inform you of your requirements to regain driving privileges in Oklahoma.

In accordance with Title 12 O.S. Section 426, "I state under p	enalty of perjury under the laws of Oklahoma that the foregoing is true and correct."	Section 4
Date: Place (location when signed):	Signature of arresting officer:	
PAPE, CHRISTOPHER	BOARD OF TESTS	
BADGE AGENCY ADDRESS INTER-AGENCY	ZIP CODE PHONE 00000 4054252460	



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis - Reference Solution

ALPHA INSTRUMENT
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-002591
07/09/2025 12:20

Test	g/210L	Time
Air Blank	0.000	12:20
Cal Check	0.019	12:21
Air Blank	0.000	12:21
Cal Check	0.019	12:22
Air Blank	0.000	12:23
Cal Check	0.019	12:23
Air Blank	0.000	12:24
Cal: Check	0.019	12:25
Air Blank	0.000	12:25
Cal Check	0.020	12:26
Air Blank	0.000	12:26
Cal Check	0,019	12:27
Air Blank	0.000	12:28
Cal Check	0.019	12:28
Air Blank	0.000	12:29
Cal Check	0.020	12:30
Air Blank	0.000	12:30
Cal Check	0.019	12:31
Air Blank	0.000	12:31
Cal Check	0.020	12:32
Air Blank	0.000	12:33
Cal Check Stats		
Average	0.0193	
Std Dev	0.0004	
Rel Std Dev(%)	2.5026	

Manufacturer: Guth	Manufactured Date: 1/7/2025	Concentration: 0.020
Lot No.: 25010	Expiration Date: 1/7/2027	Equipment Model: Guth 2100
prescribed by the State Director of Tests pursu to 40 OAC 25-1-2.1, Approved breath-alcohol	measurement equipment and reference solutions. Ed	ference solution meets the requirements for use pursuant uipment used appears on the current or supplemented by the National Highway Traffic Safety Administration.
Name Christopher Pape	Permit No.: 230	



RCVD 5/30/25 CP

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 25010 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on January 8, 2025, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.0235% (w/vol) ethyl alcohol. The expiration date for this lot number is January 7, 2027 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.020 g/210L+/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN03122113 whose values are traceable to NIST.



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis - Reference Solution

ALPHA INSTRUMENT 2
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-004923
07/01/2025 11:06

Test	g/210L	Time
Air Blank	0.000	11:07
Cal Check	0.037	11:07
Air Blank	0.000	11:08
Cal Check	0 037	11:09
Air Blank	0.000	11:09
Cal Check	0.038	11:10
Air Blank	0.000	11:10
Cal Check	0.038	11:11
Air Blank	0.000	11:12
Cal Check	0.037	11:12
Air Blank	0.000	11:13
Cal Check	0.038	11:14
Air Blank	0.000	11:14
Cal Check	0.038	11:15
Air Blank	0.000	11:15
Cal Check	0.038	11:16
Air Blank	0.000	11:17
Cal Check	0.038	11:17
Air Blank	0.000	11:18
Cal Check	0.038	11:19
Air Blank	0.000	11:19
Cal Check Stats		
Average	0.0377	
Std Dev	0.0004	
Rel Std Dev(%)	1.2812	

Manufacturer: Guth	Manufactured Date: 3/11/2025	Concentration: 0.040
Lot No.: 25090	Expiration Date: 3/11/2027	Equipment Model: Guth 2100
prescribed by the State Director of Tests pursu to 40 OAC 25-1-2.1. Approved breath-alcohol	measurement equipment and reference solutions. Equ for Breath Alcohol Testers (72 FR 34747), published by	rence solution meets the requirements for use pursuant pment used appears on the current or supplemented the National Highway Traffic Safety Administration.
holden 14	7/1/2025	5
Name Christopher Pape	Permit No.: 2302	0



RCV1) 5/30/25

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 25090 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on March 12, 2025, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.0480% (w/vol) ethyl alcohol. The expiration date for this lot number is March 11, 2027 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.040 g/210L+/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Paulcy, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN06032102 whose values are traceable to NIST.



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis - Reference Solution

ALPHA INSTRUMENT
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-002591
09/30/2025 11:20

Test	g/210L	Time
Air Blank	0.000	11:21
Cal Check	0.099	11:21
Air Blank ·	0.000	11:22
Cal Check	0.099	11:23
Air Blank	0.000	11:23
Cal Check	0.099	11:24
Air Blank	0.000	11:25
Cal Check	0.099	11:25
Air Blank	0.000	11:26
Cal Check	0.099	11:27
Air Blank	0.000	11:27
Cal Check	0.099	11:28
Air Blank	0.000	11:28
Cal Check	0.099	11:29
Air Blank	0.000	11:30
Cal Check	0.099	11:30
Air Blank	0.000	11:31
Cal Check	0.099	11:32
Air Blank	0.000	11:32
·Cal Check	0.099	11:33
Air Blank	0.000	11:33
Cal Check Stats		
Average	0.0990	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	
	1070/07/2010/00/00/00/00/00	

Manufacturer: Guth	Manufactured Date: 3/18/2025	Concentration: 0.100
-ot No.: 25100	Expiration Date: 3/18/2027	Equipment Model: Guth 12V500
prescribed by the State Director of Tests o 40 OAC 25-1-2.1. Approved breath-alc	a Infrared Spectroscopy utilizing the Intoxilyzer 8000 and pursuant to 40 OAC 30-1-3. The above NIST traceable reports of the pursuant to 40 OAC 30-1-3. The above NIST traceable reports of the pursuant to 40 OAC 30-1-3. The above NIST traceable reports of the pursuant traceable	eference solution meets the requirements for use pursuant quipment used appears on the current or supplemented by the National Highway Traffic Safety Administration.

RCVI) 25

Certified Alcohol Reference Solution for Simulator

d by

Random Samples of Lot Number 25100 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on March 19, 2025, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1211% (w/vol) ethyl alcohol. The expiration date for this lot number is March 18, 2027 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN04042408 whose values are traceable to NIST.



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis - Reference Solution

ALPHA INSTRUMENT 2 Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-004923 06/06/2025 11:23

Test	g/210L	Time
Air Blank	0.000	11:24
Cal Check	0.200	11:24
Air Blank	0.000	11:25
Cal Check	0.199	11:26
Air Blank	0.000	11:26
Cal Check	0.199	11:27
Air Blank	0.000	11:27
Cal Check	0.198	11:28
Air Blank	0.000	11:29
Cal Check	0.199	11:29
Air Blank	0.000	11:30
Cal Check	0.199	11:31
Air Blank	0.000	11:31
Cal. Check	0.199	11:32
Air Blank	0.000	11:33
Cal Check	0.197	11:33
Air Blank	0.000	11:34
Cal Check	0.199	11:35
Air Blank	0.000	11:35
Cal Check	0.198	11:36
Air Blank	0.000	11:36
Cal Check Stats		
Average	0.1987	
Std Dev	0.0008	
Rel Std Dev(%)	0.4143	

Manufacturer: Guth	Manufactured Date: 1/14/2025	Concentration: 0.200	
Lot No.: 25020	Expiration Date: 1/14/2027	Equipment Model: Guth 2100	
prescribed by the State Director of Tests to 40 OAC 25-1-2.1, Approved breath-al	via Infrared Spectroscopy utilizing the Intoxilyzer 8000 and a pursuant to 40 OAC 30-1-3. The above NIST traceable cohol measurement equipment and reference solutions. In Units for Breath Alcohol Testers (72 FR 34747), published Date:	reference solution meets the requirements for use pursuant Equipment used appears on the current or supplemented I by the National Highway Traffic Safety Administration.	
Name Christopher Pape	Permit No.: 23	Permit No.: 23020	



RCVI) 5/30/25

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 25020 of
Alcohol Reference Solution for Simulator were analyzed by
gas chromatography on January 15, 2025, using a Perkin Elmer Gas
Chromatograph Autosystem XL S/N: 610N9030209, and found to contain
0.2409% (w/vol) ethyl alcohol. The expiration date for this lot
number is January 14, 2027 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.200 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN03132302 whose values are traceable to NIST.

Location Returned From Coal Co	<u>Unit Serial #</u> 80-003358	Date: 7/17 11/17/25	
Inventory Process Unit checked in to BOT New Unit Checked out Very Cold Cylinder #: Old Cylinder #: Old Cylinder #: Old Cylinder disposed in inventory Remove old cylinder from unit (After Lab Evaluation) Old cylinder archived			
Check in Procedure Connect unit to ethernet connection 172.019.004. 130 Download any stored records 3 Change unit location to Serial number Update Unit Date and time Previous time: 12/2/25 1524 New Time: 12/2/25 1524 Execute a new Forms load on the unit	Notes from the Field	W. C.	
Lab Evaluation Breath Hose Hold Vacuum?: Test Ch	namber Hold Vacuum?:		
3μm 12937 Chamber Temp 47.0 9μm 13681 Breath Hose Temp 42	Snap Test ok? RFI Detected?	9	
Peripheral Tests Speaker (F5) Power LED Green (F6) Pump (F10) V Atmospheric Sensor Pressure 970 Correction Factor 1.00	Receipt paper Flow Sensor Tank Sensor	Is regulator close to Tank Sensor Value? Yes / No	
Sent to CMI (If needed) Reason Sent Date Sent Invoice Sent Date Returned			
Repairs needed: Cal 995 Connector Broke dusing evaluation	Repairs Made: replaced part		
		*	
	Ready for Bench Check?		
Inventory Process ☐ Bench Check Techn ☐ Linked Cylinder to Unit ☐ Mock Subject Test C ☐ Changed Unit Status ☐ Certificate of Calibration ☐ Bench Check Report Completed ☐ Service Overview Completed	Completed CMI Workorder CMI Certificate of Calibration ompleted CMI Service Evaluation Form	1	