

Dr. Jarrad Wagner, Ph.D. F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

ATTESTATION

I, the undersigned *Joshua Smith* (Director/Records Custodian), attest under oath pursuant to 12 O.S. §2902 this record was made, at or near the time of the occurrence of the matters set forth by or from information transmitted by a person having knowledge of those matters; was kept in the course of the regularly conducted business activity of the Board; and was made pursuant to the regularly conducted activity. The duplicate record provided is listed below:

1) The attached "Service Overview" is a true and accurate copy of maintenance performed on the Intoxilyzer 8000, serial number 80-003361, in compliance with the agency's Maintenance Policy and Procedures.

Signed

Joshua Smith

Name of Position: Director/Records Custodian

JORDAN FAITH THOMPSON NOTARY PUBLIC STATE OF OKLAHOMA COMMISSION NO. 25010889 EXPIRES 09-16-2029 State of Oklahoma

County of Oklahoma

Signed or attested before me on 10107 2025

by Joshua Smith.

(Signature of notarial officer)

Service Overview

3600 N. Martin Luther King Bldg #9 Oklahoma City, OK 73111 405-425-2460

BENCH CHECK DATE: 7/1/2025

SERIAL: 80-003361

DOCUMENTATION WITHIN SERVICE PACKET:

- BOT Certificate of Maintenance, Calibration and Operation
- Bench Check
- CMI Certificate of Calibration
- CMI Evaluation Form
- CMI Work Order
- BOT Dry Gas Certificate of Analysis
- Ilmo Specialty Gases Certificate of Analysis
- Mock Subject Test
- BOT Certificate of Analysis Reference Solution (4)
- Guth Laboratories Certificate of Analysis (4)



Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

CERTIFICATE of MAINTENANCE, CALIBRATION, and OPERATION

This certifies the Intoxilyzer 8000 referenced below, manufactured by CMI, Inc., a subsidiary of MPD, Inc. of Owensboro, Kentucky, is maintained by the Board of Tests for Alcohol and Drug Influence in accordance with 40 O.A.C. 30-1-3. The instrument reported below is hereby deemed to be in compliance with Board rules via linearity accuracy examination and operational/functional procedure evaluation also referred to as a "Bench Check". The Intoxilyzer 8000 appears on the National Highway Traffic Safety Administration's Conforming Products List of Evidential Breath Alcohol Measurement Devices (82 FR 50940 et seq) as an approved device. In accordance with 40 O.A.C. 25-1-2, this instrument is approved and certified for use in the State of Oklahoma as an automated evidential breath-alcohol measurement device by the Board of Tests pursuant to 47 O.S. § 759.

All analyses conducted by this instrument are considered precise, valid, and in accordance with 40 O.A.C. 30-1-3(b) Operating Procedure(s). A printed breath test operator affidavit containing the seal of the Board by this instrument confirms a calibration check successfully verified the instrument's calibration at the time of test. The printed breath test affidavit also serves as proof that safeguards in place to protect the breath samples, including but not limited to Ambient Air Interference, Radio Frequency Interference, Residual Mouth Alcohol, Deficient Samples, etc. were not triggered and the test is a true and accurate scientific measurement of the person's alcohol concentration.

- Maintenance was performed on the Intoxilyzer 8000, Serial No. 80-003361 on 7/1/2025 in accordance with 40 O.A.C. 30-1-3 and internal BOT maintenance policy and procedures.
- A dry gas canister appearing on the National Highway Traffic Safety Administration's Conforming Products List of Calibrating Units for Breath Alcohol Testers (77 FR 64588 et seq) and approved in accordance with 40 O.A.C. 25-1-3 and BOT maintenance policy was installed.
 - Manufacturer: ILMO Unit: 208 ppm / 105 L
 - Cylinder Lot No. Canister No.: 13225080A3-006
 - BOT analysis average test result: 0.081
 - Cylinder expiration date: 6/5/2027

3. In accordance with 40 O.A.C. 25-1-2.1, the approved reference standards and equipment were used:

Simulator Model			
Guth 2100	Guth 2100	Guth 2100	Guth 2100
Serial #			
DR3748	DR3694	DR2455	DR3028
Concentration			
.02	.04	.10	.20
LOT#	- A	770/07	C. Carlan
24150	25090	24390	25020
Manufactured Date		A	
4/16/2024	3/11/2025	10/29/2024	1/14/2025
Expiration Date			, united at
4/16/2026	3/11/2027	10/29/2026	1/14/2027
Solution Commission Date			
7/1/2025	7/1/2025	7/1/2025	7/1/2025

Completed by BOT Personnel:

Christopher Pape

Breath-Alcohol Testing Program Administrator

Title

Certified by:

Signature

Joshua Smith Name State Director Title

7/01/2025 UMENT TYPE		STRUMENT MODEL	80	1.779		NSTRUMENT SERIAL NUMBER	
ntoxilyzer		I-8000				3361	
		TES	T RESUL	TS			
Diagnost PASS							
Mouth Alcoho	37 17 2 2 3						
RFI Tes							
PASS		-	Duplica	ate Breath	Sampling Re	sults	
Abort Te		Val	ue	Test 1	Test 2	Control Test	
Improper Sam PASS		.0	2	0.018	0.018	0.019	
Interferent Det	ect Test	. 0	4	0.036	0.036	0.038	
PASS		.1	.0	0.095	0.096	0.098	
No Sample Giv		. 2	0	0.191	0.190	0.194	
Insufficient San				PA	SS		
Air Blank Contam							
Range Exceed							
		CERTIFYI	ING TECH	INICIAN		T. 20. T. S.	
I swear under penalt	ty of perjury	that in regards to	o the ab	ove listed	instrument, I	have complied	
the Intoxliyzer Bend	ch Checking	Procedure appro	ved by	the State o	f Oklahoma.	and a serio	

0000023020

PAPE, CHRISTOPHER

b)))

))))

))))

))))

))))

(iii)

DD))

Certificate of Calibration

This is to certify the calibration of INOXIIVZET ® serial number 80-003361, manufactured by CMI Inc., a subsidiary of MPD, Inc. of Owensboro, Kentucky, was tested and found to conform to the National Highway Traffic Safety Administration (NHTSA) Standard for Devices to Measure Breath Alcohol (Federal Register, Vol.58 No.179, pp 48705-48710, Sept. 17, 1993) for accuracy and precision. Reference materials are traceable through the National Institute of Standards and Technology (NIST) to the International System of Units (SI).

Date June. 14, 2025 Signed Christophen Mendall

INC.

((0)

(((0)

 $\|(u)\|_{L^{2}(\Omega)}$

((0)

((((

((0)

(((0))

(0)

((((i)

((0

(((0,

((()

 $\|\|_{L^{2}}$

(((0)

316 East 9th Street Owensboro, KY 42303 USA

Part No. 650517 Rev.A



This form MUST be completed and enclosed with item to be serviced.

*** Hazardous Material Warning! – DO NOT return gas cylinder with instrument! ***

Contact information: Name Chris Pape	Dhamas / 405					
Name Chris Pape	Dhama, 1, 405					
A Secretary of the Control of the Co	Pnone: (400) 425-2460				
Email: christopher.pape@bot.ok.gov	Customer # 731362	_ (contact Customer Service)				
Your Billing Address Agency Name: Board of Tests	Agency Name: Board	Your Shipping Address Agency Name: Board of Tests				
Address: 3600 N ML King, Bldg 9	Address: 3600 N ML					
	73111 City: Oklahoma City	State: OK Zip:73111				
Credit Card/PO #:	Contact Person: Chri					
Name on Card:		Contact Phone: (405) 425-2460				
Expiration Date: CVV	Email: christopher.pa	pe@bot.ok.gov				
Note: For instruments not under warranty \$43.50 (fuel cell), plus return shipping cos	그렇지 이 이 방송하다 이 경기에 있었다. 이 시간에 되었다. 그런 그렇게 되었다. 그런 그렇게 되었다. 그런	ared and I-200) or				
☑ I authorize all repairs up to \$500 to						
After repair, add "Certificate of Calibra	tion" for: 🔲 S-D2, S-D5, and I-	-800 (\$19)				
Authorized By:	1-200, I-240, I-300,	, I-400, and I-600 (\$39)				
2,490,900,000,000,000	reath Testing Admin Ship	item to:				
Name (Please Print)	Title CMI,	Inc. Service Dept.				



316 E 9th Street / Owensboro KY 42303 / USA Phone: 866-835-0690 Fax: 270-685-6268 **SERVICE WORK ORDER**

409098

Bill To:

OK BOT

DATE: 9/27/2024

PO Box 36307

Oklahoma City, OK 73136-2307

USA

Customer #: 731362 Phone: 405-425-2424 Fax: 405-425-2490 fax Ship To:

OK BOT PO Box 36307

Oklahoma City OK 73105-3204 USA

Phone: 405-425-2424

MODEL #: 002480OK

Serial Number: 80-003361

BILL CODE: Out of Warranty

EXTRA PARTS RCVD:

BH, KB, BATTERY, GAS SHELF

DESCRIPTION OF PROBLEM

UNIT WILL NOT POST

WORK PERFORMED:

Replaced parts listed; calibrated/final tested; 8133.16 software

PARTS USED

Seq. No.	Part	Description	Quantity
110	441169	COVER DUST,5/8in x 1/2in	1.00 EA
120	470145	CAP,PLASTIC,.25IDx.50LG,RED	1.00 EA
130	470154	CAP,PLSTC,.406IDx.438562,RED	1.00 EA
140	470218	SEAL, TAMPER EVIDENT, SERVICE	1.00 EA
. 150	690040	BATTERY,3VOLT,LITHIUM,CR2032	1.00 EA
160	230115	IC,SMD,DATAFLASH,4MB,TSOP-28	1.00 EA

LABOR / TESTING

Misc Code	Description	Hours
LABR	Service Repair Labor	
LCAL	Service - Calibration Adjust	0.50
ĻFT	Service - Cal / Final Test	1.00

Technician Name: Dewayne Varvel

SRV4090980001

Page: 1 of 3



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis - Dry Gas

ALPHA INSTRUMENT
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-002591
06/09/2025 10:06

Test	g/210L	Time
Air Blank	0.000	10:07
Cal Check	0.081	10:07
Air Blank	0.000	10:07
Cal Check	0.080	10:08
Air Blank	0.000	10:08
Cal Check	0.081	10:09
Air Blank	0.000	10:09
Cal Check	0.081	10:09
Air Blank	0.000	10:10
Cal Check	0.081	10:10
Air Blank	0.000	10:11
Cal Check	0.081	10:11
Air Blank	0.000	10:12
Cal Check	0.081	10:12
Air Blank	0.000	10:12
Cal Check	0.081	10:13
Air Blank	0.000	10:13
Cal Check	0.081	10:14
Air Blank	0.000	10:14
Cal Check	0.081	10:15
Air Blank	0.000	10:15
Cal Check Stats	3	
Average	0.0809	
Std Dev	0.0003	
Rel Std Dev(%)	0.3908	

0.08 g/210L = 208 PPM ÷ 2605*	Ethanol Gas (208 ppm/105 L)
Manufacturer Cert. No.: 17293	Cylinder Size: 105L
Cylinder No.: 006	Expiration Date: 06/05/2027
	Manufacturer Cert. No.: 17293

Anal. Toxicol. (1996) 20(6): pp. 484-491. [g/210L = PPM +2605].

Signature 06/09/2025

Name Christopher Pape Permit No.: 23020



7 Eastgate Dr. • P.O. Box 790 • Jacksonville, IL 62651-0790 217-245-2183 • Fax: 217-243-7634 • www.ilmoproducts.com

Certificate of Analysis

Certificate ID:

17293

Part #:

BAC105L080T

Cylinder Size:

105L

Lot Number:

13225080A3

Expiration:

6/5/2027

RCVD 5/28/2025 CP

0.080 BAC (For the calibration of instruments used to determine breath alcohol concentration)

Contents:

105 Liters @ 1000 psig 70°F (21°C)

Analytical

Analytical

Component:

Accuracy

Concentration:

(U, k=2): abs

Method:

Ethanol Nitrogen 208 ppm Balance

Reported

[5.2 ppm]

+/-0.002 BAC(G/210L) NDIR

CMI Inc.

Distributed by:

316 East Ninth Street Owensboro, KY 42303 Phone 866-835-0690 www.alcoholtest.com

*Reference Standard:

NIST Traceable Certified Reference Material - 260.6 µmol/mol Ethanol in Nitrogen - Serial No. GN0015020 Lot No. 00424E1 Store in dry area, away from sources of heat, ignition and direct sunlight. Do not allow storage area to exceed 52 °C (125 °F).

Specialty Gas Lab Tech

Issuance Date

Accreditation #61895

The calibration results within this certificate were obtained at the facility listed above using equipment and standards capable of producing analytical results traceable to NIST, and apply only to the items contained on this certificate. ILMO Products Company makes no warranty or representation as to the suitability of the use of any information provided for any particular purpose. The information use is at the sole discretion and risk of the user. Liability shall be limited to established replacement cost of this material or service. This certificate applies only to the items described and shall not be reproduced other than in full, without written approval from the issuing facility.

10/20

OFFICER'S IMPAIRED DRIVING AFFIDAVIT

								Section
TOTAL PROPERTY OF THE PROPERTY	RREST TIME 1500	COUNTY#	CITY	"	CITATIO	N#		
ARREST LOCATION BOT		CI	TY			OK	LAHC	MA
NORTH, GAYLON LEVER	NE	mark to be a fine	DATE OF	09/19	61 N		HEIGHT 51) 165
3600 MARTIN LUTHER	KING A	OKLAHOMA CIT	'Y		OK		73:	
T083454871		06/30/2		OK OK	A	RESTR	ICTIONS	ENDORSEMENTS M
VEHICLE MAKE MOCK	MOCK		TAG M	OCK		OK	N	N N

On the above date, time, and location, the above named person was arrested, and I had reasonable grounds to believe the person was driving or in actual physical control of a motor vehicle upon a public road, highway, street, tumpike, public parking lot, or other public place, or upon a private road, street, alley or lane which provides access to one or more single or multi-family dwelling within this state while under the influence of alcohol and/or other intoxicating substances as prohibited by law.

(Describe driving behavior or circumstances): MOCK SUBJECT TEST

Describe person's condition (odor, actions, etc.):

POST REPAIR MAINTENANCE, INSTALLED NEW CYLINDER

THE PERSON W	AS INFORMED OF THE IMPLIED CO	NSENT TEST REQUEST	
BREATH-ALCOHOL ANALYSIS RECORD, REPORT OF THE PERSON NAMED II	N SECTION 1 ABOVE, and LOG OF TE	ST AND MAINTENANCE RECORD.	Section 2
STANDARD Dry	O 0 3 3 6 1 BOARD TARGET VALUE 0 . 0 8 0	OF TESTS PERMIT NUMBER 23020	
MANUFACTURER ILMO	13225080A3	07/01/2025	DATE CYLINDER EXPIRES 06/05/2027
Cal Check 0.07 15:42 Air Blank 0.00 15:43	period was observed to the time of te I administered said E compliance with the	deprivation/observation prior to and continuing st administration. Breath-Alcohol Test in applicable operating Board of Tests.	Test Result BrAC g/210L 0.00

What happens next?

Section 3

- Upon receipt and review of the Officer's Impaired Driving Affidavit and any associated test results, if applicable, the State of Oklahoma may revoke your Oklahoma driving privileges.
- If you hold a commercial driver license or were driving a commercial motor vehicle, the State of Oklahoma may disqualify your privilege to operate a commercial motor vehicle.
 You will continue to have driving privileges, if otherwise eligible, unless or until an Order of Revocation and/or Disqualification is issued by the State of Oklahoma.
- The State of Oklahoma will send you the Order of Revocation and/or Disqualification to the address you have on file with Service Oklahoma.
 To confirm or change your address on file, contact Service Oklahoma.
- 5. The Order and Notice of Revocation will inform you of your requirements to regain driving privileges in Oklahoma.

in accordance with Titl	e 12 O.S. Section 426, "I state under	penalty or perj	ury under the laws of Oklahon	ia triat trie iore	going is true and correct.	Section 4
Date:	Place (location when signed):	resting officer:				
PAPE, CHR	STOPHER		BOARD OF TH	ESTS		
0000000484	AGENCY ADDRESS INTER-AGENCY		2 10 10 20 20	ZIP CODE 00000	PHONE 4054252460	
OTHER WITNESSES						Section :
1 NAME		TITLE	ADDRESS		PHONE	
2 NAME		TITLE	ADDRESS		PHONE	

11/20

OFFICER'S IMPAIRED DRIVING AFFIDAVIT

		*********					Section
The state of the s	SOO	55	CITY#	CITATION	V#		
ARREST LOCATION BOT		CI	TY		OK	LAHC	AM
SUBJECT NAME NORTH, GAYLON LEVERN	E		08/09/	1961 M		HEIGHT	165
ADDRESS 3600 MARTIN LUTHER K	ING A	OKLAHOMA CIT	Y	OK		731	
T083454871		06/30/2	019 OK	A	RESTRI	CTIONS	ENDORSEMENTS M
PEHICLE MAKE MOCK	MOCK		MOCK		OK	N	N PLACARD

On the above date, time, and location, the above named person was arrested, and I had reasonable grounds to believe the person was driving or in actual physical control of a motor vehicle upon a public road, highway, street, tumpike, public parking lot, or other public place, or upon a private road, street, alley or lane which provides access to one or more single or multi-family dwelling within this state while under the influence of alcohol and/or other intoxicating substances as prohibited by law.
(Describe driving behavior or circumstances):

MOCK SUBJECT TEST

Describe person's condition (odor, actions, etc.):

POST REPAIR MAINTENANCE, INSTALLED NEW CYLINDER

		THE PERSON W	AS INFORMED OF	THE IMPLIED CO	ONSENT TEST REQUEST	
BREATH-ALCOHOL ANALYSIS	S RECORD, REPORT OF THE P	PERSON NAMED II	SECTION 1 ABOV	E, and LOG OF 1	EST AND MAINTENANCE RECORD.	Section 2
I-8000		SERIAL N	UMBER 003361	BOARD	OF TESTS	
STANDARD Dry			0.080	JÉ	PERMIT NUMBER 23020	
MANUFACTURER ILMO			13225	080A3	07/01/2025	06/05/2027
1515 0	7/01/2025				r) deprivation/observation prior to and continuing	
RESULT TABLE Test	g/210L	Time	W. C.		est administration.	OF THE STATE
Air Blank Subject Test Breath Volume Air Blank Wait Air Blank Subject Test Breath Volume Air Blank Cal Check Air Blank	0.00 0.00 2.261 LITERS 0.00 0.00 0.00 2.140 LITERS 0.00 0.07	0.00 15:38 0.00 15:38 2.261 LITERS 0.00 15:39 15:41 0.00 15:41 0.00 15:42 2.140 LITERS 0.00 15:42 0.00 15:42		tered said ace with the edure of the educe o	Test Result BrAC g/210L	
1 () ()			BADGE # 0000000 AGENCY BOARD C	TALLET TA	0000023020	0.00

What happens next?

Section 3

- Upon receipt and review of the Officer's Impaired Driving Affidavit and any associated test results, if applicable, the State of Oklahoma may revoke your Oklahoma driving privileges.
- If you hold a commercial driver license or were driving a commercial motor vehicle, the State of Oklahoma may disqualify your privilege to operate a commercial motor vehicle.
 You will continue to have driving privileges, if otherwise eligible, unless or until an Order of Revocation and/or Disqualification is issued by the State of Oklahoma.
- The State of Oklahoma will send you the Order of Revocation and/or Disqualification to the address you have on file with Service Oklahoma.
 To confirm or change your address on file, contact Service Oklahoma.
- 5. The Order and Notice of Revocation will inform you of your requirements to regain driving privileges in Oklahoma.

In accordance with Title 12 O.S. Section 426, "I state under penalty of perjury under the laws of Oklahoma that the foregoing is true and correct."		
Date: Place (location when signed):_	Signature of arresting officer:	
PAPE, CHRISTOPHER	BOARD OF TESTS	
BADGE AGENCY ADDRESS TNTER - AGENCY	ZIP CODE PHONE 4054252460	



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis - Reference Solution

ALPHA INSTRUMENT	
Intoxilyzer - Alcoho	ol Analyzer
Model 8000	SN 80-002591
01/15/2025	09:13

Test	g/210L	Time
Air Blank	0,000	09:14
Cal Check	0.021	09:14
Air Blank	0.000	09:15
Cal Check	0.021	09:16
Air Blank	0.000	09:16
Cal Check	0.021	09:17
Air Blank	0.000	09:17
Cal Check	0.021	09:18
Air Blank	0.000	09:19
Cal Check	0.020	09:19
Air Blank	0.000	09:20
Cal Check	0.021	09:21
Air Blank	0.000	09:21
Cal Check	0.021	09:22
Air Blank	0.000	09:22
Cal Check	0.021	09:23
Air Blank	0.000	09:24
Cal Check	0.021	09:24
Air Blank	0.000	09:25
Cal Check	0.020	09:26
Air Blank	0.000	09:26
Cal Check Stats		
Average	0.0208	
Std Dev	0.0004	
Rel Std Dev(%)	2.0270	

Manufacturer: Guth	Manufactured Date: 4/16/2024	Concentration: 0.020
Lot No.: 24150	Expiration Date: 4/16/2026	Equipment Model: Guth 2100
prescribed by the State Director of Tests proceed to 40 OAC 25-1-2.1. Approved breath-alco	Infrared Spectroscopy utilizing the Intoxilyzer 8000 and oursuant to 40 OAC 30-1-3. The above NIST traceable re shol measurement equipment and reference solutions. Ed with for Breath Alcohol Testers (72 FR 34747), published	oference solution meets the requirements for use pursuant quipment used appears on the current or supplemented



RCVD 12/3/24 CP

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 24080 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on February 14, 2024, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.2419% (w/vol) ethyl alcohol. The expiration date for this lot number is February 13, 2026 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.200 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN03132302 whose values are traceable to NIST.



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis - Reference Solution

ALPHA INSTRUMENT 2 Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-004923 07/01/2025 11:06

Test	g/210L	Time
Air Blank	0.000	11:07
Cal Check	0.037	11:07
Air Blank	0.000	11:08
Cal Check	0037	11:09
Air Blank	0.000	11:09
Cal Check	0.038	11:10
Air Blank	0.000	11:10
Cal Check	0.038	11:11
Air Blank	0.000	11:12
Cal Check	0.037	11:12
Air Blank	0.000	11:13
Cal Check	0.038	11:14
Air Blank	0.000	11:14
Cal Check	0.038	11:15
Air Blank	0.000	11:15
Cal Check	0.038	11:16
Air Blank	0.000	11:17
Cal Check	0.038	11:17
Air Blank	0.000	11:18
Cal Check	0.038	11:19
Air Blank	0.000	11:19
Cal Check Stats		
Average	0.0377	
Std Dev	0.0004	
Rel Std Dev(%)	1.2812	

Manufacturer: Guth	Manufactured Date:	3/11/2025	Concentration: 0.040
Lot No.: 25090	Expiration Date: 3/	11/2027	Equipment Model: Guth 2100
prescribed by the State Director of Test to 40 OAC 25-1-2.1. Approved breath-a	s pursuant to 40 OAC 30-1-3. The cohol measurement equipment an	above NIST traceable refe d reference solutions. Equ 2 FR 34747), published by Date: 7/1/202	
Name Christopher Pape		Permit No.: 2302	20



RCVV) 5/30/25

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 25090 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on March 12, 2025, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.0480% (w/vol) ethyl alcohol. The expiration date for this lot number is March 11, 2027 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.040 g/210L+/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN06032102 whose values are traceable to NIST.



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis - Reference Solution

ALPHA INSTRUMENT Intoxilyzer - Alcohol Analyzer Model 8000

SN 80-002591

01/15/2025

09:30

Test · ·	g/210L		Time
Air Blank	0.000		09:31
Cal Check	0.100		09:31
Air Blank	0.000		09:32
Cal Check	0.100		09:33
Air Blank	0.000		09:33
Cal Check	0.099		09:34
Air Blank	0.000		09:35
Cal Check	0.099		09:35
Air Blank	0.000		09:36
Cal Check	0.099		09:37
Air Blank	0.000		09:37
Cal Check	0.099	1.6	09:38
Air Blank	0.000		09:38
Cal Check	0.100		09:39
Air Blank	0.000		09:40
Cal Check	0.099		09:40
Air Blank	0.000		09:41
Cal Check	0.099		09:42
Air Blank	0.000		09:42
Cal Check	0.100		09:43
Air Blank	0.000		09:43
Cal Check Stats			
Average	0.0994		
Std Dev	0.0005		
Rel Std Dev(%)	0.5195		

Manufacturer: Guth	Manufactured Date: 10/29/2024	Concentration: 0.100
Lot No.: 24390	Expiration Date: 10/29/2026	Equipment Model: Guth 2100
prescribed by the State Director of Tests ports 40 OAC 25-1-2.1, Approved breath-alcol	Infrared Spectroscopy utilizing the Intoxilyzer 8000 and w ursuant to 40 OAC 30-1-3. The above NIST traceable ref hol measurement equipment and reference solutions. Equ its for Breath Alcohol Testers (72 FR 34747), published b	erence solution meets the requirements for use pursuant uipment used appears on the current or supplemented
Signature / August /	Date: 1/15/20	25
Name Christopher Pape	Permit No.: 2302	20

RCVP 12/3/24

Certified Alcohol Reference Solution for Simulator

cp

Random Samples of Lot Number 24390 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on October 30, 2024, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1205% (w/vol) ethyl alcohol. The expiration date for this lot number is October 29, 2026 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN03072301 whose values are traceable to NIST.



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis - Reference Solution

ALPHA INSTRUMENT 2 Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-004923 06/06/2025 11:23

Test	g/210L	Time
Air Blank	0.000	11:24
Cal Check	0.200	11:24
Air Blank	0.000	11:25
Cal Check	0.199	11:26
Air Blank	0.000	11:26
Cal Check	0.199	11:27
Air Blank	0.000	11:27
Cal Check	0.198	11:28
Air Blank	0.000	11:29
Cal Check	0.199	11:29
Air Blank	0.000	11:30
Cal Check	0.199	11:31
Air Blank	0.000	11:31
Cal Check	0.199	11:32
Air Blank	0.000	11:33
Cal Check	0.197	11:33
Air Blank	0.000	11:34
Cal Check	0.199	11:35
Air Blank	0.000	11:35
Cal Check	0.198	11:36
Air Blank	0.000	11:36
Cal Check Stats		
Average	0.1987	
Std Dev	0.0008	
Rel Std Dev(%)	0.4143	

Manufacturer: Guth	Manufactured Date: 1/14/2025	Concentration: 0.200	
Lot No.: 25020	Expiration Date: 1/14/2027	Equipment Model: Guth 2100	
prescribed by the State Director of Tes to 40 OAC 25-1-2.1. Approved breath-	via Infrared Spectroscopy utilizing the Intoxilyzer 8000 and vits pursuant to 40 OAC 30-1-3. The above NIST traceable realcohol measurement equipment and reference solutions. Economics for Breath Alcohol Testers (72 FR 34747), published	rerence solution meets the requirements for use pursuant quipment used appears on the current or supplemented	
Signature Million	Date: 6/6/202	Date: 6/6/2025	
Name Christopher Pape	Permit No.: 230	20	



RCVI) 5/30/25

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 25020 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on January 15, 2025, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.2409% (w/vol) ethyl alcohol. The expiration date for this lot number is January 14, 2027 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.200 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN03132302 whose values are traceable to NIST.

Location Returned From Lawton PD 2	<u>Unit Serial #</u> 80-003361	Date: 8/13/24
Inventory Process	Old Cylinder #: 0152408	0A1-049
Check in Procedure Connect unit to ethernet connection 172.019.0 Download any stored records Change unit location to Serial number Update Unit Date and time Previous time: New Time:	Notes from the Field	Issues in the Field
Lab Evaluation Breath Hose Hold Vacuum?:	Test Chamber Hold Vacuum?:	
DVM 3μm Chamber Temp 9μm Breath Hose Temp Peripheral Tests	Snap Test ok? RFI Detected?	
Speaker (F5) Power LED Gre Gas Regulator (F11) Pump Atmospheric Sensor		Is regulator close to Tank Sensor Value?
	n Sent Unit will not	
Repairs needed:	Repairs Made:	
☐ Linked Cylinder to Unit ☐ Mock Sub☐ Changed Unit Status ☐ Certificate	oject Test Completed	ete if unit sent to CMI