

Dr. Jarrad Wagner, Ph.D. F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

ATTESTATION

I, the undersigned <u>Joshua Smith</u> (Director/Records Custodian), attest under oath pursuant to 12 O.S. §2902 (11)(a) this record was made, at or near the time of the occurrence of the matters set forth by or from information transmitted by a person having knowledge of those matters; was kept in the course of the regularly conducted business activity of the Board; and was made pursuant to the regularly conducted activity. The duplicate record provided is listed below:

1) The attached "Service Overview" is a true and accurate copy of maintenance performed on the Intoxilyzer 8000, serial number 80-004423, in compliance with the agency's Maintenance Policy and Procedures.

This packet consists of 17 pages.

Signed _____

Name of Position: Director/Records Custodian

Date of Attestation: 07/12/2022

Service Overview



3600 N. Martin Luther King Bldg #9 Oklahoma City, OK 73111 405-425-2460

BENCH CHECK DATE: 07/12/2022

SERIAL: 80-004423

DOCUMENTATION WITHIN SERVICE PACKET:

- BOT Certificate of Maintenance, Calibration and Operation
- Bench Check
- . BOT Dry Gas Certificate of Analysis
- Ilmo Specialty Gases Certificate of Analysis
- Mock Subject Test
- BOT Certificate of Analysis on Certified Solutions (4)
- Guth Laboratories Certificate of Analysis (4)



Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

CERTIFICATE of MAINTENANCE, CALIBRATION, and OPERATION

This certifies the Intoxilyzer 8000 referenced below, manufactured by CMI, Inc., a subsidiary of MPD, Inc. of Owensboro, Kentucky, is maintained by the Board of Tests for Alcohol and Drug Influence in accordance with 40 O.A.C. 30-1-3. The instrument reported below is hereby deemed to be in compliance with Board rules via linearity accuracy examination and operational/functional procedure evaluation also referred to as a "Bench Check". The Intoxilyzer 8000 appears on the National Highway Traffic Safety Administration's Conforming Products List of Evidential Breath Alcohol Measurement Devices (82 FR 50940 et seq) as an approved device. In accordance with 40 O.A.C. 25-1-2, this instrument is approved and certified for use in the State of Oklahoma as an automated evidential breath-alcohol measurement device by the Board of Tests pursuant to 47 O.S. § 759.

All analyses conducted by this instrument are considered precise, valid, and in accordance with 40 O.A.C. 30-1-3(b) Operating Procedure(s). A printed breath test operator affidavit containing the seal of the Board by this instrument confirms a calibration check successfully verified the instrument's calibration at the time of test. The printed breath test affidavit also serves as proof that safeguards in place to protect the breath samples, including but not limited to Ambient Air Interference, Radio Frequency Interference, Residual Mouth Alcohol, Deficient Samples, etc. were not triggered and the test is a true and accurate scientific measurement of the person's alcohol concentration.

- Maintenance was performed on the Intoxilyzer 8000, Serial No. <u>80-004423</u> on <u>07/12/2022</u> in accordance with 40 O.A.C. 30-1-3 and internal BOT maintenance policy and procedures.
- A dry gas canister appearing on the National Highway Traffic Safety Administration's Conforming Products List of Calibrating Units for Breath Alcohol Testers (77 FR 64588 et seq) and approved in accordance with 40 O.A.C. 25-1-3 and BOT maintenance policy was installed.
 - Manufacturer: ILMO Unit: 208 ppm / 105 L
 - Cylinder Lot No. Canister No.: 07622080A2-039
 - BOT analysis average test result: 0.080
 Cylinder expiration date: 05/05/2024
- In accordance with 40 O.A.C. 25-1-2.1, the approved reference standards and equipment were used:

Simulator Model			
Guth 2100	Guth 2100	Guth 2100	Guth 2100
Serial #		- W-W-6-	
DR2834	DR3694	DR2455	DR3470
Concentration			4 600
.02	.04	.10	.20
LOT#			
21120	21070	22080	22050
Manufactured Date	- C C C C C C C C.		
4/7/2021	3/1/2021	3/7/2022	2/7/2022
Expiration Date			i i i i i i i i i i i i i i i i i i i
4/7/2023	3/1/2023	3/7/2022	2/7/2024
Solution Commission Date			
06/27/2022	06/27/2022	06/27/2022	07/06/2022
177 (177)			The state of the s

27.2971/2	06/27/2022	06/27/2022	06/27/2022	
	06/21/2022	06/2/12022	06/27/2022	
Completed by BC	T Personnel:			
Lity Co	Christopher Pape	Proofb Alcohol Tost	ing Program Administrator	
Signature/Name	Christopher Pape	Title	ing Program Administrator	

Certified by:

Signature Joshua Smith

State Director Title

STATE OF OKLAHOMA - BOARD OF TESTS FOR ALCOHOL AND DRUG INFLUENCE CMI INTOXILYZER BENCH CHECK PROCEDURE REPORT

Test Test Pass Duplicate Breath Sampling Results	07/12/2022	10:3		80-00442	23	
Diagnostics	strument type Intoxilyzer		I - 8000			
Mouth Alcohol Test			TEST	RESULTS		
No Sample Test	PASS Mouth Alcoho PASS RFI Test	ol Test t		Dundingto Brooth	Sometice D	
Improper Sample Test	Abort Te	st		Duplicate Breath	Sampling Re	esults
PASS	PASS		Valu	e Test 1	Test 2	Control Test
Interferent Detect Test				0.010	0.017	0.010
PASS No Sample Given Test			. 02	0.018	0.017	0.019
10 0.098 0.096 0.099			. 04	0.037	0.037	0.038
PASS Insufficient Sample Test PASS Air Blank Contamination Test PASS Range Exceeded Test			.10	0.098	0.096	0.099
PASS Air Blank Contamination Test PASS Range Exceeded Test		en Test	.20	0.196	0.196	0.197
PASS Range Exceeded Test		ple Test	,	P	ASS	
		nation Test				
		ed Test				
CERTIFYING TECHNICIAN			CERTIFYIN	IG TECHNICIAN		
swear under penalty of perjury that in regards to the above listed instrument, I have com the Intoxliyzer Bench Checking Procedure approved by the State of Oklahoma.						

CMI Form OKBE1 01/2019

0000023020

CHRISTOPHER

PAPE,



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis - Dry Gas

ALPHA INSTRUMENT

Intoxilyzer - Alcohol Analyzer

Model 8000

SN 80-002591

07/05/2022

12:13

Test	g/210L	Time
Air Blank	0.000	12:14
Cal Check	0.080	12:14
Air Blank	0.000	12:15
Cal Check	0.080	12:15
Air Blank	0.000	12:16
Cal Check	0.080	12:16
Air Blank	0.000	12:16
Cal Check	0.080	12:17
Air Blank	0.000	12:17
Cal Check	0.080	12:18
Air Blank	0.000	12:18
Cal Check	0.080	12:18
Air Blank	0.000	12:19
Cal Check	0.080	12:19
Air Blank	0.000	12:20
Cal Check	0.080	12:20
Air Blank	0.000	12:21
Cal Check	0.080	12:21
Air Blank	0.000	12:21
Cal Check	0.080	12:22
Air Blank	0.000	12:22
Cal Check Stats		
Average	0.0800	
Std Dev	0.0000	
Rel Std Dev(%)	0.0000	

Manufacturer: ILMO	0.08 g/210L = 208 PPM ÷ 2605* Manufacturer Cert. No.: 14727	Cylinder Size: 105L
Lot No.: 07622080A2	Cylinder No.: 39	Expiration Date: 05/05/2024
		reference standard meets the requirements for use as a
calibrating unit in accordance with 47 O.S. testers published by the U.S. Department Anal. Toxicol. (1996)-20(6): pp. 484-491.	. 759, 40 OAC 25-1-3, and appears on the current confo of Transportation in the Federal Register [77 FR 64588	reference standard meets the requirements for use as a rming products list of calibrating units for breath alcohol et seq]. * See Dubowski, K.M. and Essary N., Journal of
calibrating unit in accordance with 47 O.S. testers published by the U.S. Department	. 759, 40 OAC 25-1-3, and appears on the current confo of Transportation in the Federal Register [77 FR 64588	rming products list of calibrating units for breath alcohol et seq]. * See Dubowski, K.M. and Essary N., <i>Journal of</i>



7 Eastgate Dr. • P.O. Box 790 • Jacksonville, IL 62651-0790 217-245-2183 • Fax: 217-243-7634 • www.ilmoproducts.com

Certificate of Analysis

Certificate ID:

14727

Part #:

BAC105L080T

Cylinder Size:

105L

Lot Number:

07622080A2

Expiration:

5/5/2024

12/12/22 RCVP

0.080 BAC (For the calibration of instruments used to determine breath alcohol concentration)

Contents:

105 Liters @ 1000 psig 70°F (21°C)

Analytical

Reported

Accuracy

Analytical

Component:

Concentration:

(U, k=2):

Method:

Ethanol Nitrogen

208 ppm Balance +/-0.002 BAC(G/210L) NDIR

[5.2 ppm]

Distributed by:

CMI Inc.

316 East Ninth Street Owensboro, KY 42303 Phone 866-835-0690 www.alcoholtest.com

*Traceable to:

Certified Reference Material - 251.8 µmol/mol

Ethanol in Nitrogen - Serial No. CC198518 Lot No. 111318E1

Store in dry area, away from sources of heat, ignition and direct sunlight. Do not allow storage area to exceed 52 °C (125 °F).

Specialty Gas Lab Tech

04-28-2022

Issuance Date

Accreditation #61895

The calibration results within this certificate were obtained using equipment and standards capable of producing analytical results traceable to NIST, and apply only to the items contained on this certificate. ILMO Products Company makes no warranty or representation as to the suitability of the use of any information provided for any particular purpose. The information use is at the sole discretion and risk of the user. Liability shall be limited to established replacement cost of this material or service.

7/17

OFFICER'S AFFIDAVIT AND NOTICE OF

REVOCATION/DISQUALIFICATION Section 1 ARREST TIME COUNTY II 07/12/2022 1115 55 ARREST LOCATION BOT OKLAHOMA SUBJECT NAME DATE OF BIRTH HEIGHT WEIGHT 08/09/1961 NORTH, GAYLON LEVERNE M 510 165 ADDRESS ZIP CODE 3600 MARTIN LUTHER KING A OKLAHOMA CITY OK 73133 EXPIRATION DATE DRIVER LICENSE / RESTRICTIONS ENDORSEMENTS STATE CLASS T083454871 06/30/2019 OK A M VEHICLE MAKE MODEL PLACARD MOCK MOCK MOCK OK N N On the above date, time, and location, the above named person was arrested, and I had reasonable grounds to believe the person was driving or in actual physical control of a motor vehicle upon a public road, highway, street, turnpike, public parking lot, or other public place, or upon a private road, street, alley or lane which provides access to one or more single or multi-family dwelling within this state while under the influence of alcohol and/or other intoxicating substances as prohibited by law (Describe driving behavior or circumstances): MOCK SUBJECT TEST Describe person's condition (odor, actions, etc.): ROUTINE MAINTENANCE, INSTALLED NEW CYLINDER THE PERSON WAS INFORMED OF THE IMPLIED CONSENT TEST REQUEST

BREATH-ALCOHOL ANALY	SIS RECORD, REPORT OF THE	PERSON NAMED	IN SECTION 1 ABO	E, and LOG OF T	EST AND	MAINTENANCE RECORD.	Section 3
I-8000			NUMBER	SPECIALIST	OF	mpama	
		80	-004423	BOARD		TESTS	
STANDARD Dry			0.080		100,000	T NUMBER 020	
MANUFACTURER ILMO			07622	080A2	DATE O	/12/2022	DATE CYLINDER EXPIRES 05/05/2024
0BSERVATION START TO	EST DATE 07/12/2022				r) depr	ivation/observation	
RESULT TABLE Test	g/210L	Time		the time of to		o and continuing ninistration.	OF THE STATE
Air Blank Subject Test	0.00	12:09 12:10	Charles and the control of	and when he was been		-Alcohol Test in	
Breath Volum	e 1.824 LITERS	12:10	Daniel And Lide Common to a second			icable operating ard of Tests.	
Wait	0.00	12:10	, /	cuare of th	0	nu or resis.	
Air Blank	0.00	12:12	1/1/10	Ma //	10		
Subject Test Breath Volum	0.00 e 2.156 LITERS	12:12	SIGNATURE C	FOPERATOR	-	~	1907
Air Blank	0.00	12:13	NAME				
Cal Check	0.07	12:13	The state of the s	ID T GEODIN			Test Result
Air Blank	0.00	12:14	PAPE, C. BADGE # 000000	HRISTOPHE 1484	PERMIT	00023020	BrAC g/210L 0 . 0 0
			BOARD O	F TESTS			

EFFECTIVE THIRTY (30) DAYS FROM DATE SERVED. Notice of Revocation is given that your driver's license is revoked or denied for a period of one hundred eighty (180) Section 3 days or more, if you refused or failed the state's test(s). In addition, your commercial driving privilege in this state will be disqualified if you refused or failed the state's test(s) (operaling a CMV or CDL holder) OR a test result of .04 or more while operating a commercial motor vehicle. You may appeal any Departmental action against your driver license as a result of this arrest by filing a petition in the District Court of the County in which you were arrested. Your petition must be filed within thirty (30) days of the date listed below.

You may be eligible to participate in the Impaired Driver Accountability Program (IDAP) administered by the Department of Public Safety (DPS). Completion of IDAP may prevent a revocation as a result of this arrest from appearing on your driving record. Participation in IDAP may reduce the amount of fees you will be required to pay to the State. You must request IDAP within thirty (30) days of the date listed below. Mall the request for IDAP to the Department of Public Safety at P.O. Box 11415, Oklahoma City, DK 73136.

SERVING OF NOTICE: I personally hand delivered a copy of the above containing the Notice of Revocation and the Notice of Department Impaired Driving Accountability Program to the person named in Section 1 above.

Date served:	Signature of officer:			Badge #:			
In accordance with Title 1	2 O.S. Section 426, "I state under per	nalty of perj	ury under the laws of Oklahor	na that the fore	going is true and correct."	Section 4	
Date:	Place (location when signed):		Signature	of arresting of	ficer:		
PAPE, CHRI	STOPHER		BOARD OF T	ESTS			
0000000484	AGENCY ADDRESS INTER-AGENCY			OOOE	PHONE 4054252460		
OTHER WITNESSES:						Section 5	
1 NAME		TITLE	ADDRESS	7	PHONE		
2 NAME		TITLE	ADDRESS		PHONE		

8/17

OFFICER'S AFFIDAVIT AND NOTICE OF

RREST DATE		ARREST TIME		ICOUNTY		ICATIO	ICITATI	ON #		Sect
07/12/2022		1115		55		7000	1 3/10/			
RREST LOCATION BOT		4111			CITY			COUN	LAHOI	ΛA
JBJECT NAME		1434				ATE OF BIRTH		X	HEIGHT	WEIGHT
	LON LEVE	RNE	Loiru			08/09/19		M	510	165
DRESS 3600 MARTIN	LUTHER	KING A	OKL	AHOMA C	TTY		OK		731	
IVER LICENSE #	- / -	11000	1,27	EXPIRATION D	DATE	STATE	CLASS			NDORSEME
T083454871		MODEL		06/30	1/201	19 OK	A	ICTATE.	low.	M
10CK		MOC	K			MOCK		OK	N F	N
the above date, time, and loc tway, street, turnpike, public prence of alcohol and/or other scribe driving behavior or circ MOCK SUBJECT	parking lot, or other put intoxicating substance: :umstances):	blic place, or upon a priv	and I had reasonable gro vate road, street, alley or t	unds to believe the ane which provide:	e person wa s access to	s driving or in actual pl one or more single or	hysical conti multi-family	ol of a motor v dwelling within	ehicle upon a this state wh	public road, le under the
cribe person's condition (odo	r. actions, etc.):									
EATH-ALCOHOL ANALYSIS UIPMENT	RECORD, REPORT (OF THE PERSON NAM	AL NUMBER	E, and LOG OF TE	EST AND N	AINTENANCE RECO	RD.			Secti
ANDARD		81	0 - 004423 TARGET VALUE	BOARD	PERMIT	TESTS				
ry			0.080		230	20				
NUFACTURER LMO			07622	00070	DATE CY	LINDER INSTALLED		DATE CYLIND	ER EXPIRES	4
COLUMN TO SERVICE AND ADDRESS OF THE PARTY ADDRESS OF THE PARTY AND ADD	DATE				-	12/2022		05/05	/202	4
.130 07 SULT TABLE	7/12/202: g/210L	2 Time	period w		prior to	ration/observati and continuing nistration.		18	THE ST	ATE
Air Blank Subject Test Breath Volume Air Blank Wait Air Blank Subject Test Breath Volume	0.00 0.00 1.824 LITE 0.00 0.00 0.00 2.156 LITE	12:10 12:12 12:12 12:12	compliand	ce with the	applie	Alcohol Tes cable operati d of Tests.	0.300	SHEAT SEAL	1907	
ir Blank al Check	0.00	12:13 12:13	NAME	A consideration			_		Test Res	alt
Air Blank	0.00	12:14	BADGE # 0000000 AGENCY	RISTOPHE 484	PERMIT A	0023020		В	rAC g/2	10L
			BOARD OF	TESTS				_		
FECTIVE THIRTY (30) DAY s or more, if you refused or erating a CMV or CDL holds a result of this arrest by filing a may be eligible to partic ocation as a result of this uest IDAP within thirty (36 RVING OF NOTICE: countability Program to te served: accordance with Title 1	failed the state's tester) OR a test result of g a petition in the Disipate in the Impaired arrest from appear in days of the date lill personally hand on the person named 2 O.S. Section 42	(s). In addition, your r .04 or more while op- trict Court of the Cour- d Driver Accountabil ing on your driving r isted below. Mail the delivered a copy ed in Section 1 ab Signature of offic 26, "I state under p	ocation is given that yo commercial driving privi- erating a commercial mity in which you were an ity Program (IDAP) ad- ecord. Participation is erequest for IDAP to to of the above contain ove.	ur driver's license lege in this state otor vehicle. You rested. Your peli ministered by th n IDAP may redu he Department of ning the Notice	will be discount of the period of the among the among the among the among Public See of Revenues of Okla	ualified if you refuse all any Departmental be filed within thirty (3 ment of Public Safet bount of fees you wisafety at P.O. Box 1 ocation and the Norma that the for	d or failed to action aga action aga action aga action aga action	he state's testinst your drive the date lister completion o red to pay to homa City, O Department Badge #:_	(s) Ir license If below. If IDAP may the State. K 73136.	You must
s or more, if you refused or arrating a CMV or CDL holds result of this arrest by filing may be eligible to participation as a result of this usest IDAP within thirty (30 RVING OF NOTICE: countability Program to e served:	failed the state's tester) OR a test result of g a petition in the Disipate in the Impaired arrest from appear in days of the date lill personally hand on the person named 2 O.S. Section 42	(s). In addition, your I .04 or more while op- trict Court of the Coun- d Driver Accountabiling on your driving ra- isted below. Mail the delivered a copy- ed in Section 1 ab- Signature of office	ocation is given that yo commercial driving privi- erating a commercial mity in which you were an ity Program (IDAP) ad- ecord. Participation is erequest for IDAP to to of the above contain ove.	ur driver's license lege in this state otor vehicle. You rested. Your peli ministered by th n IDAP may redu he Department of ning the Notice	will be disc may appe lion must be the Departmice the am of Public se e of Reve s of Okla	ualified if you refuse al any Departmental be filed within thirty (2 nent of Public Safet ount of fees you w lafety at P.O. Box 1 ocation and the N	d or failed to action aga action aga action aga action aga action	he state's testinst your drive the date lister completion o red to pay to homa City, O Department Badge #:_	(s) Ir license If below. If IDAP may the State. K 73136.	prevent a You must Driving



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis - Reference Solution

ALPHA INSTRUMENT 2 Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-004923 03/08/2022 11:33

Test	g/210L	Time
Air Blank	0.000	11:34
Cal Check	0.022	11:34
Air Blank	0.000	11:35
Cal Check	0.021	11:36
Air Blank	0.000	11:36
Cal Check	0.021	11:37
Air Blank	0.000	11:37
Cal Check	0.021	11:38
Air Blank	0.000	11:38
Cal Check	0.021	11:39
Air Blank	0.000	11:40
Cal Check	0.021	11:40
Air Blank	0.000	11:41
Cal Check	0.020	11:42
Air Blank	0.000	11:42
Cal Check	0.021	11:43
Air Blank	0.000	11:43
Cal Check	0.021	11:44
Air Blank	0.000	11:45
Cal Check	0.020	11:45
Air Blank	0.000	11:46
Cal Check Stats		
Average	0.0209	
Std Dev	0.0005	
Rel Std Dev(%)	2.7158	

Manufacturer: Guth	Manufactured Date: 04/07/2021	Concentration: 0.02	
Lot No.: 21120	Expiration Date: 04/07/2023	Equipment Model: Guth 2100	
to 40 OAC 25-1-2.1. Approved breath- Conforming Products List of Calibrating	ts pursuant to 40 OAC 30-1-3. The above NIST traceable ref alcohol measurement equipment and reference solutions. Eq g Units for Breath Alcohol Testers (72 FR 34747), published b	uipment used appears on the current or supplemented y the National Highway Traffic Safety Administration.	
Name Charles Office	03/08/2		
Name Christopher Pape	Permit No.: 2302	Permit No.: 23020	



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis - Reference Solution

ALPHA INSTRUMENT 2 Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-004923 03/08/2022 11:53

Test	g/210L	Time
Air Blank	0.000	11:53
Cal Check	0.041	11:54
Air Blank	0.000	11:54
Cal Check	0.041	11:55
Air Blank	0.000	11:56
Cal Check	0.042	11:56
Air Blank	0.000	11:57
Cal Check	0.042	11:58
Air Blank	0.000	11:58
Cal Check	0.042	11:59
Air Blank	0.000	11:59
Cal Check	0.041	12:00
Air Blank	0.000	12:01
Cal Check	0.041	12:01
Air Blank	0.000	12:02
Cal Check	0.041	12:02
Air Blank	0.000	12:03
Cal Check	0.041	12:04
Air Blank	0.000	12:04
Cal Check	0.042	12:05
Air Blank	0.000	12:05
Cal Check Stats		
Average .	0.0414	
Std Dev	0.0005	
Rel Std Dev(%)	1.2473	

Manufacturer: Guth	Manufactured Date: 03/01/2021	Concentration: 0.04	
Lot No.: 21070	Expiration Date: 03/01/2023	Equipment Model: Guth 2100	
prescribed by the State Director of Tests to 40 OAC 25-1-2.1. Approved breath-ald Conforming Products List of Calibrating I	la Infrared Spectroscopy utilizing the Intoxilyzer 8000 and w pursuant to 40 OAC 30-1-3. The above NIST traceable rel cohol measurement equipment and reference solutions. Eq Jnits for Breath Alcohol Testers (72 FR 34747), published b	ference solution meets the requirements for use pursuant	
Signature	Date: 03/08/2	022	
Name Christopher Pape	Permit No.: 230	Permit No.: 23020	



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis - Reference Solution

ALPHA INSTRUMENT Intoxilyzer - Alcohol Analyzer Model 8000

SN 80-002591

06/27/2022

11:40

Test	g/210L	Time
Air Blank	0.000	11:40
Cal Check	0.100	11:41
Air Blank	0.000	11:41
Cal Check	0.100	11:42
Air Blank	0.000	11:43.
Cal Check	0.100	11:43
Air Blank	0.000	11:44
Cal Check	0.100	11:45
Air Blank	0.000	11:45
Cal Check	0.100	11:46
Air Blank	0.000	11:47
Cal Check	0.101	11:47
Air Blank	0.000	11:48
Cal Check	0.101	11:49
Air Blank	0.000	11:49
Cal Check	0.101	11:50
Air Blank	0.000	11:50
Cal Check	0.100	11:51
Air Blank	0.000	11:52
Cal Check	0.100	11:52
Air Blank	0.000	11:53
Cal Check Stats		
Average	0.1003	
Std Dev	0.0004	
Rel Std Dev(%)	0.4815	

Manufacturer: Guth	Manufactured Date: 03/07/2022	Concentration: 0.10	
Lot No.: 22080	Expiration Date: 03/07/2024	Equipment Model: Guth 2100	
to 40 OAC 25-1-2.1. Approved breath-alco	oursuant to 40 OAC 30-1-3. The above NIST traceable re ohol measurement equipment and reference solutions. Eq nits for Breath Alcohol Testers (72 FR 34747), published t	Ulpment used appears on the current or supplemented	
Signature A /	Date:	Date: 06/27/2022	
Conforming Products List of Calibrating Un Signature Name Christopher Pape	06/27/2	022	



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis - Reference Solution

ALPHA INSTRUMENT

Intoxilyzer - Alcohol Analyzer

Model 8000

SN 80-002591

06/27/2022

11:11

Test	g/210L	Time
Air Blank	0.000	17:11
Cal Check	0.199	11:12
Air Blank	0.000	11:12
Cal Check	0.200	11:13
Air Blank	0.000	11:14
Cal Check	0.200	11:14
Air Blank	0.000	11:15
Cal Check	0.200	11:16
Air Blank	0.000	11:16
Cal Check	0.201	11:17
Air Blank	0.000	11:17
Cal Check	0.200	11:18
Air Blank	0.000	11:19
Cal Check	0.201	11:19
Air Blank	0.000	11:20
Cal Check	0.200	11:21
Air Blank	0.000	11:21
Cal Check	0.200	11:22
Air Blank	0.000	11:23
Cal Check	0.201	11:23
Air Blank	0.000	11:24
Cal Check Stats		
Average	0.2002	
Std Dev	0.0006	
Rel Std Dev(%)	0.3159	

NIST Traceable Breath-A	lcohol Equivalent (g/210L +/- 3%) Ratio of ethanol and	water free of test interfering substances.
Manufacturer: Guth	Manufactured Date: 02/07/2022	Concentration: 0.20
Lot No.: 22050	Expiration Date: 02/07/2024	Equipment Model: Guth 2100
prescribed by the State Director of Tests to 40 OAC 25-1-2.1. Approved breath-ald	ia Infrared Spectroscopy utilizing the Intoxilyzer 8000 and pursuant to 40 OAC 30-1-3. The above NIST traceable r cohol measurement equipment and reference solutions. E Units for Breath Alcohol Testers (72 FR 34747), published	eference solution meets the requirements for use pursuant
Signature	Date: 06/27/	
Name Christopher Pape	Permit No.: 230	020



Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 21120 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on April 16, 2021, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.0237% (w/vol) ethyl alcohol. The expiration date for this lot number is April 7, 2023 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.020 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN08031603 whose values are traceable to NIST.



RCV'D 8/26/21 CP

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 21070 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on March 2, 2021, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.0488% (w/vol) ethyl alcohol. The expiration date for this lot number is March 1, 2023 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.040 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN08211802 whose values are traceable to NIST.



5/9/2022

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 22080 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on March 9, 2022, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1214% (w/vol) ethyl alcohol. The expiration date for this lot number is March 7, 2024 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN03052002 whose values are traceable to NIST.



Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 22050 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on February 9, 2022, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.2437% (w/vol) ethyl alcohol. The expiration date for this lot number is February 7, 2024 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.200 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN06231703 whose values are traceable to NIST.

Location Returned From	Unit Serial #	Date: 17/17
Stigler PD	80-004423	7/7/22
Unit checked in to BOT New Unit Checked out Value Value	Old Cylinder #: 1402008 ster from Unit	inventory
Check in Procedure Connect unit to ethernet connection 172.019.0 Download any stored records Change unit location to Serial number Update Unit Date and time Previous time: 7/8/22 0940 New Time: 7/8/22 0949 Execute a new Forms load on the unit	Notes from the Field Routine Swap	Issues in the Field
Lab Evaluation Breath Hose Hold Vacuum?: DVM 3μm 13/14 Speaker (F5) Power LED Green Gas Regulator (F1.1) Pump (F1.1) Atmospheric Sensor Pressure 966 Correction Facto Sent to CMI (If needed) Date Sent —	RFI Detected? Red (F7) Receipt paper Flow Sensor Resistance Resistance	Is regulator close to Tank Sensor Value 956 Yes / No
Invoice Sent Date Returned Repairs needed:	Repairs Made:	
Inventory Process	Ready for Bench Check?	if unit sout to CAN
☐ Linked Cylinder to Unit ☐ Mock Subject ☐ Changed Unit Status ☐ Certificate o	t Test Completed GMI Workor Calibration and Operation CMI Certification	e if unit sent to CMI der ate of Calibration Evaluation Form Intoxilyzer Check-in sheet Version 2.0