

Dr. Jarrad Wagner, Ph.D. F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

ATTESTATION

I, the undersigned *Joshua Smith* (Director/Records Custodian), attest under oath pursuant to 12 O.S. §2902 (11)(a) this record was made, at or near the time of the occurrence of the matters set forth by or from information transmitted by a person having knowledge of those matters; was kept in the course of the regularly conducted business activity of the Board; and was made pursuant to the regularly conducted activity. The duplicate record provided is listed below:

1) The attached "Service Overview" is a true and accurate copy of maintenance performed on the Intoxilyzer 8000, serial number 80-004876, in compliance with the agency's Maintenance Policy and Procedures.


This packet consists of 20 pages.

Signed _____

Name of Position: Director/Records Custodian

Date of Attestation: 06/15/2022

Service Overview



3600 N. Martin Luther King Bldg #9 Oklahoma City, OK 73111 405-425-2460

BENCH CHECK DATE: 06/15/2022

SERIAL: 80-004876

DOCUMENTATION WITHIN SERVICE PACKET:

- BOT Certificate of Maintenance, Calibration and Operation
- Bench Check
- CMI Certificate of Calibration
- CMI Evaluation Form
- CMI Work Order
- BOT Dry Gas Certificate of Analysis
- Ilmo Specialty Gases Certificate of Analysis
- Mock Subject Test
- BOT Certificate of Analysis Reference Solution (4)
- Guth Laboratories Certificate of Analysis (4)



Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

CERTIFICATE of MAINTENANCE, CALIBRATION, and OPERATION

This certifies the Intoxilyzer 8000 referenced below, manufactured by CMI, Inc., a subsidiary of MPD, Inc. of Owensboro, Kentucky, is maintained by the Board of Tests for Alcohol and Drug Influence in accordance with 40 O.A.C. 30-1-3. The instrument reported below is hereby deemed to be in compliance with Board rules via linearity accuracy examination and operational/functional procedure evaluation also referred to as a "Bench Check". The Intoxilyzer 8000 appears on the National Highway Traffic Safety Administration's Conforming Products List of Evidential Breath Alcohol Measurement Devices (82 FR 50940 et seg) as an approved device. In accordance with 40 O.A.C. 25-1-2, this instrument is approved and certified for use in the State of Oklahoma as an automated evidential breath-alcohol measurement device by the Board of Tests pursuant to 47 O.S. § 759.

All analyses conducted by this instrument are considered precise, valid, and in accordance with 40 O.A.C. 30-1-3(b) Operating Procedure(s). A printed breath test operator affidavit containing the seal of the Board by this instrument confirms a calibration check successfully verified the instrument's calibration at the time of test. The printed breath test affidavit also serves as proof that safeguards in place to protect the breath samples, including but not limited to Ambient Air Interference, Radio Frequency Interference, Residual Mouth Alcohol, Deficient Samples, etc. were not triggered and the test is a true and accurate scientific measurement of the person's alcohol concentration.

- Maintenance was performed on the Intoxilyzer 8000, Serial No. 80-004876 on 06/15/2022 in accordance with 40 O.A.C. 30-1-3 and internal BOT maintenance policy and procedures.
- A dry gas canister appearing on the National Highway Traffic Safety Administration's Conforming Products List of Calibrating Units for Breath Alcohol Testers (77 FR 64588 et seq) and approved in accordance with 40 O.A.C. 25-1-3 and BOT maintenance policy was installed.
 - Manufacturer: ILMO Unit: 208 ppm / 105 L
 - Cylinder Lot No. Canister No.: 06022080A2-052
 - BOT analysis average test result: 0.081
- Cylinder expiration date: 04/05/2024

3.	In accordance with 40 O.A.C. 25-1-	I, the approved reference st	tandards and equipment we	re usea:
	Simulator Model			
	O. 45 0400	0.45.0400	0.45.0400	- 0

Simulator Model			
Guth 2100	Guth 2100	Guth 2100	Guth 2100
Serial #			
DR3685	DR3694	DR3594	DR3754
Concentration			
.02	.04	.10	.20
LOT#	7777		
21120	21070	21190	20510
Manufactured Date			
4/7/2021	3/1/2021	6/8/2021	12/3/2020
Expiration Date			
4/7/2023	3/1/2023	6/8/2023	12/3/2022
Solution Commission Date		- A - A - A - A - A - A - A - A - A - A	
06/08/2022	06/08/2022	06/08/2022	06/08/2022

1 200	06/08/2022	06/08/2022	06/08/2022
Completed by BC	OT Personnel: Christopher Pape	Breath-Alcohol Test	ing Program Administrator
Signature/Mame		Title	

Joshua Smith 9ignature Name

ertified

State Director Title

STATE OF OKLAHOMA - BOARD OF TESTS FOR ALCOHOL AND DRUG INFLUENCE CMI INTOXILYZER BENCH CHECK PROCEDURE REPORT

NCH CHECK DATE 06/15/2022	10:30:5		INSTRUMENT LOCATION 80-00487	6	
TRUMENT TYPE Intoxilyzer	INSTRUMENT MO			80-00	
ncoxilyzer	1 000		ESULTS	00 00	1070
121.55631					
Diagnostic PASS	5				
Mouth Alcohol	Test				
PASS					
RFI Test					
PASS		127		Shill all all	235
Abort Test	7	Du	plicate Breath	Sampling Re	sults
PASS		Value	Test 1	Test 2	Control Test
Improper Sample	e Test			0.010	
PASS		.02	0.020	0.019	0.020
Interferent Detec	t Test	.04	0.041	0.040	0.040
PASS		.10	0.102	0.101	0.102
No Sample Give	n Test	.20	0.204	0.202	0.199
Insufficient Samp	le Test		PA	SS	
Air Blank Contamina	ation Test				
	d Test				

CMI Form OKBE1 01/2019

PAPE,

CHRISTOPHER

0000023020

Certificate of Calibration

This is to certify the calibration of INONINZET® serial number 80-004876, manufactured by CMI Inc., a subsidiary of MPD, Inc. of Owensboro, Kentucky, was tested and found to conform to the National Highway Traffic Safety Administration (NHTSA) Standard for Devices to Measure Breath Alcohol (Federal Register, Vol.58, No.179, pp 48705-48710, Sept. 17, 1993) for accuracy and precision. Reference materials are traceable through the National Institute of Standards and Technology (NIST) to the International System of Units (SI).

Date June 3, 2022

Signed 3

Technician

INC.

316 East 9th Street Owensboro, KY 42303 USA

Part No. 650517 Rev.A

CMI INC. Service Evaluation Form

This form MUST be completed and enclosed with instrument to be serviced. Failure to complete and return this form may cause delays in service.

(Note: please	ship items in their	original shipping cont	ainer(s) or a similar protective box.)
. Contact informa	ation: Custome	r Number 731362	(contact Customer Service)
Name: Board	d of Tests	Phone: (405) 425 2460
Fax: (405) 425	2490	Email: christopher.	pape@bot.ok.gov
2. Bill to Address:	Board of Tests	Ship to Addr	ess: Board of Tests
	3600 N ML King, B	ldg 9	3600 N ML King, Bldg 9
	OKC, OK 73111		OKC, OK 73111
. Serial Number:	80-004876	Instrument M	odel: Intoxilyzer 8000
values are lov ***Hazardou I Authorize Rep	v. Replace Flow se s Material Warning! pairs Up To: \(\sum_All \)	nsor. Green LED is we - <u>DO NOT</u> return gas cy ☐\$250 ☑\$500 ☐\$7	nospheric sensor and Flow Sensor ak. Inders with instrument!*** Other \$
	7 - CENTRAL - MANAGEM - M	2017.01.11.11.11.11.11.11	
uthorized By:			Ship to:
Chris Pape		Breath Testing Admin	
lame (Please Prir	111)	itle	
1/2/1/	110		CMI, Inc
Wille 1		9/21/2021	Attn: Service Dept.
signature //		9/21/2021 Date	Attn: Service Dept. 316 East Ninth Street
Ca	II: 405-425-2460 or e	Date mail_	Attn: Service Dept.
 Control of the control of the control		Date mail_	Attn: Service Dept. 316 East Ninth Street
Ca No, please sen Note: An e	II: 405-425-2460 or e d estimate before re stimate will be faxed	Date mail epairs are made. d before performing any	Attn: Service Dept. 316 East Ninth Street



Phone: 866-835-0690 Fax: 270-685-6268

SERVICE WORK ORDER

406352

DATE: 9/27/2021

Bill To:

Oklahoma Board Of Tests Blvd Adcu Annex Bldg #9 3600 Martin Luther King Oklahoma City, OK 73136-3000 USA

Customer #: 731362 Phone: 405-425-2424 Fax: 405-425-2490 fax Ship To:

Oklahoma Board Of Tests Blvd Adcu Annex Bldg #9 3600 Martin Luther King Oklahoma City OK 73136-3000 USA

Phone: 405-425-2424

MODEL #: 002480OK

Serial Number: 80-004876

BILL CODE: Out of Warranty

A 4!4 . .

EXTRA PARTS RCVD:

KB(W/COVER), BATTERY B11-2525, HAS SHELF W/REG, BH

DESCRIPTION OF PROBLEM

VALUE LOWER THAN EXPECTED, CURRENTLY AT 11044, ATMOSPHERIC SENSOR AND FLOW SENSOR VALUES ARE LOW. REPLACE FLOW SENSOR. GREEN LED IS WEAK.

WORK PERFORMED:

Display was dim; LED was burnt out; Pressure sensor/erratic; Solenoid leaked; Contamination found in recirculation port; Contamination found in cell; Tygon tubing/Dirty; System Board battery was low; Fan Filter/Dirty; Contamination found on source; 2-D Barcode reader did not light up; Replaced all parts listed; Calibrated and Final Tested; 8133.15 Software;

PARTS USED . ..

Seq. No.	Part	Description	Quantity
70	650517	CERTIFICATE OF CALIBRATION	2.00 EA
80	441169	COVER DUST,5/8in x 1/2in	1.00 EA
90	470145	CAP,PLASTIC,.25IDx.50LG,RED	3.00 EA
100	470154	CAP,PLSTC,.406IDx.438562,RED	1.00 EA
110	470218	SEAL,TAMPER EVIDENT,SERVICE	1.00 EA
111	320643	HDR,14PIN,2ROW,.100CC,4-SIDED	1.00 EA
112	140112	DISPLAY,2x20 VFD 9MM	1.00 EA
115	021355	LED ASSY,POWER,I8000	1.00 EA
120	210120	SENSOR,SMD,PRESSURE	1.00 EA
140	021357	SOLENOID ASSY,18000	1.00 EA
145	441014	FTG,RECIRC. OUT	1.00 EA
146	441109	FTG,QUICK DISC,M,1/8FL 1/8NPT	1.00 EA
147	441044	FTG,RT ANGLE,1/8BARD X 1032	1.00 EA
150	690040	BATTERY,3VOLT,LITHIUM,CR2032	1.00 EA
160	470202	FILTER,FOAM,1.85x2.20	1.00 EA
161	530030	TUBING,TYGON,.375IDx.562OD	2.00 IN
162	021383FLS	CELL SUB-ASSY,SVC	1.00 EA
165	020982	SOURCE ASSY	1.00 EA
170	210119	SENSOR,SMD,PRESSURE	1.00 EA
180	140119	READER,2D BARCODE,SMART FOCUS	1.00 EA

LABOR / TESTING

Misc Code	Description	Hours
LABR	Service Repair Labor	
LCAL	Service - Calibration Adjust	0.50
LFT	Service - Cal / Final Test	1.00

Technician Name: Zach Estes

SRV4063520001

Page: 1 of 3



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis - Dry Gas

ALPHA INSTRUMENT Intoxilyzer - Alcohol Analyzer Model 8000

04/25/2022

SN 80-002591 13:08

Test	g/210L	Time
Air Blank	0.000	13:08
Cal Check	0.080	13:08
Air Blank	0.000	13:09
Cal Check	0.080	13:09
Air Blank	0.000	13:10
Cal Check	0.080	13:10
Air Blank	0.000	13:11
Cal Check	0.082	13:11
Air Blank	0.000	13:11
Cal Check	0.081	13:12
Air Blank	0.000	13:12
Cal Check	0.081	13:13
Air Blank	0.000	13:13
Cal Check	0.081	13:13
Air Blank	0.000	13:14
Cal Check	0.081	13:14
Air Blank	0.000	13:15
Cal Check	0.080	13:15
Air Blank	0.000	13:16
Cal Check	0.081	13:16
Air Blank	0.000	13:16
Cal Check Stats		
Average	0.0807	
Std Dev	0.0006	
Rel Std Dev(%)	0.8363	

	h-Alcohol Equivalent (g/210L) Ratio of Nitrogen and 0.08 g/210L = 208 PPM ÷ 2605*	Ethanol Gas (208 ppm/105 L)						
Manufacturer: ILMO	Manufacturer Cert. No.: 14689	Cylinder Size: 105L						
Lot No.: 06022080A2	Cylinder No.: 52	Expiration Date: 04/05/2024						
The above 105 L cylinder has been analyzed via Infrared Spectroscopy utilizing the Intoxilyzer 8000 and was performed in accordance with the procedures prescribed by the State Director of Tests pursuant to 40 OAC 30-1-3. The above NIST traceable reference standard meets the requirements for use as a calibrating unit in accordance with 47 O.S. 759, 40 OAC 25-1-3, and appears on the current conforming products list of calibrating units for breath alcohol testers published by the U.S. Department of Transportation in the Federal Register [77 FR 64588 et seq]. * See Dubowski, K.M. and Essary N., <i>Journal of Anal. Toxicol.</i> (1996) 20(6): pp. 484-491. [g/210L = PPM +2605].								
Clausetine		The state of the s						
Name Christopher Pape	Date: 04/25/	/2022						



7 Eastgate Dr. • P.O. Box 790 • Jacksonville, IL 62651-0790 217-245-2183 • Fax: 217-243-7634 • www.ilmoproducts.com

Certificate of Analysis

Certificate ID:

14689

Part #:

BAC105L080T

Cylinder Size:

105L

Lot Number:

06022080A2

Expiration:

4/5/2024

4/11/2022

0.080 BAC (For the calibration of instruments used to determine breath alcohol concentration)

Contents:

105 Liters @ 1000 psig 70°F (21°C)

Analytical

Reported

Accuracy

Analytical

Component:

Concentration:

(U, k=2):

Method:

Ethanol

208 ppm

+/-0.002 BAC(G/210L) NDIR [5.2 ppm]

Nitrogen

Balance

CMI Inc.

Distributed by:

316 East Ninth Street Owensboro, KY 42303 Phone 866-835-0690 www.alcoholtest.com

*Traceable to:

Certified Reference Material - 251.8 µmol/mol

Ethanol in Nitrogen - Serial No. CC198518 Lot No. 111318E1

Store in dry area, away from sources of heat, ignition and direct sunlight. Do not allow storage area to exceed 52 °C (125 °F).

Specialty Gas Lab Tech

03-28-2022

Issuance Date



The calibration results/within this certificate were obtained using equipment and standards capable of producing analytical results traceable to NIST, and apply only to the items contained on this certificate. ILMO Products Company makes no warranty or representation as to the suitability of the use of any information provided for any particular purpose. The information use is at the sole discretion and risk of the user. Liability shall be limited to established replacement cost of this material or service.

OFFICER'S AFFIDAVIT AND NOTICE OF 10/20 REVOCATI

ON/DISQUALIFICATION	
"아이 HANDY YEAR HOUNTY TO IN IT TO THE FOLK HERE HE	Section I

ARREST DATE	A	ST TIME	TOT WID	COUNTY	PARTY C	CITY#	444	TION#		Section
06/15/2022		00		55		51111	- 1			
ARREST LOCATION BOT					CITY				UNTY OKLAHO	OMΔ
UBJECT NAME						ATE OF BIRTH		EX	HEIGHT	the state of the s
BREATHTEST,	DELTA					01/01/1		M	510	
DDRESS 3600 MLK			OKT	AHOMA (VTT		OI		73	
RIVER LICENSE #			OLL	EXPIRATION	DATE	STATE	CLASS		STRICTIONS	ENDORSEMENT
B083463042				10/31	1/202		D		7,000	EDEX
EHICLE MAKE MOCK		MOCK				MOCK		OK	N CMV	N N
the above date, time, and loca	ation, the above named perso			unds to believe the	e person wa		al physical con			
ghway, street, turnpike, public p luence of atcohol and/or other i escribe driving behavior or circt MOCK SUBJECT	ntoxicating substances as pro imstances): TEST		road, street, alley or t	ane which provide	s access (c	one or more single	or multi-tamily	dwelling with	nin this state v	vnile under the
scribe person's condition (odor		TNOWALLE	D NEW CVI	TNDED						
POST REPAIR M	AINTENANCE,	INSTALLE	D NEW CYL	INDER						
		THE PERSON V	VAS INFORMED OF	THE IMPLIED CO	NSENT TE	ST REQUEST				
EATH-ALCOHOL ANALYSIS	RECORD, REPORT OF THE		IN SECTION 1 ABOV	E, and LOG OF T	EST AND N	MAINTENANCE RE	CORD.			Section
I - 8000			-004876	BOARD	OF '	TESTS				
ANDARD			TARGET VALUE		PERMIT	NUMBER				
Dry			0.080		230	20				
NUFACTURER LLMO			06022	08022	O.6	15/202	2	O4/	05/202	24
SERVATION START TEST	DATE		_					04/0	33/202	4
.115 06	/15/2022					ation/observ and continui				
SULT TABLE			No. of the control of	the time of to	ACCO. 12.7 P. 100		ng		THE S	TATE
Test	g/210L	Time	101	ine time of te	est admi	nistration.		1	** * *	10
			I administe	arad said I	Pronth	Alcohol T	oct in	12/	***	1 * * * /S
Air Blank Subject Test	0.00	11:48 11:49	1 4 4 4 4 5 1 W. W. W. W. W. W. W.				SWITH BAR	188	ON	A SEZ E
Breath Volume	2.023 LITERS	11:49				cable opera		E-1 40		** A
Air Blank	0.00	11:49	proce	edure of th	e Boar	d of Tests		1 4	**	*** 5
V ait		11:51	1/1/-	1	10			12	*	1.5
Air Blank	0.00	11:51	La sala	WING	PC			-	1000	
Subject Test Breath Volume	0.00 2.003 LITERS	11:52	SIGNATURE OF	OPERATOR	-				190	
Air Blank	0.00	11:52	and the second							
Cal Check	0.07	11:53	NAME	D T GMODILI	10			III T	Test Re	esult
Air Blank	0.00	11:53	PAPE, CH	RISTOPHE	PERMIT	,			BrAC g	210L
			0000000	484		002302	0			
			AGENCY	h det debet	1 3 7/3				0.0	00
	the state of the s		BOARD OF	F TESTS				_		
FECTIVE THIRTY (30) DAY: so or more, if you refused or erating a CMV or CDL holde a result of this arrest by filing u may be eligible to partici ocation as a result of this ruest IDAP within thirty (30 RVING OF NOTICE: I countability Program to	failed the state's test(s). In r) OR a test result of .04 or a petition in the District Co- pate in the Impaired Driva arrest from appearing on) days of the date listed by personally hand deliv	addition, your cor more while opera surt of the County or Accountability your driving rec- celow. Mail the re- ered a copy of	nmercial driving privi ting a commercial m in which you were ar Program (IDAP) ad ord. Participation in equest for IDAP to t the above contai	lege in this state otor vehicle. You rested. Your pet Iministered by the IDAP may red the Department	will be dis a may appo- ition must he Departa ace the an of Public	qualified if you ref eal any Departme be filed within thir ment of Public S nount of fees you Safety at P.O. Bo	used or failed ntal action ag ty (30) days o afety (DPS). u will be requ ox 11415, Okl	the state's ainst your d f the date lis Completion aired to pay ahoma City	test(s) river license sted below. n of IDAP ma to the State n, OK 73136.	. You must
te served:	Sign	ature of officer						Badge	#:	
accordance with Title 12	2 O.S. Section 426, "I	state under pe	nalty of perjury	under the law	s of Okla	shoma that the	foregoing	s true and	d correct."	Section
te:	Place (location wh	en signed):	100 O C 34 N	1511	Signa	ature of arresti	ng officer:_			
APE, CHRI	STOPHER			BOARI	OF	TESTS				
OGE	AGENCY ADDRESS INTER-AGEI	icv		DOARCE	, OI	ZIP CODE 0000	PHON	5425	2460	
000000484	INIEK-AGEI	NC I				10000	0 40	3425	2400	Section
HER WITNESSES:			TITLE	ADDRESS			In	HONE		Section
NAME			TITALE.	ADDITE 33				TOTAL		
NAME			TITLE	ADDRESS	-		P	HONE		
2			77.7	2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2						

OFFICER'S AFFIDAVIT AND NOTICE OF REVOCATION/DISQUALIFICATION

11/20

Section I COUNTY # ARREST DATE CITATION # 06/15/2022 1100 55 ARREST LOCATION BOT OKLAHOMA SUBJECT NAME DATE OF BIRTH HEIGHT 01/01/1985 BREATHTEST, DELTA M 510 200 ZIP CODE 3600 MLK OKLAHOMA CITY 73111 OK CLASS ENDORSEMENTS DRIVER LICENSE # EXPIRATION DATE RESTRICTIONS B083463042 10/31/2022 OK D VEHICLE MAKE MODEL PLACARD MOCK MOCK MOCK N OK N On the above date, time, and location, the above named person was arrested, and I had reasonable grounds to believe the person was driving or in actual physical control of a motor vehicle upon a public road, highway, street, turnpike, public parking lot, or other public place, or upon a private road, street, alley or lane which provides access to one or more single or multi-family dwelling within this state while under the influence of alcohol and/or other intoxicating substances as prohibited by law. (Describe driving behavior or circumstances): MOCK SUBJECT TEST Describe person's condition (odor, actions, etc.): POST REPAIR MAINTENANCE, INSTALLED NEW CYLINDER THE PERSON WAS INFORMED OF THE IMPLIED CONSENT TEST REQUEST BREATH-ALCOHOL ANALYSIS RECORD, REPORT OF THE PERSON NAMED IN SECTION 1 ABOVE, and LOG OF TEST AND MAINTENANCE RECORD. Section 2 SERIAL NUMBER SPECIALIST I-8000 80-004876 BOARD OF TESTS PERMIT NUMBER STANDARD ARGET VALUE Dry MANUFACTURER 0.080 23020 DATE CYLINDER EXPIRES DATE CYLINDER INSTALLED OT # ILMO 06022080A2 06/15/2022 04/05/2024 DRSERVATION START TEST DATE A 15-minute (or greater) deprivation/observation 1115 06/15/2022 period was observed prior to and continuing RESULT TABLE to the time of test administration. Time Test g/210L I administered said Breath-Alcohol Test in Air Blank 0.00 11:48 Subject Test 0.00 11:49 compliance with the applicable operating Breath Volume 2.023 LITERS procedure of the Board of Tests. Air Blank 0.00 11:49 11:51 Wait Air Blank 0.00 11:51 Subject Test 0.00 11:52 Breath Volume 2.003 LITERS 0.00 Air Blank 11:52 VAME Cal Check 0.07 11:53 Test Result PAPE, CHRISTOPHER Air Blank 0.00 11:53 BrAC g/210L 0000000484 0000023020 0.00 AGENCY BOARD OF TESTS EFFECTIVE THIRTY (30) DAYS FROM DATE SERVED. Notice of Revocation is given that your driver's license is revoked or denied for a period of one hundred eighty (180) Section 3 days or more, if you refused or failed the state's test(s). In addition, your commercial driving privilege in this state will be disqualified if you refused or failed the state's test(s) (operating a CMV or CDL holder) OR a test result of .04 or more while operating a commercial motor vehicle. You may appeal any Departmental action against your driver license as a result of this arrest by filing a petition in the District Court of the County in which you were arrested. Your petition must be filed within thirty (30) days of the date listed below. You may be eligible to participate in the Impaired Driver Accountability Program (IDAP) administered by the Department of Public Safety (DPS). Completion of IDAP may prevent a revocation as a result of this arrest from appearing on your driving record. Participation in IDAP may reduce the amount of fees you will be required to pay to the State. request IDAP within thirty (30) days of the date listed below. Mail the request for IDAP to the Department of Public Safety at P.O. Box 11415, Oklahoma City, OK 73136. SERVING OF NOTICE: I personally hand delivered a copy of the above containing the Notice of Revocation and the Notice of Department Impaired Driving Accountability Program to the person named in Section 1 above. Signature of officer: Section 4 In accordance with Title 12 O.S. Section 426, "I state under penalty of perjury under the laws of Oklahoma that the foregoing is true and correct." Place (location when signed): Signature of arresting officer: Date: AGENCY NAME PAPE CHRISTOPHER BOARD OF TESTS AGENCY ADDRESS 0000000484 00000 4054252460 INTER-AGENCY



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis - Reference Solution

ALPHA INSTRUMENT 2 Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-004923 03/08/2022 11:33

Test	g/210L	Time
Air Blank	0.000	11:34
Cal Check	0.022	11:34
Air Blank	0.000	11:35
Cal Check	0.021	11:36
Air Blank	0.000	11:36
Cal Check	0.021	11:37
Air Blank	0.000	11:37
Cal Check	0.021	11:38
Air Blank	0.000	11:38
Cal Check	0.021	11:39
Air Blank	0.000	11:40
Cal Check	0.021	11:40
Air Blank	0.000	11:41
Cal Check	0.020	11:42
Air Blank	0.000	11:42
Cal Check	0.021	11:43
Air Blank	0.000	11:43
Cal Check	0.021	11:44
Air Blank	0.000	11:45
Cal Check	0.020	11:45
Air Blank	0.000	11:46
Cal Check Stats		
Average	0.0209	
Std Dev	0.0005	
Rel Std Dev(%)	2.7158	

Manufacturer: Guth	Manufactured Date: 04/07/2021	Concentration: 0.02	
Lot No.: 21120	Expiration Date: 04/07/2023	Equipment Model: Guth 2100	
prescribed by the State Director of Test to 40 OAC 25-1-2.1. Approved breath-a Conforming Products List of Calibrating	via Infrared Spectroscopy utilizing the Intoxilyzer 8000 and w s pursuant to 40 OAC 30-1-3. The above NIST traceable ref lcohol measurement equipment and reference solutions. Equ Units for Breath Alcohol Testers (72 FR 34747), published b	erence solution meets the requirements for use pursuant sipment used appears on the current or supplemented	
Signature AMIL OF	Date: 03/08/2	Date: 03/08/2022	
Name Christopher Pape	Permit No.: 2303	Permit No.: 23020	



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis - Reference Solution

ALPHA INSTRUMENT 2 Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-004923 03/08/2022 11:53

Test	g/210L	Time
Air Blank	0.000	11:53
Cal Check	0.041	11:54
Air Blank	0.000	11:54
Cal Check	0.041	11:55
Air Blank	0.000	11:56
Cal Check	0.042	11:56
Air Blank	0.000	11:57
Cal Check	0.042	11:58
Air Blank	0.000	11:58
Cal Check	0.042	11:59
Air Blank	0.000	11:59
Cal Check	0.041	12:00
Air Blank	0.000	12:01
Cal Check	0.041	12:01
Air Blank	0.000	12:02
Cal Check	0.041	12:02
Air Blank	0.000	12:03
Cal Check	0.041	12:04
Air Blank	0.000	12:04
Cal Check	0.042	12:05
Air Blank	0.000	12:05
Cal Check Stats		
Average .	0.0414	
Std Dev	0.0005	
Rel Std Dev(%)	1.2473	

Manufacturer: Guth	Manufactured Date: 03/01/2021	Concentration: 0.04	
Lot No.: 21070	Expiration Date: 03/01/2023	Equipment Model: Guth 2100	
to 40 OAC 25-1-2.1. Approved breath-a Conforming Products List of Calibrating	s pursuant to 40 OAC 30-1-3. The above NIST traceable ref lcohol measurement equipment and reference solutions. Eq Units for Breath Alcohol Testers (72 FR 34747), published b	ulpment used appears on the current or supplemented	
Signature	Date: 03/08/2		
Name Christopher Pape	Permit No.: 0000	Permit No.: 23020	



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis - Reference Solution

ALPHA INSTRUMENT 2 Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-004923 03/08/2022 12:11

(5)		
Test	g/210L	Time
Air Blank	0.000	12:11
Cal Check	0.101	12:12
Air Blank	0.000	12:13
Cal Check	0.101	12:13
Air Blank	0.000	12:14
Cal Check	0.101	12:15
Air Blank	0.000	12:15
Cal Check	0.102	12:16
Air Blank	0.000	12:16
Cal Check	0.101	12:17
Air Blank	0.000	12:18
Cal Check	0.101	12:18
Air Blank	0.000	12:19
Cal Check	0.101	12:20
Air Blank	0.000	12:20
Cal Check	0.101	12:21
Air Blank	0.000	12:21
Cal Check	0.101	12:22
Air Blank	0.000	12:23
Cal Check	0.101	12:23
Air Blank	0.000	12:24
Cal Check Stats		
Average	0.1011	
Std Dev	0.0003	
Rel Std Dev(%)	0.3127	

Manufacturer: Guth	Manufactured Date: 06/08/2021	Concentration: 0.10	
Lot No.: 21190	Expiration Date: 06/08/2023	Equipment Model: Guth 2100	
prescribed by the State Director of Tests to 40 OAC 25-1-2.1. Approved breath-alc Conforming Products List of Calibrating U	a Infrared Spectroscopy utilizing the Intoxilyzer 8000 and w pursuant to 40 OAC 30-1-3. The above NIST traceable rel ohol measurement equipment and reference solutions. Eq Inits for Breath Alcohol Testers (72 FR 34747), published b	ference solution meets the requirements for use pursuan uipment used appears on the current or supplemented by the National Highway Traffic Safety Administration.	
Signature Chapter A		Date: 03/08/2022	
Name Christopher Pape	Permit No.: 220	Permit No.: 23020	



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

Certificate of Analysis - Reference Solution

ALPHA INSTRUMENT 2 Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-004923 03/08/2022 12:28

Test	g/210L	Time
Air Blank	0.000	12:29
Cal Check	0.201	12:29
Air Blank	0.000	12:30
Cal Check	0.200	12:31
Air Blank	0.000	12:31
Cal Check	0.201	12:32
Air Blank	0.000	12:33
Cal Check	0.201	12:33
Air Blank	0.000	12:34
Cal Check	0.201	12:35
Air Blank	0.000	12:35
Cal Check	0.201	12:36
Air Blank	0.000	12:36
Cal Check	0.201	12:37
Air Blank	0.000 -	12:38
Cal Check	0.200	12:38
Air Blank	0.000	12:39
Cal Check	0.201	12:40
Air Blank	0.000	12:40
Cal Check	0.200	12:41
Air Blank	0.000	12:41
Cal Check Stats		
Average	0.2007	
Std Dev	0.0004	
Rel Std Dev(%)	0.2406	

Manufacturer: Guth	Manufactured Date: 12/03/2020	Concentration: 0.20
Lot No.: 20510	Expiration Date: 12/03/2022	Equipment Model: Guth 2100
to 40 OAC 25-1-2.1. Approved breath-a Conforming Products List of Calibrating	ts pursuant to 40 OAC 30-1-3. The above NIST traceable r alcohol measurement equipment and reference solutions. E Junits for Breath Alcohol Testers (72 FR 34747), published	guipment used appears on the current or supplemented
Signature Assignment of the Si	Date: 03/08/	2022
Name Christopher Pape		



RCVD 8/26/21

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 21120 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on April 16, 2021, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.0237% (w/vol) ethyl alcohol. The expiration date for this lot number is April 7, 2023 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.020 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN08031603 whose values are traceable to NIST.



RCV'D 8/26/21 CP

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 21070 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on March 2, 2021, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.0488% (w/vol) ethyl alcohol. The expiration date for this lot number is March 1, 2023 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.040 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN08211802 whose values are traceable to NIST.

RCV D 8/26/21

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 21190 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on June 9, 2021, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1206% (w/vol) ethyl alcohol. The expiration date for this lot number is June 8, 2023 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN03052002 whose values are traceable to NIST.



RCV'D \$/26/21 CP

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 20510 of
Alcohol Reference Solution for Simulator were analyzed by
gas chromatography on December 7, 2020, using a Perkin Elmer Gas
Chromatograph Autosystem XL S/N: 610N9030209, and found to contain
0.2419% (w/vol) ethyl alcohol. The expiration date for this lot
number is December 3, 2022 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.200 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN06231703 whose values are traceable to NIST.

kilyzer Check-in sheet Version 2.0	luto	BOT Internal Use Only	10-1	20/20
	☐ CMI Service Evaluation Form			Bench Check Report Complete
	CMI Certificate of Calibration	n and Operation	Certificate of Calibratio	Changed Unit Status
	CMI Workorder	pleted	Mock Subject Test Com	☐ Linked Cylinder to Unit
	To complete if unit sent to CMI	n Report	☐ Bench Check Technicia	Inventory Process
	Bench Check?	Ready for		
,				
				·
			1, VM Pa-	10000 TWO
	<u>ade:</u>	M szisgeA		Repairs needed:
resvid s	on: Neplece Flow	1 30100 (2)	714 EE,	Date Returned 6/15/
7 7 7 3	ow. Replace Flow	1 0 100	57	Invoice Sent
	ouss Than expected	1 JOIN W	15 77	Date Sent
그를 즐기 때문에 가를 받는다 되었다. 기계, 그 기기 기기 기기를 보고 있다.	LOCOUND IN H 7.		Reason Sent	Sent to CMI (If needed)
		. 1		
ON / (Sal)	48P = 136	onstsiseA	Correction Factor	Pressure 886
Tank Sensor Value?	Sor Tank Sensor	Flow Sens	/// (Atmospheric Sensor
ls regulator close to		·		·
	aper	Receipt p	Pump (F10)	Gas Regulator (F11)
	Fan (F9)	Red (F7)	Power LED Green (F6)	Speaker (F5)
		7	20 A C	
e e		η		<u>Peripheral Tests</u>
	7FI Detected?		Breath Hose Temp	LC) url6
1	iyo asau daya	. , -	Ch ame Teso Hatrens	67 1/11 winds
	Snap Test ok?	,	Chamber Temp	hhall mule
			•	DAM
	ان: ا	nber Hold Vacuun	ıezr cusı	Breath Hose Hold Vacuum?:
			/	
				Lab Evaluation
			าเบท อนา นด	Execute a new Forms load o
			1011	_ /
			LTH 14	5/4/8 :9miT w9N
			8/hl 181	Previous time:
			•	emit bns etsel tinU etsbqU
				السر
	ble remotelγ.	l l	ial number	To Ser Inotation to Ser
eldenu bne noitoenno	out of tolerance. Unit on one-way o	Cal check	sbr	Download any stored reco
·	tine Swap	nov [] -	onnection 172.019.004. 135	2 Jounect unit to etnember c
			351 100 010 051 4544444	
	<u>bləi4 ma</u>	Notes fro		Check in Procedure
	aluation) Old cylinder archived	sv∃ der Lab Ev	Memove old cylinder from	V New Unit Checked out
	uger disposed in inventory	uit 🗖 Old Cylli	Umlink Gas Canister from U	V Unit checked in to BOT
	nder #: 05620080A2-019	OI9 CVI!		
	DEEDOODA DAO	,	,	Inventory Process
17/81/8		9८8700-08		McCurtain Co.
Date:		4 leirae 1 jin L		Location Returned From
.0450	İ	4 10:20 4:01		AMAN LONNING MOHONO