

STATE OF OKLAHOMA BOARD OF TESTS FOR ALCOHOL AND DRUG INFLUENCE

Dr. Jarrad Wagner, Ph.D. F-ABFT Chairman Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

ATTESTATION

I, the undersigned *Joshua Smith* (Director/Records Custodian), attest under oath pursuant to 12 O.S. §2902 (11)(a) this record was made, at or near the time of the occurrence of the matters set forth by or from information transmitted by a person having knowledge of those matters; was kept in the course of the regularly conducted business activity of the Board; and was made pursuant to the regularly conducted activity. The duplicate record provided is listed below:

 The attached "Service Overview" is a true and accurate copy of maintenance performed on the Intoxilyzer 8000, serial number , in compliance with the agency's Maintenance Policy and Procedures.

This packet consists of pages. Signed

Name of Position: Director/Records Custodian

Date of Attestation:

Service Overview

Board of Tests for Alcohol and Drug Influence Enhancing public safety by emaining the accuracy and acientific reliability of blood and breath tasting

> 3600 N. Martin Luther King Bldg #9 Oklahoma City, OK 73111 405-425-2460

BENCH CHECK DATE: 03-26-2021

SERIAL: 80-005042

DOCUMENTATION WITHIN SERVICE PACKET:

- BOT Certificate of Calibration and Operation
- CMI Certificate of Calibration 03-18-2021
- Bench Check
- Bench Check Technician's Report
- BOT Dry Gas Certificate of Analysis
- <u>CMI Service Evaluation Form</u>
- CMI Work Order
- Mock Subject Test
- BOT Certificates of Analysis on Certified Solutions (4)



STATE OF OKLAHOMA BOARD OF TESTS FOR ALCOHOL AND DRUG INFLUENCE

Administrative Office at: Board of Tests for Alcohol and Drug Influence Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307

Certificate of Calibration and Operation

This is to certify that the calibration of INTOXILYZER 8000, serial number <u>80-005042</u>, manufactured by CMI Inc., a subsidiary of MPD, Inc. of Owensboro, Kentucky, was tested and found to be in compliance with National Highway Traffic Safety Administration Standard for Devices to Measure Breath Alcohol (F.R., No. 179 48705-48710 Sept. 17, 1993) and with International Recommendation OIML R 126 Edition 1998 (E). Calibration solutions are traceable to NIST standard material 1828.

In addition the operational procedures and functions of this Intoxilyzer have also been tested and found to be in compliance with the Oklahoma Board of Tests for Alcohol and Drug Influence - Rules and Regulations OAC Title 40: 25-1-2 and 30-1-3 and is hereby certified for use as an evidential breath testing device within this state.

CERTIFICATION TYPE	DATE
Manufacturer Calibration	02-17-10
Board of Tests Bench Check	03-23-10
Board of Tests Bench Check	12-12-11
Board of Tests Bench Check	11-20-13
Board of Tests Bench Check	09-11-15
Board of Tests Bench Check	06-21-17
Board of Tests Bench Check	05-06-19
Manufacturer Calibration	03-18-21
Board of Tests Bench Check	03-26-21



Presentation of this form certifies that the Intoxilyzer, SN: <u>80-005042</u>, functioned properly at the time of the breath test and hereby certifies the breath test result as valid.

Date:

Re: Request for Log of Test on Intoxilyzer 8000

Name:

Dear Sir:

The Log of Test information from a breath test conducted on a Intoxilyzer 8000 is printed in Section 2 of the Affidavit and is no longer retained as a separate form as with the Intoxilyzer 5000. The 8000 uses a nitrogen/ethanol dry gas for the external control test. The dry gas cylinder has the capacity to perform approximately 325 tests and does not require maintenance until its expiration date or can no longer provide the required pressure to perform a control test. The cylinder has an expiration date of two years from date filled at the manufacturer and is programmed in the 8000 upon installation, both dates are printed in Section 2 of the Affidavit. The Intoxilyzer is programmed to disable itself when the date of expiration is reached or the pressure monitor indicates 25 lbs psi.

Due to the reliability and stability of dry gas, maintenance on the 8000 is not required every 30 days or 25 tests, therefore maintenance is no longer performed by officers in the field. The Board of Tests for Alcohol and Drug Influence certifies each Intoxilyzer 8000 for field use by performing a bench check. The bench check is designed to check the operation and procedures programmed in the 8000 to ensure it is functioning properly for use in Oklahoma. A bench check is performed every time a new dry gas cylinder is installed or the instrument has returned from the manufacturer for repair. Once the instrument passes the bench check, all that is required to set it up for field use is to plug it in and turn it on and does not require a permit to set up. The Board of Tests has the ability to monitor each instrument once it is set up in the field.

Each Intoxilyzer has a database of valid permit numbers that it checks prior to allowing a breath test to be performed. The database is updated on a regular basis, if a permit number is invalid that officer will not be able to conduct a breath test on the 8000. If the 8000 prints an affidavit the breath test operator listed in Section 2 had a valid permit at the time of the test and the breath test is valid.

If you have any further questions please feel free to call me.

Sincerely,

poher Blackelun

J. Robert Blakeburn (OBA11059) State Director of Tests for Alcohol and Drug Influence

Certificate of Calibration

This is to certify the calibration of Intoxilyzer ® serial number 80-005042, manufactured by CMI Inc., a subsidiary of MPD, Inc. of Owensboro, Kentucky, was tested and found to conform to the National Highway Traffic Safety Administration (NHTSA) Standard for Devices to Measure Breath Alcohol (Federal Register, Vol.58, No.179, pp 48705-48710, Sept. 17, 1993) for accuracy and precision. Reference materials are traceable through the National Institute of Standards and Technology (NIST) to the International System of Units (SI).

Date Mar18, 2021 Signed Alice Jo Technician lor




316 East 9th Street Owensboro, KY 42303 USA

Part No. 650517 Rev.A

6/16

STATE OF OKLAHOMA - BOARD OF TESTS FOR ALCOHOL AND DRUG INFLUENCE CMI INTOXILYZER BENCH CHECK PROCEDURE REPORT

CH CHECK DATE 3/26/2021	10:42			0 - 00504	2	
RUMENT TYPE ntoxilyzer					INSTRUMENT SE	
		TES	T RESI	JLTS		
Diagnost						
PASS						
Mouth Alcoh						
PASS						
RFI Tes						
PASS			Dupli	cate Breath	Sampling Re	sults
Abort Te						Control
PASS		Val	lue	Test 1	Test 2	Test
Improper Sam				100000		
PASS		. 0)2	0.019	0.018	0.019
Interferent Det		. 0)4	0.039	0.038	0.040
PASS		.1	.0	0.100	0.098	0.099
No Sample Giv PASS		.2	20	0.198	0.195	0.197
Insufficient San PASS	nple Test			PA	SS	
Air Blank Contam PASS						
Range Exceed						
		a las constants de la la		CHNICIAN		
swear under penalt he Intoxliyzer Benc	y of perjury the Checking P	hat in regards t	o the a	bove listed	instrument, I	have complie
URE	in checking i	rocedure appre	over 0	y the state c	DÁTE	
Mita	N.	230 11				2/01/0
V Myllaft	1	Lende	/			120/2
					PERMIT #	



3600 N. Martin Luther King Bldg #9 Oklahoma City, OK 73111 405-425-2460

20160

03/18/2020

03/18/2022

3/22/2021

BENCH CHECK TECHNICIAN'S REPOR

Date: 03/26/2021

Start Time: <u>10:40</u>

End Time: 12:45

INSTRUMENT

20070

02/13/2020

02/13/2022

3/22/2021

Intoxilyzer Model		Serial #				
8000		80-005042				
GAS CANISTER						
LOT #		EXP Date				
00521080A2-007		02/05/2	023			
REFERENCE	Simulate	or Model				
Guth 2100	Guth 2100	Guth 2100	Guth 2100			
	Se	erial #	the second s			
DR2834	DR3700	DR2455 DR34				
	Conc	entration				
02	.04	.10	.20			

LOT #

Manufactured Date

Expiration Date

Solution Commission Date

20190

04/06/2020

04/06/2022

3/22/2021

REASON FOR BENCH CHECK

☑ POST REPAIR ☑ CYLINDER REPLACEMENT

TROUBLESHOOTING ROUTINE MAINTENANCE

20060

02/10/2020

02/10/2022

3/22/2021

COMMENTS

Pressure sensor and Source replaced by CMI.

Technician Signature / Permit #

Christopher Pape #23020 Printed Name

22/2021	0.080					8/16	
f Analysis	Labelled target value (g/210L)		A CONTRACTOR OF A CONTRACTOR A CONTR				
07	0.081						
er#	Average test result				DI DUD INCTDUMENT	-	
HSTopher Pal	re cliquite D'al	23020			ALPHA INSTRUMENT Intoxilyzer - Alc		
chnician Name and Signat		the second second		Contraction of the local division of the loc	Model 8000		SN 80-00492
ordance with BOT Boll	cy and Procedure Statement BF	BIR BU	ecialty g	ases	02/22/2021		11:1
nced dry gas cylinder is	suitable for use as an external re	ference with the datasily zerrad	ove	onville, IL 62651-0790			
		217-245-2192 . Eav	2.60x /90 · Jacks	onville, IL 62651-0790	Test	0/2101	Tim
			-243-7634 .	www.ilmoproducts.com	Test	g/210L	Tim
					Air Blank	0.000	11:08
			2000 C		Cal Check	0.082	11:09
		Certi	ficate of	Analysis	Air Blank	0.000	11:09
				11101/515	Cal Check	0.081	11:1
	/	·	10		Air Blank	0.000	11:1
Certificate	ID: 13601		NCILL		Cal Check	0.081	11:1
- 영거야 않으면 생			KUV		Air Blank	0.000	11:1
Part #:	BAC105L08	ЗОТ	17. 1	21	Cal Check	0.081	11:1
Cylinder Si	ze: 105L		- /19/2	OOL .	Air Blank	0.000	11:1
		- /	1 1 10		Cal Check	0.081	11:1
I - A DI - I	00504000						
Lot Numbe			RCV12/19/2		Air Blank	0.000	11:1
Lot Numbe Expiration:	er: 005210804 2/5/2023		10		Air Blank Cal Check	0.000 0.081	11:1 11:1
			CP.		Air Blank Cal Check Air Blank	0.000 0.081 0.000	11:1 11:1 11:1
			CP.		Air Blank Cal Check Air Blank Cal Check	0.000 0.081 0.000 0.081	11:1 11:1 11:1 11:1
			CP.		Air Blank Cal Check Air Blank Cal Check Air Blank	0.000 0.081 0.000 0.081 0.000	11:17 11:17 11:17 11:19 11:19
Expiration:	2/5/2023		CP.		Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check	0.000 0.081 0.000 0.081 0.000 0.000	11:1 11:1 11:1 11:1 11:1 11:1
Expiration:	2/5/2023		CP.		Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank	0.000 0.081 0.000 0.081 0.000 0.081 0.000	11:1: 11:1: 11:1: 11:1: 11:1: 11:1: 11:1: 11:1:
Expiration:			CP.		Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check	0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081	11:1 11:1 11:1 11:1 11:1 11:1 11:1 11:
Expiration:	2/5/2023		CP.		Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank	0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000	11:11 11:12 11:12 11:14 11:14 11:14 11:14 11:17 11:17 11:17
Expiration: 0.080 BAC (F	2/5/2023	uments used to determ	CP.		Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check	0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081	11:11 11:11 11:12 11:12 11:12 11:14 11:14 11:17 11:17 11:17 11:17 11:17
Expiration:	2/5/2023	uments used to determ	CP.		Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank	0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000	11:11 11:12 11:12 11:14 11:14 11:14 11:14 11:17 11:17 11:17
Expiration: 0.080 BAC (F	2/5/2023	uments used to determ	CP.		Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Stats	0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081	11:1 11:1 11:1 11:1 11:1 11:1 11:1 11:
Expiration: 0.080 BAC (Fr Contents:	2/5/2023	uments used to determ) psig 70°F (21°C) Analytical	ine breath alcohol c		Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank	0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081	11:1 11:1 11:1 11:1 11:1 11:1 11:1 11:
Expiration: 0.080 BAC (F	2/5/2023 or the calibration of instr 105 Liters @ 1000	uments used to determ) psig 70°F (21°C) Analytical Accuracy	ine breath alcohol c	oncentration)	Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Stats Auerage	0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000	11:1 11:1 11:1 11:1 11:1 11:1 11:1 11:
Expiration: 0.080 BAC (F Contents: Component:	2/5/2023 or the calibration of instr 105 Liters @ 1000 Reported Concentration:	uments used to determ) psig 70°F (21°C) Analytical Accuracy (U, k=2):	ine breath alcohol c Analytical Method:	oncentration)	Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Stor Blank Cal Check Stats Average Std Dev	0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000	11:1 11:1 11:1 11:1 11:1 11:1 11:1 11:
Expiration: 0.080 BAC (F Contents: Component: Ethanol	2/5/2023 or the calibration of instr 105 Liters @ 1000 Reported Concentration: 208 ppm	uments used to determ D psig 70°F (21°C) Analytical Accuracy (U, k=2): +/-0.002 BAC(G/210L)	ine breath alcohol c Analytical Method:	oncentration)	Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Stor Blank Cal Check Stats Average Std Dev	0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000	11:1 11:1 11:1 11:1 11:1 11:1 11:1 11:
Expiration: 0.080 BAC (F Contents: Component:	2/5/2023 or the calibration of instr 105 Liters @ 1000 Reported Concentration:	uments used to determ) psig 70°F (21°C) Analytical Accuracy (U, k=2):	ine breath alcohol c Analytical Method:	oncentration)	Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Stor Blank Cal Check Stats Average Std Dev	0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000	11:1 11:1 11:1 11:1 11:1 11:1 11:1 11:
Expiration: 0.080 BAC (F Contents: Component: Ethanol	2/5/2023 or the calibration of instr 105 Liters @ 1000 Reported Concentration: 208 ppm	uments used to determ D psig 70°F (21°C) Analytical Accuracy (U, k=2): +/-0.002 BAC(G/210L)	ine breath alcohol c Analytical Method:	Distributed by: CMI Inc. 316 East Ninth Str	Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Stor Blank Cal Check Stats Average Std Dev	0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000	11:1 11:1 11:1 11:1 11:1 11:1 11:1 11:
Expiration: 0.080 BAC (F Contents: Component: Ethanol	2/5/2023 or the calibration of instr 105 Liters @ 1000 Reported Concentration: 208 ppm	uments used to determ D psig 70°F (21°C) Analytical Accuracy (U, k=2): +/-0.002 BAC(G/210L)	ine breath alcohol c Analytical Method:	Distributed by: CMI Inc. 316 East Ninth Str Owensboro, KY 4	Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Stor Blank Cal Check Stats Average Std Dev	0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000	11:1 11:1 11:1 11:1 11:1 11:1 11:1 11:
Expiration: 0.080 BAC (F Contents: Component: Ethanol	2/5/2023 or the calibration of instr 105 Liters @ 1000 Reported Concentration: 208 ppm	uments used to determ D psig 70°F (21°C) Analytical Accuracy (U, k=2): +/-0.002 BAC(G/210L)	ine breath alcohol c Analytical Method:	Distributed by: CMI Inc. 316 East Ninth Str Owensboro, KY 4. Phone 866-835-06'	Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Stor Blank Cal Check Stats Average Std Dev	0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000	11:1 11:1 11:1 11:1 11:1 11:1 11:1 11:
Expiration: 0.080 BAC (F Contents: Component: Ethanol	2/5/2023 or the calibration of instr 105 Liters @ 1000 Reported Concentration: 208 ppm	uments used to determ D psig 70°F (21°C) Analytical Accuracy (U, k=2): +/-0.002 BAC(G/210L)	ine breath alcohol c Analytical Method:	Distributed by: CMI Inc. 316 East Ninth Str Owensboro, KY 4	Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Stor Blank Cal Check Stats Average Std Dev	0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000	11:1 11:1 11:1 11:1 11:1 11:1 11:1 11:
Expiration: 0.080 BAC (F Contents: Component: Ethanol	2/5/2023 or the calibration of instr 105 Liters @ 1000 Reported Concentration: 208 ppm	uments used to determ D psig 70°F (21°C) Analytical Accuracy (U, k=2): +/-0.002 BAC(G/210L)	ine breath alcohol c Analytical Method:	Distributed by: CMI Inc. 316 East Ninth Str Owensboro, KY 4. Phone 866-835-06'	Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Air Blank Cal Check Stats Auerage Std Deu Rel Std Deu(%)	0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000 0.081 0.000	11:1 11:1 11:1 11:1 11:1 11:1 11:1 11:

*Traceable to: Certified Reference Material - 262.4 μmol/mol Ethanol in Nitrogen - Serial No. GN0015026 Lot No. 050319ΕΙΙ

as 200

Specialty Gas Lab Tech

100

Store in dry area, away from sources of heat, ignition and direct sunlight. Do not allow storage area to exceed 52 °C (125 °F).

01-25-2021 Issuance Date



The calibration results within this certificate were obtained using equipment and standards capable of producing analytical results traceable to NIST, and apply only to the items contained on this certificate. ILMO Products Company makes no warranty or representation as to the suitability of the use of any information provided for any particular purpose. The information use is at the sole discretion and risk of the user. Liability shall be limited to established replacement cost of this material or service.

ISO/IEC 17025:2017 Accredited Laboratory

Service Evaluation Form This form MUST be completed and enclosed with instrument to be serviced. Failure to complete and return this form may cause delays in service. (Note: please ship items in their original shipping container(s) or a similar protective box.) (contact Customer Service) 1. Contact information: Customer Number 731362 Phone: (405) 425 2460 Name: **Board of Tests** christopher.pape@bot.ok.gov Email: Fax: (405) 425 2490 2. Bill to Address: Board of Tests Ship to Address: Board of Tests 3600 N ML King, Bldg 9 3600 N ML King, Bldg 9 OKC, OK 73111 OKC, OK 73111 3. Serial Number: 80-005042 Instrument Model: Intoxilyzer 500 / 8000 Detailed Description of Problem: Unit's source failed snap test during evaluation. Unit's Flow sensor is at 122. Unit's Atmospheric sensor has a high CF at 1.14 ***Hazardous Material Warning! - DO NOT return gas cylinders with instrument!*** X I Authorize Repairs Up To: All \$250 \$500 \$750 Other \$_____ Purchase Order Number (attach a copy of P.O. if applicable) _ Ship to: Authorized By: Chris Pape Breath Testing Admin Name (Please Print) Title CMI, Inc Attn: Service Dept. 10/19/2020 316 East Ninth Street Signature Date Owensboro, KY 42303 Call: 405-425-2460 or email No, please send estimate before repairs are made.

CMI INC.

Note: An estimate will be faxed <u>before</u> performing any repairs and may cause delays in service. An evaluation fee (\$79.00 or actual costs) will apply to estimates that are not repaired.





316 E 9th Street / Owensboro KY 42303 / USA Phone: 866-835-0690 Fax: 270-685-6268

Ship To:

USA

Oklahoma Board Of Tests

Blvd Adcu Annex Bldg #9

Oklahoma City OK 73136-3000

3600 Martin Luther King

Phone: 405-425-2424

SERVICE WORK ORDER

BILL CODE: Out of Warranty

405441

Bill To:

DATE: 10/26/2020

Oklahoma Board Of Tests Blvd Adcu Annex Bldg #9 3600 Martin Luther King Oklahoma City, OK 73136-3000 USA Customer #: 731362 Phone: 405-425-2424 Fax: 405-425-2490 fax

Serial Number: 80-005042

MODEL #: 002480OK EXTRA PARTS RCVD:

N/A

DESCRIPTION OF PROBLEM

UNIT'S SOURCE FAILED SNAP TEST DURING EVALUATION. UNIT'S FLOW SENSOR IS AT 122. UNIT'S ATMOSPHERIC SENSOR HAS A HIGH CF AT 1.14

WORK PERFORMED:

source/ erratic; flow sensor/ erratic; tygon tubing/ dirty; O-ring/ cracked; replaced parts listed; calibrated/ final tested;

PARTS USED

Seq. No.	Part	Description	Quantity		
50	210120	SENSOR,SMD,PRESSURE	1.00 EA		
60	471201	O-RING,BUNA-N,-008,3/16inx5/16	1.00 EA		
80	650517	CERTIFICATE OF CALIBRATION	2.00 EA		
90	441169	COVER DUST,5/8in x 1/2in	1.00 EA		
100	470145	CAP, PLASTIC, 25IDx.50LG, RED	3.00 EA		
110	470154	CAP,PLSTC,.406IDx.438562,RED	1.00 EA		
120	470218	SEAL, TAMPER EVIDENT, SERVICE	1.00 EA		
130	530030	TUBING, TYGON, .375 IDx. 562 OD	2.00 IN		
140	020982FL	SOURCE ASSY	1.00 EA		

LABOR / TESTING		
Misc Code	Description	Hours
LABR	Service Repair Labor	
LCAL	Service - Calibration Adjust	0.50
LFT	Service - Cal / Final Test	1.00

Technician Name: Dewayne Varvel



OFFICER'S AFFIDAVIT AND NOTICE OF REVOCATION/DISQUALIFICATION

11/16

Badge #:

	RREST TIME		COUNTY #	C	ΓY #	C	ITATION	#		
ARREST LOCATION			CITI	Y				OK	Y LAHC	MA
SUBJECT NAME NORTH, GAYLON LEVER	NE	A.C			/09/1	961	SEX M	W	HEIGHT	Contraction of the second
ADDRESS 3600 MARTIN LUTHER	KING A	OKLAH	OMA CITY	ζ			OK		ZIP COL 73	DE 133
DRIVER LICENSE # T083454871			PIRATION DATE	19	STATE OK	CLAS	S	RESTR	ICTIONS	ENDORSEMENT
VEHICLE MAKE MOCK	MODEL				G # MOCK			OK	CMV N	PLACARD N

On the above date, time, and location, the above named person was arrested, and I had reasonable grounds to believe the person was driving or in actual physical control of a motor vehicle upon a public road, highway, street, turnpike, public parking lot, or other public place, or upon a private road, street, alley or lane which provides access to one or more single or multi-family dwelling within this state while under the influence of alcohol and/or other intoxicating substances as prohibited by law. (Describe driving behavior or circumstances):

MOCK SUBJECT TEST

Describe person's condition (odor, actions, etc.):

POST REPAIR MAINTENANCE, INSTALLED NEW CYLINDER

THE DEDCON	WAS INFORMED OF	THE IMPLIED	CONCENT	TECT DECH	FOT
THE FERSON	WAS INFORMED OF	I HE IMPLIED	CONSERT	1 COL REQU	6.51

UIPMENT [-8000			NUMBER SPECIALIST		
andard Dry			0.080	PERMIT NUMBER	
		wa nin temperat	00521080A2	DATE CYLINDER INSTALLED 03/26/2021	DATE CYLINDER EXPIRES
215 C	ST DATE 03/26/2021	G		ater) deprivation/observation	1
SULT TABLE Test	g/210L	Time		ed prior to and continuing of test administration.	OF THE STATE
Air Blank Subject Test Breath Volume Air Blank Nait Air Blank Subject Test Breath Volume Air Blank Sal Check	0.00 0.00 0.00 1.527 LITERS 0.00 0.08	12:38 12:38 12:38 12:40 12:40 12:43 12:44 12:44	compliance with		Test Result
Air Blank	0.00	12:44	BADGE # 0000000484 AGENCY BOARD OF TES	PERMIT # 0000023020	BrAC g/210L 0.00

SERVING OF NOTICE: I personally hand delivered a copy of the above containing the Notice of Revocation and the Notice of Department Impaired Driving Accountability Program to the person named in Section 1 above.

Date served:	03/	26	2021	
Date served:	057	201	AUAL	Signature of office

Date: 03/26/2021 Place (location when sig	ned):	Signature of arrest	ing officer:
PAPE, CHRISTOPHER		BOARD OF TESTS	
Addee Agency address INTER-AGENCY			CONTRACTOR OF A DECIMAL OF A
THER WITNESSES:	1.512		Sectio
I NAME	TITLE	ADDRESS	PHONE
2 NAME	TITLE	ADDRESS	PHONE

OFFICER'S AFFIDAVIT AND NOTICE OF REVOCATION/DISQUALIFICATION

12/16

Badge #:

ARREST DATE ARREST TI 03/26/2021 1200		COUNTY #		CITY #	CI	TATION #			
ARREST LOCATION BOT			CITY		-		OK	LAHO	AMA
NORTH, GAYLON LEVERNE				E OF BIRTH 8/09/1	961	SEX M	W	HEIGHT	
ADDRESS 3600 MARTIN LUTHER KING	A	OKLAHOMA CI	TY			OK		ZIP COI	DE 133
DRIVER LICENSE # T083454871		EXPIRATION DATE 06/30/		OK STATE	CLASS		RESTR	ICTIONS	ENDORSEMENT
VEHICLE MAKE MOCK	MOCK		_	MOCK			OK	^{CMV} N	PLACARD N

On the above date, lime, and location, the above named person was arrested, and I had reasonable grounds to believe the person was driving or in actual physical control of a motor vehicle upon a public road, highway, street, tumplke, public parking lot, or other public place, or upon a private road, street, alley or lane which provides access to one or more single or multi-family dwelling within this state while under the influence of alcohol and/or other intoxicating substances as prohibited by law. (Describe driving behavior or circumstances):

MOCK SUBJECT TEST

Describe person's condition (odor, actions, etc.):

POST REPAIR MAINTENANCE, INSTALLED NEW CYLINDER

THE PERSON Y	WAS INFO	RMED OF	THE IMPL	IED CONSENT	TEST REOUES	T

QUIPMENT			NUMBER - 005042	BOARD	OF TESTS	
TANDARD		1 * *	TARGET VALUE		PERMIT NUMBER	
Dry			0.080		23020	
ANUFACTURER ILMO			005210	080A2	DATE CYLINDER INSTALLED 03/26/2021	DATE CYLINDER EXPIRES 02/05/2023
L215 03	DATE 3/26/2021	G			r) deprivation/observation	
SULT TABLE Test	g/210L	Time	and the second se		prior to and continuing est administration.	OL THE STATE
Air Blank Subject Test Breath Volume Air Blank Wait Air Blank Subject Test Breath Volume	0.00 0.00 1.964 LITERS 0.00 0.00 0.00 1.527 LITERS	12:38 12:38 12:38 12:40 12:40 12:43	complianc	e with th dure of th	Breath-Alcohol Test in e applicable operating ne Board of Tests. 23076	
Air Blank Cal Check Air Blank	0.00 0.08 0.00	12:44 12:44 12:44	BADGE # 00000000 AGENCY	RISTOPH 484 7 TESTS	PERMIT # 0000023020	Test Result BrAC g/210L 0.00

as a result of this arrest by filing a petition in the District Court of the County in which you were arrested. Your petition must be filed within thirty (30) days of the date listed below. You may be eligible to participate in the Impaired Driver Accountability Program (IDAP) administered by the Department of Public Safety (DPS). Completion of IDAP may prevent a revocation as a result of this arrest from appearing on your driving record. Participation in IDAP may reduce the amount of fees you will be required to pay to the State. You must request IDAP within thirty (30) days of the date listed below. Mail the request for IDAP to the Department of Public Safety at P.O. Box 11415, Oklahoma City, OK 73136.

SERVING OF NOTICE: I personally hand delivered a copy of the above containing the Notice of Revocation and the Notice of Department Impaired Driving Accountability Program to the person named in Section 1 above.

Date served:	03/	26/	2021	Signature of officer:
Barb Partes		the second se	concern and the second second	orginatorie er erneer

In accordance with Title 12 O.S. Section 426, "I state under penalty of perjury under the laws of Oklahoma that the foregoing is true and correct." Section 4

Date: 03/26/2021 Place (location when signed):	Signature of arresting officer:
PAPE, CHRISTOPHER	BOARD OF TESTS
BADGE AGENCY ADDRESS INTER-AGENCY	ZIP CODE PHONE 000000 4054252460

Certificate of	Analysis – Wet Bath		
8/11/2020	0.020		
Date of Analysis	Labellad tar (Value (g/210L)		
20070	2020	Nonace and	
Lot Number 🤄	Average that regult (72401 17 1-	ABORATORIES,	INC.
christopher Pch		HARRISBURG, PA 17111- 4511	
BOT Technician Name and Signature	V	BRT Form CO.	

BRT Form 6.0

CERTIFICATE

In accordance with BRT - 2.2.0, Maintenance, the above referenced simulator solution is suitable for use as an external reference in maintenance of the intoxilyzer 8000.

ALPHA INSTRUMENT Z Intoxilyzer - Alcohol Analyzer Model 8000 08/11/2020

g/210L

0.000

0.020

0.000

0.020

0.000

0.020

0,000

0.020

0,000

0.021

0.000

0.020 .

0.000

0.020

0,000

0.020

0.000

0.020

0,000

0.019

0.000

0,0200

Test

Air Blank

Cal Check

Air Blank

Cal Check Stats Average

SN 80-004923 09:42

54:60

09:43

09:43

09:44

09:45

09:45

09:46

09:47

09:47

09:48

09:48

09:49

09:49

09:50

09:51

09:51

09:52

09:53

09:53

09:54

09:54

Certified Alcohol Reference Solution for Simulator Time

OF

Random Samples of Lot Number 20070 of Alcohol Reference Solution for Simulator were analyzed by chromatography on February 17, 2020, using a Perkin Elmer Gas gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.0242% (w/vol) ethyl alcohol. The expiration date for this lot number is February 13, 2022 at 11:59 PM.

ANALYSIS

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.020 g/210L +/- 3%.

Std Dev 0.0004 Rel Std Deu(%) 2,3569

Operator's Signature

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN08031603 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.

13/16

Certificate o	of Analysis – Wet Bath	A CONTRACTOR	
7/22/2020	0.040		
Date of Analysis	Labelled targes value (g/210L)		
20060	10.040	1992[D011	
Lot Number	Arera o vat rawit (grant)	LABORATORIES,	INC.
christopher 19	CL INDRINGTON	STREET . HARRISBURG, PA 17111- 4511	• TELEPHONE: 717-564-5470
BOT Technician Name and Signatu		BRT Form 6,0	

In accordance with BRT – 2.2.0, Maintenance, the above referenced simulator solution is suitable for use as an external reference in maintenance of the Intoxilyzer 8000.

CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

RCVD 6/1/20 CP

ALPHA UNIT SUB		
Intoxilyzer - Alcohol	Analyzer	
Nodel 8000	SN	80-00340
07/22/2020		14:5

g/210L

0.000

0.040 0.000

0.039

0.000

0.039

0.000

0.039

0.000

0.039

0.000

Test

Air Blank

Cal Check

Air Blank

Time

14:57

14:58

14:59

14:59

15:00

15:01

15:01

15:02

15:02

15:03

15:04

Random Samples of Lot Number 20060 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on February 11, 2020, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.0483% (w/vol) ethyl alcohol. The expiration date for this lot number is February 10, 2022 at 11:59 PM.

When used in a calibrated Simulator, operating at $34^{\circ}C$ +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.040 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

aul

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability; Testing was conducted using Cerilliant Reference Standard lot number FN12181501 whose values are traceable to NIST. All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.

Cal Check 0.039 15:04 Air Blank 0.000 15:05 Cal Check 0.040 15:06 Air Blank 0.000 15:06 Cal Check 0.040 15:07 Air Blank 0.000 15:07 D.040 15:08 Cal Check Air Blank 0.000 15:09 0.040 15:09 Cal Check Air Blank 0.000 15:10 Cal Check Stats Average 0.0395 0.0005 Std Dev Rel Std Dev(%) 1,3342

Signature

ABORATORIES, INC.
TOET BRE BELLE AND RISE URG, PA 17111- 4511 @ TELEPHONE: 717-564-5470
simulator solution is
/zer 8000.

CERTIFICATE OF ANALYSIS

RLPHA INSTRUMENT Intoxilyzer - Alcohol Analyzer Model 8000 5N 80-002591 Certified Alcohol Reference Solution for Simulator 06/16/2020 16:18

TACA.	g/210L	Tine
AIr Blank	0,000	16:18
Cal Check	0.099	16:19
Air Blank	0.000	16:20
Cal Check	0.098	16:20
Air Blank	0.000	16:21
Cal Check	0.097	16:22
Air Blank	0.000	16:22
Cal Check	0.098	16:23
Air Blank	0.000	16:24
Cal Check	0.098	16:24
Air Blank	0.000	16:25
Cal Check	0.098	16:26
Air Blank	0.000	16:26
Cal Check	0.098	16:27
Air Blank	0.000	16:27
Cal Check	0.097	16:28
Air Blank	0.000	16:29
Cal Check	0.098	16:29
AIr Blank	0.000	16:30
:al Check	0.098	16:31
lir Blank	D. DOO	16:31
'al Check Stats		
Average	0.0979	
Std Dev	0.0005	
Rel Std Deu(%)	0.5798	

Random Samples of Lot Number 20190 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on April 8, 2020, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1206% (w/vol) ethyl alcohol. The expiration date for this lot number is April 6, 2022 at 11:59 PM.

When used in a calibrated Simulator, operating at $34^{\circ}C$ +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

ignature Operator

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability: Testing was conducted using Cerilliant Reference Standard lot number FN02271802 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.

15/16

Lawied tont value (g/210L)		
1 Ma 0.197	-110X.com	
Av ar leat sur mar Hin - K	ABORATORIES,	INC.
	TREET O HARRISBURG, PA 17111- 4511	
Signature Chal303	BRT Form 6.0	

In accordance with BRT – 2.2.0, *Molntenance*, the above referenced simulator solution is suitable for use as an external reference in maintenance of the Intoxilyzer 8000.

00200

Certificate of Analysis - Wet Bath

07/2020

20160

Christopher BOT Technician Na

Lot Nu

CERTIFICATE OF ANALYSIS

12020

Certified Alcohol Reference Solution for Simulator

ox11yzer - f e1 8000 07/2020		5N 80-003400 15;25
	g/210L	Tine
Blank	0.000	15:26
heck	0.196	15:26
lank	0,000	15:27
heck	0.196	15:28
lank	0.000	15:28
heck	0.197	15:29
lank	0.000	15:29
heck	0.196	15:30
lank	0.000	15:31
leck ·	0,198	15:31
ank	0.000	. 15:32
leck	0.198	15:33
lank	0.000	15:33
leck	0,197	15:34
lank	0.000 .	15:35
neck.	0.198	15:35
lank	0.000	15:36
leck	0.198	15:37
lank	0.000	15:37
neck.	0.197	15:38
ank	0.000	15:38
eck Stats		
ige ,	0.1971	
Deu	0.0008	
td Deu(2)	0.4442	

Random Samples of Lot Number 20160 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on March 20, 2020, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.2427% (w/vol) ethyl alcohol/ The expiration date for this lot number is March 18, 2022 at 11:59 PM.

When used in a calibrated Simulator, operating at $34^{\circ}C$ +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.200 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN08101505 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.

16/16