

## STATE OF OKLAHOMA BOARD OF TESTS FOR ALCOHOL AND DRUG INFLUENCE

Dr. Jarrad Wagner, Ph.D. F-ABFT Chairman Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

## ATTESTATION

I, the undersigned *Joshua Smith* (Director/Records Custodian), attest under oath pursuant to 12 O.S. §2902 (11)(a) this record was made, at or near the time of the occurrence of the matters set forth by or from information transmitted by a person having knowledge of those matters; was kept in the course of the regularly conducted business activity of the Board; and was made pursuant to the regularly conducted activity. The duplicate record provided is listed below:

1) The attached "Service Overview" is a true and accurate copy of maintenance performed on the Intoxilyzer 8000, serial number 80-005351, in compliance with the agency's Maintenance Policy and Procedures.

This packet consists of 13 pages.

Signed \_\_\_\_\_

Name of Position: Director/Records Custodian

Date of Attestation: 05/24/2022



Board of Yests for Alcohol and Drug Influence Entrancing public safety by ensuing the accuracy and sciuntific metability of blood and breath resting



3600 N. Martin Luther King Bldg #9 Oklahoma City, OK 73111 405-425-2460

#### **BENCH CHECK DATE:** <u>04-04-2019</u>

SERIAL: 80-005351

## DOCUMENTATION WITHIN SERVICE PACKET:

- BOT Certificate of Calibration and Operation
- Bench Check
- Bench Check Technician's Report
- BOT Dry Gas Certificate of Analysis
- Mock Subject Test
- BOT Certificates of Analysis on Certified Solutions (4)



## STATE OF OKLAHOMA BOARD OF TESTS FOR ALCOHOL AND DRUG INFLUENCE

Administrative Office at: Board of Tests for Alcohol and Drug Influence Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307

# **Certificate of Calibration and Operation**

This is to certify that the calibration of **INTOXILYZER 8000, serial number** <u>80-005351</u>, manufactured by CMI Inc., a subsidiary of MPD, Inc. of Owensboro, Kentucky, was tested and found to be in compliance with National Highway Traffic Safety Administration Standard for Devices to Measure Breath Alcohol (F.R., No. 179 48705-48710 Sept. 17, 1993) and with International Recommendation OIML R 126 Edition 1998 (E). Calibration solutions are traceable to NIST standard material 1828.

In addition the operational procedures and functions of this Intoxilyzer have also been tested and found to be in compliance with the Oklahoma Board of Tests for Alcohol and Drug Influence - Rules and Regulations OAC Title 40: 25-1-2 and 30-1-3 and is hereby certified for use as an evidential breath testing device within this state.

DATE
08-10-11
08-23-11
08-06-13
06-25-15
07-17-15
05-02-17
04-04-19



Presentation of this form certifies that the Intoxilyzer, SN: <u>80-005351</u>, functioned properly at the time of the breath test and hereby certifies the breath test result as valid.

Re: Request for Log of Test on Intoxilyzer 8000

Name:

Dear Sir:

The Log of Test information from a breath test conducted on a Intoxilyzer 8000 is printed in Section 2 of the Affidavit and is no longer retained as a separate form as with the Intoxilyzer 5000. The 8000 uses a nitrogen/ethanol dry gas for the external control test. The dry gas cylinder has the capacity to perform approximately 325 tests and does not require maintenance until its expiration date or can no longer provide the required pressure to perform a control test. The cylinder has an expiration date of two years from date filled at the manufacturer and is programmed in the 8000 upon installation, both dates are printed in Section 2 of the Affidavit. The Intoxilyzer is programmed to disable itself when the date of expiration is reached or the pressure monitor indicates 25 lbs psi.

Due to the reliability and stability of dry gas, maintenance on the 8000 is not required every 30 days or 25 tests, therefore maintenance is no longer performed by officers in the field. The Board of Tests for Alcohol and Drug Influence certifies each Intoxilyzer 8000 for field use by performing a bench check. The bench check is designed to check the operation and procedures programmed in the 8000 to ensure it is functioning properly for use in Oklahoma. A bench check is performed every time a new dry gas cylinder is installed or the instrument has returned from the manufacturer for repair. Once the instrument passes the bench check, all that is required to set it up for field use is to plug it in and turn it on and does not require a permit to set up. The Board of Tests has the ability to monitor each instrument once it is set up in the field.

Each Intoxilyzer has a database of valid permit numbers that it checks prior to allowing a breath test to be performed. The database is updated on a regular basis, if a permit number is invalid that officer will not be able to conduct a breath test on the 8000. If the 8000 prints an affidavit the breath test operator listed in Section 2 had a valid permit at the time of the test and the breath test is valid.

If you have any further questions please feel free to call me.

Sincerely,

pohert Blakeberry

J. Robert Blakeburn (OBA11059) State Director of Tests for Alcohol and Drug Influence

## STATE OF OKLAHOMA - DEPARTMENT OF PUBLIC SAFETY CMI INTOXILYZER BENCH CHECK PROCEDURE REPORT

NCH CHECK DATE 04/04/2019	TIME PROCEDURE BI	E E 1 (1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	STRUMENT LOCATION $80 - 00535$	1	
TRUMENT TYPE	UMENT TYPE INSTRUM			INSTRUMENT SE	
		TEST RES	SULTS		
Diagnosti PASS					
Mouth Alcoho PASS					
RFI Test	t				
PASS		Dur	licate Breath	Sampling Re	sults
Abort Te	st	- Dup			
PASS		Value	Test 1	Test 2	Control Test
Improper Samp PASS	철수 가슴 하는 것 같아.	. 02	0.017	0.018	0.019
Interferent Dete	ect Test	.04	0.038	0.037	0.038
PASS		.10	0.099	0.097	0.099
No Sample Giv PASS	en Test	.20	0.199	0.198	0.202
Insufficient Sam	nple Test		PA	SS	
PASS Air Blank Contami PASS	nation Test				
Range Exceede	ed Test				
T 1 1	ć i du	CERTIFYING TI			
I swear under penalty the Intoxliyzer Bencl					
				DATE	
CK					04/04/201
E					



3600 N. Martin Luther King Bldg #9 Oklahoma City, OK 73111 405-425-2460

#### **BENCH CHECK TECHNICIAN'S** REPORT

Date: 04 / 04 / 19 Start Time: 0900 End Time: 1030

INSTRUMENT Intoxilyzer Model Serial # 8000 80-005351 GAS CANISTER LOT # EXP Date 00919980AZ 03-05-2021 1007 REFERENCE Simulator Model Guth 2100 Guth 2100 Guth 2100 Guth 2100 

03-11-19	03-11-19	03-11- 19	03-11-19	
	Solution Con	nmission Date		
01-09-2020	12-06-2019	02-26-2020	10-09-2019	
	Expirat	ion Date		
01-09-2018	12-06-2017	02-26-2018	10-09-2017	
	Manufact	tured Date	and the second	
18020	17410	18070	17340	
	LC	<u>DT #</u>		
.02	.04	.10	.20	
	Conce	ntration		
DR3591	DR3753	DR3753 DR3594		
	Sei	rial #	the second second second	
Guti 2100	Guti 2100	Guil 2100	Gutt 2100	

## **REASON FOR BENCH CHECK**

POST REPAIR CYLINDER REPLACEMENT TROUBLESHOOTING ROUTINE MAINTENANCE

COMMENTS

Technician Signature / Permi sh.

Printed Name

6/13

Certificate	of Analysis – Dry Gas	
04-01-19	,080	
Date of Analysis	Labelled target value (g/210L)	
1007	,001	
Cylinder#	Average test result	
JOSH SMITH		
BOT Technician Name and Signatu	re	specialty gases
In accordance with BOT Polic referenced dry gas cylinder is s	y and Procedure Statement BRT-2.1.0, Mu uitable for use as an external reference and	

0.080 BAC (For the calibration of instruments used to determine breath alcohol concentration)

Analytical

Accuracy.

(U, k=2):

[5.2 ppm]

+/-0.002 BAC(G/210L) NDIR

105 Liters @ 1000 psig 70°F (21°C)

ALPHA INSTRUMENT Intoxilyzer - Alcohol Analyzer SN 80-002591 Model 8000 09:56 04/01/2019

Test	g/210L	Time
Air Blank	0.000	09:56
Cal Check	0.080	09:57
Air Blank	0.000	09:57
Cal Check	0.080	09:58
Air Blank	0.000	09:58
Cal Check	0.081	09:58
Air Blank	0.000	09:59
Cal Check	0.081	09:59
Air Blank	0.000	10:00
Cal Check	0.081	10:00
Air Blank	0.000	10:01
Cal Check	0.081	10:01
Air Blank	0.000	10:02
Cal Check	0.080	10:02
Air Blank	0.000	10:02
Cal Check	0.080	10:03
Air Blank	0.000	10:03
Cal Check	0.080	10:04
Air Blank	0.000	10:04
Cal Check	0.081	10:04
Air Blank	0.000	10:05
Cal Check Stats		
Average	0.0805	
Std Dev	0.0005	
Rel Std Deu(%)	0.6547	



\*NIST Traceable Reference Material Cylinder No. CC274523 / Job No. 09160306 Certified 362.2 µmol/mol Ethanol in Nitrogen

11722

105L

Reported

208 ppm

Balance

Concentration:

BAC105L080T

00919080A2

3/5/2021

Certificate ID:

Cylinder Size:

Lot Number:

Expiration:

Contents:

Component:

Ethanol.

Nitrogen

Part #:

Specialty Gas Lab Tech

Store in dry area, away from sources of heat, ignition and direct sunlight. Do not allow storage area to exceed 52 °C (125 °F).

Distributed by:

CMI Inc.

01-24-Date

217-245-2183 • Fax: 217-243-7634 • www.ilmoproducts.com

**Certificate of Analysis** 

RECVD - 03/05/2

Analytical

Method:



The calibration results within this certificate were obtained using equipment and standards capable of producing analytical results traceable to NIST, and apply only to the items contained on this certificate. ILMO Products Company makes no warranty or representation as to the suitability of the use of any information provided for any particular purpose. The information use is at the sole discretion and risk of the user. Liability shall be limited to established replacement cost of this material or service.

ISO/IEC 17025:2005 Accredited Laboratory

## OFFICER'S AFFIDAVIT AND NOTICE OF REVOCATION/DISQUALIFICATION

8/13

Badge #:

			C								Section
ARREST DATE	ARREST TIME	······	COUNTY #	CI	TY #	C	ITATION #				
04/04/2019	0930		55				40				
			CI		7			COUNT			
MOCK			1	MOCH	ξ			MO	CK		
SUBJECT NAME					OF BIRTH		SEX	RACE	HEIGHT	1	WEIGHT
MOCK, SUBJECT TEST				08	/09/	1961	M	W	510	)	165
ADDRESS		CITY				S	TATE		ZIP COD	DE	
3600 MARTIN LUTHER	KING A	OKLAH	OMA CIT	Y		1	OK		731	133	
DRIVER LICENSE #			PIRATION DATE		STATE	CLASS	3	RESTRI	CTIONS	ENDO	RSEMENTS
T083454871		(	06/30/2	019	OK	A				M	
VEHICLE MAKE	MODEL			TA	G #		ST	ATE	CMV	PLAC	ARD
NONE	NONE	3		]	NONE			ΟK	N	N	
On the above date, time, and location, the above name highway, street, turnpike, public parking lot, or other pul											
influence of alcohol and/or other intoxicating substance						,				nino ai	

(Describe driving behavior or circumstances):

MOCK SUBJECT TEST. NEW GAS CANISTER INSTALL.

Describe person's condition (odor, actions, etc.):

ROUTINE MAINTENANCE.

BREATH-ALCOHOL ANALYSIS	SRECORD, REPORT OF TH	IE PERSON NAMED	IN SECTION 1 ABOVE, and LOG OF TEST AND MAINTENANCE RECORD.	Section 2
			NUMBER SPECIALIST	
I-8000		80 -	-005351   BOARD OF TESTS	
STANDARD			TARGET VALUE PERMIT NUMBER	
Dry	λ.		0.080 16837	
MANUFACTURER			LOT # DATE CYLINDER INSTALLED	DATE CYLINDER EXPIRES
ILMO			00919080A2 04/04/2019	03/05/2021
	TDATE 4/04/2019	MOUTHPIECE G	A 15-minute (or greater) deprivation/observation	
RESULT TABLE			period was observed prior to and continuing	THE STAN
Test	g/210L	Time	to the time of test administration.	
Air Blank	0.00	10:14	I administered said Breath-Alcohol Test ir	
Subject Test Breath Volume	0.00 2.027 LITERS	10:15	compliance with the applicable operating	
Air Blank	0.00	10:15	procedure of the Board of Tests.	
Wait		10:17	1	
Air Blank	0.00	10:17		
Subject Test	0.00	10:18		1007
Breath Volume	1.316 LITERS		SIGNATURE OF OPERATOR	1901
Air Blank	0.00	10:18		
Cal Check	0.08	10:18	NAME	Test Result
Air Blank	0.00	10:19	SMITH, JOSHUA	BrAC g/210L
			BADGE # PERMIT #	Diric g/210E
			000000435 0000016837	-   0.00
			BOARD OF TESTS	
			BOARD OF TESTS	ked or denied for a Section 3

EFFECTIVE THIRTY (30) DAYS FROM DATE SERVED. Notice of Revocation is given that your driver's license is revoked or denied for a period of one hundred eighty (180) days or more, if you refused or failed the state's test(s). In addition, your commercial driving privilege in this state will be disqualified if you refused or failed the state's test(s) (operating a CMV or CDL holder) OR a test result of .04 or more while operating a commercial motor vehicle. To appeal the driver's license revocation/disqualification, the Department of Public Safety must receive your written request for an Administrative Hearing within fifteen (15) days from the date you received notice of revocation or disqualification. Submit your written request in person at 3600 N. MLK Ave. or by mail to Legal Division, Department of Public Safety, P.O. Box 11415, OKC, OK, 73136

SERVING OF NOTICE: I personally hand-delivered a copy of the above containing the Notice of Revocation to the person named in Section 1 above.

Date served: 04/04/2019	Signature of officer:
-------------------------	-----------------------

In accordance with Title 12 O.S. Section 426, "I state under penalty of perjury under the laws of Oklahoma that the foregoing is true and correct." Section 4 Date: 04/04/2019 Place (location when signed): Signature of arresting officer: NAME AGENCY SMITH, JOSHUA BOARD OF TESTS BADGE AGENCY ADDRESS ZIP CODE PHONE 000000435 INTER-AGENCY 00000 "4054252460" OTHER WITNESSES: Section 5 NAME TITLE ADDRESS PHONE 1 NAME TITLE ADDRESS PHONE 2

## OFFICER'S AFFIDAVIT AND NOTICE OF REVOCATION/DISQUALIFICATION

9/13

			-			7				Section
ARREST DATE	ARREST TIME		COUNTY #	0	CITY #	C	ITATION #			
04/04/2019	0930		55				40			
ARREST LOCATION			-	CITY				COUNT	Y	
MOCK				MOC	K			MO	CK	
SUBJECT NAME					OF BIRTH		SEX	RACE	HEIGHT	WEIGHT
MOCK, SUBJECT TEST				08	3/09/1	961	M	W	510	) 165
ADDRESS		CITY				S	TATE		ZIP COD	DE
3600 MARTIN LUTHER	KING A	OKLAHO	DMA CI'	ΤY			OK		731	L33
DRIVER LICENSE #		EX	PIRATION DATE		STATE	CLASS	3	RESTRI	CTIONS	ENDORSEMEN
T083454871		C	06/30/2	2019	OK	A				M
VEHICLE MAKE	MODEL			1	AG #		ST/	ATE	CMV	PLACARD
NONE	NONE				NONE			ΟK	N	N
On the above date, time, and location, the above name	d person was arrested, and I had r	easonable arounds to	believe the ner	con was dr	iving or in actus	hobycical	control of a	motor	hiele upon	a public road

On the above date, time, and location, the above named person was arrested, and I had reasonable grounds to believe the person was driving or in actual physical control of a motor vehicle upon a public road, highway, street, turnpike, public parking lot, or other public place, or upon a private road, street, alley or lane which provides access to one or more single or multi-family dwelling within this state while under the influence of alcohol and/or other intoxicating substances as prohibited by law. (Describe driving behavior or circumstances):

MOCK SUBJECT TEST. NEW GAS CANISTER INSTALL.

Describe person's condition (odor, actions, etc.):

ROUTINE MAINTENANCE.

THE PERSON WAS	S INFORMED	OF THE IMPLIED	CONSENT TE	ST REQUEST
----------------	------------	----------------	------------	------------

	RECORD, REPORT OF TH	E PERSON NAMED	IN SECTION 1 ABOVE, and LOG OF	TEST AND MAINTENANCE RECORD.	Section 2
EQUIPMENT			NUMBER SPECIALIST		
I-8000		80-	-005351  BOARD	OF TESTS	
STANDARD			TARGET VALUE	PERMIT NUMBER	
Dry			0.080	16837	
MANUFACTURER			LOT #	DATE CYLINDER INSTALLED	DATE CYLINDER EXPIRES
ILMO			00919080A2	04/04/2019	03/05/2021
	TATE 1/04/2019	MOUTHPIECE G	A 15-minute (or greate	er) deprivation/observation	
RESULT TABLE	1/01/2019		period was observed	l prior to and continuing	
Test	g/210L	Time	to the time of t	est administration.	OF THE STATE
Air Blank	0.00	10:14	I administered said	Breath-Alcohol Test in	
Subject Test	0.00	10:14			
Breath Volume	2.027 LITERS	10:15	compliance with th	e applicable operating	
Air Blank	0.00	10:15	procedure of th	ne Board of Tests.	A +++
Wait		10:17	1		
Air Blank	0.00	10:17			
Subject Test	0.00	10:18			1007
Breath Volume	1.316 LITERS		SIGNATURE OF OPERATOR		1907
Air Blank	0.00	10:18			
Cal Check	0.08	10:18	NAME		Test Result
Air Blank	0.00	10:19	SMITH, JOSHUA	PERMIT #	BrAC g/210L
			0000000435	0000016837	
			AGENCY	0000010037	0.00
			BOARD OF TESTS		
				nat vour driver's license is revoke	d or denied for a Section 3

EFFECTIVE THIRTY (30) DAYS FROM DATE SERVED. Notice of Revocation is given that your driver's license is revoked or denied for a second period of one hundred eighty (180) days or more, if you refused or failed the state's test(s). In addition, your commercial driving privilege in this state will be disqualified if you refused or failed the state's test(s). In addition, your commercial driving privilege in commercial motor vehicle. To appeal the driver's license revocation/disqualification, the Department of Public Safety must receive your written request for an Administrative Hearing within fifteen (15) days from the date you received notice of revocation or disqualification. Submit your written request in person at 3600 N. MLK Ave. or by mail to Legal Division, Department of Public Safety, P.O. Box 11415, OKC, OK, 73136

SERVING OF NOTICE: I personally hand-delivered a copy of the above containing the Notice of Revocation to the person named in Section 1 above.

Date served: 04/04/2019 Signature of officer:	Badge #:
In accordance with Title 12 O.S. Section 426, "I state under penalty of perjury	under the laws of Oklahoma that the foregoing is true and correct." Section 4
Date: 04/04/2019 Place (location when signed):	Signature of arresting officer:
NAME SMITH, JOSHUA	AGENCY BOARD OF TESTS

BADGE 0000000435	AGENCY ADDRESS INTER - AGENCY	PHONE "4054252460"

	i.	Certificate of Analysis – Wet Bath		
		09-26-2018	, <b>DZBG</b> <sup>0/13</sup>	
JA LE		Date of Analysis	Labelled target value (g/210L)	
JA AN		. 18020	0.0204	100
2010-	GUTH LABORATORIE		Average test result (c/210L)	
4. 1	590 NORTH 67th STREET @ HARRISBURG, PA 1711	-4511 @ TELEPHONE: 717-	564-5470	
$\sim$		BOT Technician Name and Signature	OMATLA	BRT Form
*		In accordance with RRT - 2.1.	Maintenance the shows referenced simu	

suitable for use as an external reference in maintenance of the Intoxilyzer 8000.

of

CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 18020

0.0238% (w/vol) ethyl alcohol. The expiration date for this lot

Alcohol Reference Solution for Simulator were analyzed by

chromatography on January 11, 2018, using a Perkin Elmer Gas

Chromatograph Autosystem XL S/N: 610N9030209, and found to contain

Intoxilyzer - Model 8000 09/26/2018	Alcohol Analyze	r SN 80-0025 15:
Test	g/210L	Tim
Air Blank	.0.000	15:3
Cal Check	0.020	15:3
Air Blank	0.000	15:3
Cal Check	0.020	15:3
Air Blank	0.000	15:3
Cal Check	0.020	15:35
Air Blank	0.000	15:41
Cal Check	0.020	15:41
Air Blank	0.000	15:41
Cal Check	0.020	15:42
Air Blank	0.000	15:42
Cal Check	0.020	15:43
Air Blank	0.000	15:44
Cal Check	0.021	15:44
Air Blank	0.000	15:45
Cal Check	0.021	15:46
Rir Blank	0.000	15:46
al Check	0.021	15:47
lir Blank	0.000	15:47
al Check	0.021	15:48
ir Blank	0.000	15:49
al Check Stats		
Average	0.0204	
Std Deu	0.0005	- V.
Rel Std Deu(%)	2.5312	

number is January 9, 2020 at 11:59 .PM. When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol

The alcohol and water used in this solution were free of test interfering substances.

analysis instrument reading of 0.020 g/210L +/- 3%.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN08031603 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.

Intoxilyzer Model 8000 09/26/2018	mendade formation			
Test	g/210L	Time		
Air Blank	.0.000	15:36		
Cal Check	0.020	15:37		
Air Blank	0.000	15:37		
Cal Check	0.020	15:38		
Air Blank	0.000	15:39		
Cal Check	0.020	15:39		
Air Blank	0.000	15:40		
Cal Check	0.020	15:41		
Air Blank	0.000	15:41		
Cal Check	0.020	15:42		
Air Blank	0.000	15:42		
Cal Check	0.020	15:43		
Air Blank	0.000	15:44		
Cal Check	0.021	15:44		
Air Blank	0.000	15:45		
Cal Check	0.021	15:46		
Rir Blank	0.000	15:46		
al Check	0.021	15:47		
Ma Dlask				

gas

Operato Signature

		Certificate o	f Analysis - Wet Bath	an tiller it.
∧ ®	"	09-210-18	,040	1
15 4 4		Date of Analysis	Labelled target value (g/210L)	
EXTASC	UTH LABORATOR	DEC 112410	0.0396	1001
- 590 I		17119 Water TELEPHONE	Average test result (g/210L)	
201 -000	ACKING THEET	Inch Sundh	111-001-0110	
		BOT Technician Name and Signatur	e	BRT Form 6.0
		In accordance with BRT - 2.1	.0, Maintenance, the above referenced sim	ulator solution is
	CERTIFICATE OF	A suitable for use as externa	I reference in maintenance of the Intoxilyzer	8000.

## Certified Alcohol Reference Solution for Simulator

Intoxilyzer - Model 8000 09/26/2018		SN 80-002591 15:54
1995 - 199		
Test	g/210L	Time
Air Blank	0.000	15:54
Cal Check	0.040	15:55
Air Blank	0.000	15:56
Cal Check	0,039	15:56
Air Blank	0.000	15:57
Cal Check	0.039	15:58
Air Blank	0.000	15:58
Cal Check	0.040	15:59
Air Blank	0.000	15:59
Cal Check	0.040	16:00
Air Blank	0.000	16:01
Cal Check	0.040	16:01
Air Blank	0.000	16:02
Cal Check	0.039	16:03
Air Blank	0.000	16:03
Cal Check	0.040	16:04
Air Blank	0.080	16:04
Cal Check	0.040	16:05
Air Blank	0.000	16:06
Cal Check	0.039	16:06
Air Blank	0.000	16:07
Cal Check Stats		
Average	0.0396	
Std Dev	0.0005	
Rel Std Deu(%)	1.3040	

Operator's Signature

RLPHR INSTRUMENT

Random Samples of Lot Number 17410 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on **December 7, 2017**, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.0481% (w/vol) ethyl alcohol. The expiration date for this lot number is **December 6, 2019** at 11:59 PM.

When used in a calibrated Simulator, operating at  $34^{\circ}C$  +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.040 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN12181501 whose values are traceable to NIST. All balances are calibrated annually by an outside agency using NIST traceable weights.

All balances are calibrated annually by an outside agency using NIST traceable weights Calibration verification is done prior to each use utilizing NIST traceable weights.

Certificate of Ana	Ilysis – Wet Bath	and the second state			
09-27-2018 ate of Analysis Labell 18070	·/08 edgeet value (g/210L) 34 9.9996 103/ ************************************	BORAT	ORIES.	INC.	
DT Technician Name and Signature	590 NORTH 67th STREE			1 • TELEPHONE	: 717-564-5470
in actordance with BRT – 2.1.0, Maini suitable for use as an external reference	tenance, the above referenced similation of the intenance of the intoxilyzer	ulator solution is 8000.			

Certified Alcohol Reference Solution for Simulator

Intoxilyzer - Model 8000 09/27/2018		SN 80-00259 08:28
Test	g/210L	Tine
Air Blank	0.000	08:28
Cal Check	0.100	08:29
Rir Blank	0.000	08:29
Cal Check	0.100	08:30
Air Blank	0,000	08:31
Cal Check	0.099	08:31
Air Blank	0.000	08:32
Cal Check	0.100	08:33
Air Blank	0.000	08:33
Cal Check	0.100	08:34
Air Blank	0.000	08:34
Cal Check	0.100	08:35
Air Blank	0.000	08:36
Cal Check	0.099	08:36
Air Blank	0.000	08:37
Cal Check	0.099	08:38
ir Blank	0.000	08:38
Cal Check	0.099	08:39
tir Blank	0.000	08:39
al Check	0.100	08:40
lir Blank	0.000	08:41
al Check Stats		
Average	0.0996	
Std Dev	0,0005	
Rel Std Deu(%)	0.5184	

Operator

Signatur

Random Samples of Lot Number 18070 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on February 28, 2018, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1219% (w/vol) ethyl alcohol. The expiration date for this lot number is February 26, 2020 at 11:59 PM.

When used in a calibrated Simulator, operating at  $34^{\circ}C$  +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability: Testing was conducted using Cerilliant Reference Standard lot number FN04271602 whose values are traceable to NIST. All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.

#### GUTH LABORATORIES, INC. 590 NORTH 67th STREET • HARRISBURG, PA 17111-4511 • TELEPHONE: 717-564-5470

In accordance with BRT – 2.1.0, *Maintenance*, the above referenced simulator solution is suitable for use as an external reference in maintenance of the Intoxilyzer 8000.

Labelled target value (g/210L)

test result (g/210L)

(R)

,200

0,1968

Certificate of Analysis - Wet Bath

09-27-2018

BOT Technician Name and Signature

Date of Analysis

Lot Number

17340

## CERTIFICATE OF ANALYSIS

BRT Form 6.0

Certified Alcohol Reference Solution for Simulator

ALPHA INSTRUMENT Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-002591 09/27/2018 15:27 Test g/210L Time Rir Blank 0.000 15:27 Cal Check 0.196 15:28 Air Blank 0.000 15:29 Cal Check 0.196 15:29 Air Blank 0.000 15:30 Cal Check 0.197 15:31 Air Blank 0.000 15:31 Cal Check 0.197 15:32 Air Blank 0.000 15:33 Cal Check 0.197 15:33 Air Blank 0.000 15:34 Cal Check 0.197 15:35 Air Blank 0.000 15:35 al Check 0.197 15:36 ir Blank 0.000 15:36 al Check 0.197 15:37 lir Blank 0.000 15:38 al Check 0.197 15:38 lir Blank 0.000 15:39 al Check 0.197 15:40 tin Blank 0.000 15:40 al Check Stats Average 0.1968 Std Deu 0.0004 Rel Std Deu(%) 0,2142

Operator

lignatur

Random Samples of Lot Number 17340 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on October 10, 2017, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.2410% (w/vol) ethyl alcohol. The expiration date for this lot number is October 9, 2019 at 11:59 PM.

When used in a calibrated Simulator, operating at  $34^{\circ}C$  +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.200 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN08101505 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.