

# STATE OF OKLAHOMA BOARD OF TESTS FOR ALCOHOL AND DRUG INFLUENCE

Dr. Jarrad Wagner, Ph.D. F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

#### ATTESTATION

I, the undersigned *Joshua Smith* (Director/Records Custodian), attest under oath pursuant to 12 O.S. §2902 (11)(a) this record was made, at or near the time of the occurrence of the matters set forth by or from information transmitted by a person having knowledge of those matters; was kept in the course of the regularly conducted business activity of the Board; and was made pursuant to the regularly conducted activity. The duplicate record provided is listed below:

1) The attached "Service Overview" is a true and accurate copy of maintenance performed on the Intoxilyzer 8000, serial number , in compliance with the agency's Maintenance Policy and Procedures.

This packet consists of pages.

Signed \_\_\_\_\_

Name of Position: Director/Records Custodian

Date of Attestation:

### **Service Overview**



3600 N. Martin Luther King Bldg #9 Oklahoma City, OK 73111 405-425-2460

BENCH CHECK DATE: 08-31-2020

SERIAL: 80-003462

### **DOCUMENTATION WITHIN SERVICE PACKET:**

- **BOT Certificate of Calibration and Operation**
- Bench Check
- Bench Check Technician's Report
- **BOT Dry Gas Certificate of Analysis**
- Mock Subject Test
- BOT Certificates of Analysis on Certified Solutions (4)



## STATE OF OKLAHOMA BOARD OF TESTS FOR ALCOHOL AND DRUG INFLUENCE

Administrative Office at:
Board of Tests for Alcohol and Drug Influence
Post Office Box 36307
Oklahoma City, Oklahoma 73136-2307

### **Certificate of Calibration and Operation**

This is to certify that the calibration of INTOXILYZER 8000, serial number 80-003462, manufactured by CMI Inc., a subsidiary of MPD, Inc. of Owensboro, Kentucky, was tested and found to be in compliance with National Highway Traffic Safety Administration Standard for Devices to Measure Breath Alcohol (F.R., No. 179 48705-48710 Sept. 17, 1993) and with International Recommendation OIML R 126 Edition 1998 (E). Calibration solutions are traceable to NIST standard material 1828.

In addition the operational procedures and functions of this Intoxilyzer have also been tested and found to be in compliance with the Oklahoma Board of Tests for Alcohol and Drug Influence - Rules and Regulations OAC Title 40: 25-1-2 and 30-1-3 and is hereby certified for use as an evidential breath testing device within this state.

CERTIFICATION TYPE	DATE
Manufacturer Calibration	10-07-08
Board of Tests Bench Check	10-20-08
Board of Tests Bench Check	04-27-10
Manufacturer Calibration	07-26-11
Board of Tests Bench Check	08-22-11
Board of Tests Bench Check	07-11-13
Board of Tests Bench Check	05-26-15
Board of Tests Bench Check	01-09-17
Board of Tests Bench Check	11-21-18
Board of Tests Bench Check	08-31-20

OFFICION TUDE



Presentation of this form certifies that the Intoxilyzer, SN: 80-003462, functioned properly at the time of the breath test and hereby certifies the breath test result as valid.

Date:

Re: Request for Log of Test on Intoxilyzer 8000

Name:

#### Dear Sir:

The Log of Test information from a breath test conducted on a Intoxilyzer 8000 is printed in Section 2 the Affidavit and is no longer retained as a separate form as with the Intoxilyzer 5000. The 8000 uses a nitrogen/ethanol dry gas for the external control test. The dry gas cylinder has the capacity to perform approximately 325 tests and does not require maintenance until its expiration date or can no longer provide the required pressure to perform a control test. The cylinder has an expiration date of two years from date filled at the manufacturer and is programmed in the 8000 upon installation, both dates are printed in Section 2 of the Affidavit. The Intoxilyzer is programmed to disable itself when the date of expiration is reached or the pressure monitor indicates 25 lbs psi.

Due to the reliability and stability of dry gas, maintenance on the 8000 is not required every 30 days or 25 tests, therefore maintenance is no longer performed by officers in the field. The Board of Tests for Alcohol and Drug Influence certifies each Intoxilyzer 8000 for field use by performing a bench check. The bench check is designed to check the operation and procedures programmed in the 8000 to ensure it is functioning properly for use in Oklahoma. A bench check is performed every time a new dry gas cylinder is installed or the instrument has returned from the manufacturer for repair. instrument passes the bench check, all that is required to set it up for field use is to plug it in and turn it on and does not require a permit to set up. The Board of Tests has the ability to monitor each instrument once it is set up in the field.

Each Intoxilyzer has a database of valid permit numbers that it checks prior to allowing a breath test to be performed. The database is updated on a regular basis, if a permit number is invalid that officer will not be able to conduct a breath test on the 8000. If the 8000 prints an affidavit the breath test operator listed in Section 2 had a valid permit at the time of the test and the breath test is valid.

If you have any further questions please feel free to call me.

J. Robert Blakeburn (OBA11059) State Director of Tests for Alcohol

7. pohr & Blakelun

and Drug Influence

# STATE OF OKLAHOMA - BOARD OF TESTS FOR ALCOHOL AND DRUG INFLUENCE CMI INTOXILYZER BENCH CHECK PROCEDURE REPORT

08/31/2020	108:33:3	INSTRUMENT LOCATION 80-003462				
nstrument type Intoxilyzer	INSTRUMENT MODEL I - 8000	INSTRUMENT SERIAL NUMBER 80-003462				

#### **TEST RESULTS**

**Diagnostics** 

PASS

**Mouth Alcohol Test** 

PASS

**RFI Test** 

PASS

**Abort Test** 

PASS

**Improper Sample Test** 

PASS

**Interferent Detect Test** 

PASS

No Sample Given Test

PASS

**Insufficient Sample Test** 

PASS

Air Blank Contamination Test

PASS

Range Exceeded Test

PASS

**Duplicate Breath Sampling Results** 

	noute Broati.		
Value	Test 1	Test 2	Control Test
.02	0.021	0.020	0.022
.04	0.040	0.039	0.041
.10	0.103	0.100	0.102
.20	0.203	0.199	0.201

PASS

### **CERTIFYING TECHNICIAN**

I swear under penalty of perjury that in regards to the above listed instrument, I have complied with the Intoxliyzer Bench Checking Procedure approved by the State of Oklahoma.

SIGNATURE	DATE
MM/1/2/2020	8/31/2020
PAPE, CHRISTOPHER	PERMIT # 000023020

3600 N. Martin Luther King Bldg #9 Oklahoma City, OK 73111 405-425-2460

### BENCH CHECK TECHNICIAN'S REPORT

Start Time: 08:30 Date: 08/31/2020 End Time: 09:30 INSTRUMENT Serial # Intoxilyzer Model 80-003462 8000 **GAS CANISTER** LOT# **EXP Date** 09/05/2022 20520080A1-024 REFERENCE Simulator Model Guth 2100 Guth 2100 Guth 2100 Guth 2100 Serial # DR2834 DR3700 DR2455 DR3470 Concentration .02 .04 .10 .20 LOT# 20190 20160 19050 20060 Manufactured Date 04/06/2020 03/18/2020 02/04/2019 02/10/2020 **Expiration Date** 03/18/2022 02/04/2021 02/10/2022 04/06/2022 Solution Commission Date 8/27/2020 8/27/2020 8/27/2020 8/27/2020 **REASON FOR BENCH CHECK** □ POST REPAIR □ CYLINDER REPLACEMENT ☐ TROUBLESHOOTING ☐ ROUTINE MAINTENANCE COMMENTS

Technician Signature / Permit #

Christopher Pape #23020

Printed Name

Certificate	of Analysis – Dry Gas	
8/21/2020	0.080	The state of the s
Date of Analysis	Labelled target value (g/210L)	
024	0.081	
Cylinder#	Average test result	
Christopher Pa	oc Wetter Noc 23	o je
BOT Technician Name and Signat	ure	specialty gases

In accordance with BOT Policy and Procedure Statement ERTS 30 Weightenings of Box 990 • Jacksonville, IL 62651-0790 referenced dry gas cylinder is suitable for use as an external reference with the intoxilyzer 8000.

217-245-2183 • Fax: 217-243-7634 • www.ilmoproducts.com

### Certificate of Analysis

Certificate ID:

13202 .

Part #:

BAC105L080T

Cylinder Size:

105L

Lot Number:

20520080A1

Expiration:

9/5/2022

0.080 BAC (For the calibration of instruments used to determine breath alcohol concentration)

Contents:

105 Liters @ 1000 psig 70°F (21°C)

Analytical

Reported

Accuracy

Analytical

Component:

Concentration:

(U, k=2):

Method:

Ethanol

208 ppm

+/-0.002 BAC(G/210L) NDIR

Nitrogen

Balance

[5.2 ppm]

Distributed by:

CMI Inc.

316 East Ninth Street Owensboro, KY 42303 Phone 866-835-0690

RCV1) 8/19/2020

www.alcoholtest.com

Operator's Signature

7/13

g/210L

0.000

0.080

0.000

0.081

0.000

0.081

0.000

0.081

0.000

0.080

0.000

0.081

0.000

0.081

0.000

0.081

0.000

0.081

0.000

0.081

0.000

0.0808

0.0004

0.5218

SN 80-002591 14:01

Time

14:01

14:02

14:02

14:03

14:03

14:04

14:04

14:04

14:05

14:05

14:06

14:06

14:07

14:07

14:07

14:08

14:08

14:09

14:09

14:09

14:10

ALPHA INSTRUMENT

Model 8000

08/21/2020

Test

Air Blank

Cal Check

Air Blank

Cal Check

Air Blank

Cal Check

Air Blank

Cal Check Air Blank

Cal Check

Air Blank

Cal Check

Air Blank

Cal Check

Air Blank

Cal Check

Air Blank

Cal Check

Air Blank

Cal Check

Air Blank

Average

Std Dev

Cal Check Stats

Rel Std Deu(%)

Intoxilyzer - Alcohol Analyzer

\*Traceable to:

Certified Reference Material - 262.4 µmol/mol

Ethanol in Nitrogen - Serial No. GN0015026 Lot No. 050319E11

Store in dry area, away from sources of heat, ignition and direct sunlight. Do not allow storage area to exceed 52 °C (125 °F).

Specialty Gas Lab Tech

08-10-2

Issuance Date

PJLA
Calibration and Testing
Accreditation #61895

The calibration results within this certificate were obtained using equipment and standards capable of producing analytical results traceable to NIST, and apply only to the items contained on this certificate. ILMO Products Company makes no warranty or representation as to the sultability of the use of any information provided for any particular purpose. The information use is at the sole discretion and risk of the user. Liability shall be limited to established replacement cost of this material or service.

### OFFICER'S AFFIDAVIT AND NOTICE OF DEVOCATION/DISQUALIFICATION

8/		3
O,	•	

	KE	OCA	HOND	ISQUA	LIF	10	AII	Or	4			Sec	tion
RREST DATE 08/31/2020		ST TIME		GOUNTY A	,	CITY	" #		CITATION	#			
RREST LOCATION	109	.00		123	CITY	1				COUN		7 10 L	
BOT											LAHO		
BJECT NAME NORTH, GAYI	ON LEVERNE						09/1	96	1 SEX	RACE	510	WEIGH 16	
DRESS	LUTHER KI	toward -	CITY	АНОМА С					STATE		731		
IVER LICENSE #	LOTHER KI	NG A	OKL	EXPIRATION D	ATE		STATE	CLA	Latin Secretarion	RESTR	RICTIONS	ENDORSE	MENT
083454871				06/30	/201		OK	A			Later	M	
HICLE MAKE 10CK		MOCK				TAG	OCK			OK.	N	PLACARD	
the above date, time, and lock hway, street, turnpike, public p uence of alcohol and/or other i scribe driving behavior or circu	arking lot, or other public plac ntoxicating substances as pro	e, or upon a private	of I had reasonable gr e road, street, alley or	ounds to believe the r lane which provides	person wa access to	one or	g or in actu more single	ai physic e or mult	i-family dw	elling within	this state w	a public roa hile under th	0
scribe person's condition (odor	actions atc.)					_		-		_			_
MOCK SUBJECT		E MAINTE	NANCE, IN	NSTALL NE	W CY	LINI	DER						
		THE BERSON	WAS INFORMED OF	THE IMPLIED CO.	NEWNIT TIP	er pre	MIEET						
ATH-ALCOHOL ANALYSIS	RECORD, REPORT OF THE							CORD.	-	_		Sec	tion
UIPMENT			NUMBER	SPECIALIST	0	n	mc						
-8000		80-	-003462	BOARD	OF !	res		_		-			
NDARD TY			0.080		230								
NUFACTURER			LOT #	UKUUK S Y	DATE CY	LINDER	RINSTALL	ED			ER EXPIRE		
LMO			20520	0080A1	08/	31,	/202	0	(	19/05	5/202	2.2	
905 08	B/31/2020	MOUTHPIECE G	period v	ute (or greater was observed p	prior to	and o	continui				THE S	TA	
'est	g/210L	Time	to	the time of te	st admi	nistra	tion.			103	***		
				and and m		A 1 -	- 1- 1 T			13/2	***/#	/* * /	1
Air Blank	0.00	09:22	The second of th	tered said E							00.00	2 867	B
Subject Test Breath Volume	0.00 1.417 LITERS	09:22		nce with the					5	H **	7	1/	F
Air Blank	0.00	09:22	proc	edure of the	Boar	d of	Tests			1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1		**	8
Wait		09:24	1/1/	X V	11-	- 1	רמהו	10		12/	1	*/	3/
Air Blank	0.00	09:25	and	1/1-19	1	0	200	0		1	***	**//	
Subject Test Breath Volume	0.00 1.390 LITERS	09:25	SIGNATURE C	OF OPERATOR			_				190		
Air Blank	0.00	09:25											
Cal Check	0.08	09:26	NAME	up ramonum	-						Test Re	esult	
Air Blank	0.00	09:26	PAPE, C	HRISTOPHE	PERMIT	,,		_	- 1	E	BrAC g	210L	
			0000000	0484	000	00:	2302	0			0 (	20	
			AGENCY	n mnama	-						0.0	00	
				F TESTS			1. 14. 10.			4	. (400)	Sac	tion
FECTIVE THIRTY (30) DAY ys or more, if you refused or cerating a CMV or CDL holde a result of this arrest by filing u may be eligible to partic vocation as a result of this quest IDAP within thirty (30 ERVING OF NOTICE: ccountability Program t	failed the state's test(s). In pr) OR a test result of .04 or g a petition in the District Co ipate in the Impaired Drivi arrest from appearing on 0) days of the date listed to I personally hand delification.	addition, your co r more while opera- burt of the County er Accountability your driving rec below. Mall the r vered a copy	mmercial driving pri aling a commercial of in which you were or Program (IDAP) a cord. Participation request for IDAP to of the above co	ivilege in this state motor vehicle. You arrested. Your petiadministered by the in IDAP may reduct the Department of the	will be disonally may appoint the must be departed to the architecture of Public architecture architecture to the architecture of Public architecture architectur	qualifie eal any be filed ment o mount Safety	d if you re Departme I within this I Public S of fees yo at P.O. Be	fused or ental act rty (30) Safety (I ou will b ox 1141	r failed the don agains days of th DPS). Co be require 5, Oklaho	e state's tes st your driv e date liste mpletion o d to pay to oma City, (	et (s) per license and below. of IDAP ma o the State OK 73136.	y prevent	a
ate served: 08/3		gnature of offic	17.0							_ Badg			
accordance with Title	12 O.S. Section 426, 1	"I state under	penalty of perj	ury under the la	aws of C	Oklaho	oma that	t the fo	regoing	is true a	and corre	ct." Sec	ction
	O Place (location w	vhen signed):_		Linguis	Si	gnatu	re of arr	esting	officer:				
ME PAPE, CHRI	STOPHER			BOARD	OF	TE	STS						
000000484	AGENCY ADDRESS INTER-AGE	NCY					ZIP CODE		PHONE 405	4252	460		
HER WITNESSES:					_		T.	-				Sec	ction
NAME			TITLE	ADDRESS					PHO	NE			-
N/T				1						. 1/2			
NAME			TITLE	ADDRESS					PHO	NE			
2													

#### 9/13

# OFFICER'S AFFIDAVIT AND NOTICE OF REVOCATION/DISQUALIFICATION

1/	1	Э	)
1	_	_	

ARREST DATE		REST TIME		COUNTY	1	CIT	Υ //		CITATIO	N#			
08/31/2020 RREST LOCATION	10	900		55	CITY	1				Icou	NTY		_
BOT					Oill						KLAH	AMC	1
JBJECT NAME	To a de la trada de la companya de l				DA	TE O	FBIRTH	-	SEX			T	WEIGHT
DDRESS	ON LEVERN		CITY		17.300	)8/	09/1	96	STATE	I W	ZIP CO	DE	165
3600 MARTIN	LUTHER K	ING A	OKLA	AHOMA C JEXPIRATION D	ATE		STATE	CLA	OK ss	RES	73 TRICTIONS	133 END	3 ORSEMEN
T083454871				06/30	/201		OK	A			Mro-Vi	M	
EHICLE MAKE MOCK		MOCE MOCE				TAC	OCK			OK	N	PLAC	CARD
n the above date, time, and loca ghway, street, tumplike, public lluence of alcohol and/or other i sescribe driving behavior or circu	arking lot, or other public p ntoxicating substances as p	ace, or upon a privat	e road, street, alley or la	ane which provide	access to c	one or	more single	or mult	i-family di	welling with	in this state	while u	nder the
escribe person's condition (odor	, actions, etc.):	- 1.4 F-10		A ! 1	10.0		0.07						
MOCK SUBJECT		THE PERSON	WAS INFORMED OF T	HE IMPLIED CO	NSENT TES	T RE	QUEST	CORD.					Section
QUIPMENT			NUMBER	SPECIALIST	LESSON TO		50 / Cor						
I-8000		80	-003462	BOARD	OF I								
TANDARD			0.080		230		ER						
Dry ANUFACTURER			LOT#				RINSTALL	ED	ĪĎ	ATE CYLIN	NDER EXPIR	RES	_
ILMO			205200	080A1	08/	31	/202	0			5/20		
	DATE	MOUTHPIECE		te (or greater	-								
	3/31/2020	G		as observed							_	-	
ESULT TABLE	W4556			he time of te				6		1	THE	STA	E
Test	g/210L	Time	101	ne mile of te	or delities	11.50		_	-	19	***	\**.	189
Air Blank Subject Test Breath Volume Air Blank Wait Air Blank Subject Test Breath Volume Air Blank	0.00 0.00 1.417 LITERS 0.00 0.00 0.00 1.390 LITERS 0.00	09:22 09:22 09:22 09:24 09:25 09:25	SIGNATURE OF	ce with the	applic	abi	le oper	ating	27	GREAT BEA		0	
Cal Check	0.08	09:26	NAME	e r amabar	185						Test R	Resul	t
Air Blank	0.00	09:26	PAPE, CH	RISTOPHE	PERMIT #				-		BrAC g	g/210	OL
			0000000	484			2302	0			0	00	1
			BOARD OF	mpama							U.	U	_
FECTIVE THIRTY (30) DAY ays or more, if you refused or perating a CMV or CDL holds a result of this arrest by filing ou may be eligible to partic vocation as a result of this quest IDAP within thirty (3)	failed the state's lest(s).  br) OR a test result of .04 g a petition in the District  pate in the Impaired Dr  arrest from appearing  days of the date liste  I personally hand de	In addition, your or or more while oper Court of the County iver Accountability on your driving red below. Mail the elivered a copy	ommercial driving privi ating a commercial my in which you were an y Program (IDAP) ad cord. Participation in request for IDAP to t of the above cont	lege in this state otor vehicle. You rested. Your pet ministered by the n IDAP may red he Department	will be disq i may appe- ition must b ne Departm ice the am of Public S	ualifical any be file nent count safety	ed if you rely Department of Public Store of fees you at P.O. Bo	fused or ental act rty (30) afety (I ou will b ox 1141	failed the	e state's t nst your dr he date ils ompletion ed to pay noma City of Depar	est(s) fiver license sted below. n of IDAP m to the Star n OK 73136	nay prote. Yo	u must
ERVING OF NOTICE: ccountability Program t	All the Charles of Annals and Annal State of the Annals and		non								lge #:		
ERVING OF NOTICE: ccountability Program t ate served: 08/3	1/2020	Signature of offi		ry under the l	aws of O	klah	oma that	the fo	regoin	g is true	and corr	ect."	Sectio
ERVING OF NOTICE:	1/2020 12 O.S. Section 426		penalty of perju	ry under the l			oma that ure of arr				and corr	ect."	Section
erving of Notice: ccountability Program to ate served: 08/3 accordance with Title ate: 08/31/202	1/2020 12 O.S. Section 426 0 Place (location	, "I state under	penalty of perju	AGENCY	Sig	gnati	ure of arr				and corr	ect."	Sectio
erving of notice: countability Program to the served: 08/3 accordance with Title tate: 08/31/202	1/2020 12 O.S. Section 426	, "I state under	penalty of perju		Sig	gnati		esting			and corr	ect."	Sectio

Certificate of Analysis – Wet Bath 12-10-19 Date of Analysis targa value (g/210L) 19050

020 ABORATORIES, INC.

In accordance with BRT - 2.2.0, Maintenance, the above referenced simulator solution is

suitable for use as an external reference in maintenance of the Intoxilyzer 8000.

### CERTIFICATE OF ANALYSIS

ALPHA INSTRUMENT Z Intoxilyzer - Alcohol Analyzer Model Bood 12/10/2019

**BOT Technician Name and Signature** 

Lot Number

14:52

SN 80-000923 Certified Alcohol Reference Solution for Simulator

Test g/210L Time Air Blank 0.000 14:53 Cal Check 0.019 14:53 Air Blank 0.000 14:54 Cal Check 0.020 14:55 Air Blank 0.000 14:55 Cal Check 0.019 14:56 Air Blank 0.000 14:56 Cal Check 0.020 14:57 Rir Blank 0.000 Cal Check 0.019 14:58 Air Blank 0.000 14:59 Cal Check 0.020 15:00 Air Blank 0.000 15:00 Cal Check 0.019 15:01 Air Blank 0.000 15:01 Cal Check 0.020 15:02 -Rir Blank 0.000 15:02 Cal Check 0.020 15:03 Air Blank 0,000 Cal Check 0.020 Air Blank 0.000 Cal Check Stats 0.0196 Ruerage

Random Samples of Lot Number 19050 Alcohol Reference Solution for Simulator were analyzed by chromatography on February 5, 2019, using a Perkin Elmer Gas gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.0236% (w/vol) ethyl alcohol. The expiration date for this lot number is February 4, 2021 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.020 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

perator's Signature

0.0005

Std Deu

Rel Std Deu(%) 2,6345

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN08031603 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.

Certificate o	f Analysis – Wet Bath	
7/22/2020 Date of Analysis	0.040 Labelled targetypalue (g/210L)	
20060 Lot Number	Acera le talt result (g/23-01.)	
Christopher Pap	the same of the sa	HABORATORIES, INC.
BOT Technician Name and Signatur	.0, Maintenance, the above reference	BRT Form 6.0

suitable for use as an external reference in maintenance of the intoxilyzer 8000.

SN 80-003400

14:57

Tine

14:57

14:58

14:59

14:59

15:00

15:01

15:01

15:02

15:02

15:03

15:04

15:04

15:05

15:06

15:06

15:07

15:07

15:08

15:09

15:09

15:10

CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 20060 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on February 11, 2020, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.0483% (w/vol) ethyl alcohol. The expiration date for this lot number is February 10, 2022 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.040 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Operator's Signature

ALPHA UNIT SUB

Model 8000

07/22/2020

Test

Air Blank

Cal Check

Air Blank

Average Std Dev

Cal Check Stats

Rel Std Deu(%)

Intoxilyzer - Alcohol Analyzer

g/210L

0.000

0.040

0.000

0.039

0.000

0.039

0.000

0.039

0.000

0.039

0.000

0.039

0.000

0.040

0.000

0.040

0.000

0.040

0.000

0,040

0.000

0.0395

0.0005

1.3342

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN12181501 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.

Certificate of Analysis – Wet Bath						
06/16/	20	,100				
Date of Analysis	-	Labelled target value (g/210L)				
20190		.098				
Lot Number		August test result (g/210L)				
1						

16:18

LABORATORIES, INC.

WEET BRE FORM 5.0 RISE JRG. PA 17111- 4511 9 TELEPHONE: 717-564-5470

In accordance with BRT - 2,2.0, Maintenance, the above referenced simulator solution is suitable for use as an external reference in maintenance of the intoxilyzer 8000.

### CERTIFICATE OF ANALYSIS

ALPHA INSTRUMENT Intoxilyzer - Alcohol Analyzer Certified Alcohol Reference Solution for Simulator Model 8000 SN 80-002591

9/210L Time Air Blank 0.000 16:18 Cal Check 0.099 16:19 Air Blank 0.000 16:20 Cal Check 0.098 16:20 Air Blank 0.000 16:21 Cal Check 0.097 16:22 Air Blank 0.000 16:22 Cal Check 0.098 16:23 Air Blank 0.000 16:24 Cal Check 0.098 16:24 Air Blank 0.000 16:25 Cal Check 0.098 16:26 Air Blank 0.000 16:26 Cal Check 0.098 16:27 air Blank 0.000 16:27 Cal Check 0.097 16:28 Rir Blank 0.000 16:29 tal Check 0.098 16:29 lir Blank 0.000 16:30 :al Check 0.098 16:31

0.000

0.0979

0.0005

0.5798

16:31

06/16/2020

Test

lir Blank

Std Dev

al Check Stats Average

Rel Std Deu(%)

Random Samples of Lot Number 20190 Alcohol Reference Solution for Simulator were analyzed by gas chromatography on April 8, 2020, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1206% (w/vol) ethyl alcohol. The expiration date for this lot number is April 6, 2022 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Operator Signature

> Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN02271802 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.



In accordance with BRT - 2.2.0, Maintenance, the above referenced simulator solution is suitable for use as an external reference in maintenance of the Intoxilyzer 8000.

### CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 20160 of
Alcohol Reference Solution for Simulator were analyzed by
gas chromatography on March 20, 2020, using a Perkin Elmer Gas
Chromatograph Autosystem XL S/N: 610N9030209, and found to contain
0.2427% (w/vol) ethyl alcohol/ The expiration date for this lot
number is March 18, 2022 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.200 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Operator's Standure

LPHR UNIT SUB htoxilyzer - Alcohol Analyzer

g/210L

0.196

0.000

0.196

0.000

0.197

0.000

0.196

0.000

0.198

0.000

0.198

0.000

0.197

0.000

0.198

0.000

0.198

0.000

0.197

0.000

0.1971

0.0008

7/07/2020

ir Blank

al Check

Hr Blank

al Check

ir Blank al Check Stats

Average

Std Deu

Rel Std Deu(%) 0.4442

SN 80-003400

15:25

Time

15:26

15:26

15:27

15:28

15:29

15:30

15:31

15:31

15:32

15:33

15:33

15:34

15:35

15:35

15:36 15:37

15:37

15:38

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN08101505 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.