

Dr. Jarrad Wagner, Ph.D. F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

#### ATTESTATION

I, the undersigned *Joshua Smith* (Director/Records Custodian), attest under oath pursuant to 12 O.S. §2902 (11)(a) this record was made, at or near the time of the occurrence of the matters set forth by or from information transmitted by a person having knowledge of those matters; was kept in the course of the regularly conducted business activity of the Board; and was made pursuant to the regularly conducted activity. The duplicate record provided is listed below:

1) The attached "Service Overview" is a true and accurate copy of maintenance performed on the Intoxilyzer 8000, serial number 80-004896, in compliance with the agency's Maintenance Policy and Procedures.

Signed

Name of Position: Director/Records Custodian

Date of Attestation: 11/08/2023

# **Service Overview**

3600 N. Martin Luther King Bldg #9 Oklahoma City, OK 73111 405-425-2460

BENCH CHECK DATE: 11/8/2023

SERIAL: 80-004896

#### **DOCUMENTATION WITHIN SERVICE PACKET:**

- BOT Certificate of Maintenance, Calibration and Operation
- Bench Check
- BOT Dry Gas Certificate of Analysis
- · Ilmo Specialty Gases Certificate of Analysis
- Mock Subject Test
- BOT Certificate of Analysis on Certified Solutions (4)
- Guth Laboratories Certificate of Analysis (4)



Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

# CERTIFICATE of MAINTENANCE, CALIBRATION, and OPERATION

This certifies the Intoxilyzer 8000 referenced below, manufactured by CMI, Inc., a subsidiary of MPD, Inc. of Owensboro, Kentucky, is maintained by the Board of Tests for Alcohol and Drug Influence in accordance with 40 O.A.C. 30-1-3. The instrument reported below is hereby deemed to be in compliance with Board rules via linearity accuracy examination and operational/functional procedure evaluation also referred to as a "Bench Check". The Intoxilyzer 8000 appears on the National Highway Traffic Safety Administration's Conforming Products List of Evidential Breath Alcohol Measurement Devices (82 FR 50940 et seg) as an approved device. In accordance with 40 O.A.C. 25-1-2, this instrument is approved and certified for use in the State of Oklahoma as an automated evidential breath-alcohol measurement device by the Board of Tests pursuant to 47 O.S. § 759.

All analyses conducted by this instrument are considered precise, valid, and in accordance with 40 O.A.C. 30-1-3(b) Operating Procedure(s). A printed breath test operator affidavit containing the seal of the Board by this instrument confirms a calibration check successfully verified the instrument's calibration at the time of test. The printed breath test affidavit also serves as proof that safeguards in place to protect the breath samples, including but not limited to Ambient Air Interference, Radio Frequency Interference, Residual Mouth Alcohol, Deficient Samples, etc. were not triggered and the test is a true and accurate scientific measurement of the person's alcohol concentration.

- Maintenance was performed on the Intoxilyzer 8000, Serial No. 80-004896 on 11/8/2023 in accordance with 40 O.A.C. 30-1-3 and internal BOT maintenance policy and procedures.
- 2. A dry gas canister appearing on the National Highway Traffic Safety Administration's Conforming Products List of Calibrating Units for Breath Alcohol Testers (77 FR 64588 et seg) and approved in accordance with 40 O.A.C. 25-1-3 and BOT maintenance policy was installed.
  - Manufacturer: ILMO Unit: 208 ppm / 105 L
  - Cylinder Lot No. Canister No.: 14323080A4-038
  - BOT analysis average test result: 0.081
  - Cylinder expiration date: 6/5/2025

3. In accordance with 40 O.A.C. 25-1-2.1, the approved reference standards and equipment were used:

Simulator Model		In the second second second	and the second
Guth 2100	Guth 2100	Guth 2100	Guth 2100
Serial #			
DR2834	DR3694	DR2455	DR3470
Concentration			
.02	.04	.10	.20
LOT#			
23060	22460	23180	23340
Manufactured Date			
2/28/2023	12/28/2022	5/17/2023	9/18/2023
Expiration Date			LT Takenda Ja
2/28/2025	12/28/2024	5/17/2025	9/18/2025
Solution Commission Date		Academie	
11/6/2023	11/6/2023	11/6/2023	11/8/2023

	2/28/2025	12/28/2024	5/17/2025	9/18/2
Solutio	n Commission Date		Andrews	
	11/6/2023	11/6/2023	11/6/2023	11/8/2
LUN VALLE	A. A. L. A. A.	**		
Completed by BC	F Personnel:			
1/1/8/				
4/1/2/1/20	Christopher Pape	Breath-Alcohol Testi	ng Program Administrator	
Signature/Name	The state of the s	Title	CARL CONTRACTOR CONTRACTOR	

Certified

Joshua Smith State Director Signature

H CHECK DATE L / 08/2023	11:40:5		STRUMENT LOCATION	6	
ntoxilyzer	INSTRUMEN			80-004	
itoxilyzei	1-80	TEST RE	SULTS	80-00	1090
Diagnosti	cs				
PASS					
Mouth Alcoho	ol Test				
RFI Tes	t				
PASS		Duj	olicate Breath	Sampling Re	sults
Abort Te PASS		Value	Test 1	Test 2	Control Test
Improper Samp PASS		.02	0.018	0.018	0.019
Interferent Det	ect Test	.04	0.038	0.037	0.039
PASS		.10	0.095	0.094	0.098
No Sample Giv PASS		.20	0.192	0.190	0.196
Insufficient Sam			PA	ss	
PASS	ipie rest				
Air Blank Contami PASS	nation Test				
Range Exceed PASS					
swear under penalt he Intoxliyzer Benc			above listed	The state of the s	and the second of the second of the second of the second
URE THO XITY ZET BETTE	in Checking I rocce	aure approved	by the State C	DATE	1
Charles	Per			15	18/20
PE, CHRISTOR	PHER			PERMIT#	023020



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

### Certificate of Analysis - Dry Gas

ALPHA INSTRUMENT

Intoxilyzer - Alcohol Analyzer

Model 8000

SN 80-002591

09/26/2023

08:44

Test	g/210L	Time
Air Blank	0.000	08:44
Cal Check	0.081	08:44
Air Blank	0.000	08:45
Cal Check	0.081	08:45
Air Blank	0.000	08:46
Cal Check	0.082	08:46
Air Blank	0.000	08:47
Cal Check	0.082	08:47
Air Blank	0.000	08:47
Cal Check	0.081	08:48
Air Blank	0.000	08:48
Cal Check	0.082	08:49
Air Blank	0.000	08:49
Cal Check	0.082	08:49
Air Blank	0.000	08:50
Cal Check	0.081	08:50
Air Blank	0.000	08:51
Cal Check	0.081	08:51
Air Blank	0.000	08:52
Cal Check	0.081	08:52
Air Blank	0.000	08:52
Cal Check Stats		
Average	0.0814	
Std Dev	0.0005	
Rel Std Dev(%)	0.6343	

Known Breat	h-Alcohol Equivalent (g/210L) Ratio of Nitrogen and 0.08 g/210L = 208 PPM ÷ 2605*	
Manufacturer: ILMO	Manufacturer Cert. No.: 15762	Cylinder Size: 105L
Lot No.: 14323080A4	Cylinder No.: 038	Expiration Date: 06/05/2025
prescribed by the State Director of Tests p	ursuant to 40 OAC 30-1-3. The above NIST traceable 759, 40 OAC 25-1-3, and appears on the current confo of Transportation in the Federal Register [77 FR 64588	not only was performed in accordance with the procedures reference standard meets the requirements for use as a corming products list of calibrating units for breath alcohol et seq]. * See Dubowski, K.M. and Essary N., Journal of
Signature And Advantage of the Signature	Date: 09/26	/2023
Name Christopher Pape	Permit No.: 23	020





7 Eastgate Dr. • P.O. Box 790 • Jacksonville, IL 62651-0790 217-245-2183 • Fax: 217-243-7634 • www.ilmoproducts.com

# Certificate of Analysis

Certificate ID:

15762

Part #:

BAC105L080T

Cylinder Size:

105L

Lot Number:

14323080A4

Expiration:

6/5/2025

9/18/2023

0.080 BAC (For the calibration of instruments used to determine breath alcohol concentration)

Contents:

105 Liters @ 1000 psig 70°F (21°C)

Analytical

Reported

Analytical

Component: 😤

Accuracy

Concentration:

(U, k=2):

Method:

Ethanol Nitrogen 208 ppm Balance +/-0.002 BAC(G/210L) NDIR [5.2 ppm]

Distributed by:

CMI Inc.

316 East Ninth Street Owensboro, KY 42303 Phone 866-835-0690 www.alcoholtest.com

\*Traceable to:

Certified Reference Material - 261.0 µmol/mol

Ethanol in Nitrogen - Serial No. ND7017 Lot No. 080722E2

Store in dry area, away from sources of heat, ignition and direct sunlight. Do not allow storage area to exceed 52 °C (125 °F).

06-01-2023

Issuance Date

Accreditation #61895

The calibration results within this certificate were obtained at the facility listed above using equipment and standards capable of producing analytical results traceable to NIST, and apply only to the items contained on this certificate. ILMO Products Company makes no warranty or representation as to the suitability of the use of any information provided for any particular purpose. The information use is at the sole discretion and risk of the user. Liability shall be limited to established replacement cost of this material or service. This certificate applies only to the items described and shall not be reproduced other than in full, without written approval from the issuing facility.

#### 7/17

# OFFICER'S IMPAIRED DRIVING AFFIDAVIT

			10000000000000					Section
11/08/2023		1300	55	CITY#	CITATION	#		
ARREST LOCATION BOT			Cit	TY.		OK	LAHC	MA
BREATHTEST,	FRANK			03/15/			509	9 175
3600 MARTIN	LUTHER	KING	OKLAHOMA CIT	Y	OK	T-	73:	
H083463197			EXPIRATION DATE 10/31/2	022 OK	D	RESTR	RICTIONS	ENDORSEMENTS
MOCK		MOCK		MOCK		OK	N CMV	N N

On the above date, time, and location, the above named person was arrested, and I had reasonable grounds to believe the person was driving or in actual physical control of a motor vehicle upon a public road, highway, street, turnpike, public parking lot, or other public place, or upon a private road, street, alley or lane which provides access to one or more single or multi-family dwelling within this state while under the influence of alcohol and/or other intoxicating substances as prohibited by law. (Describe driving behavior or circumstances):

MOCK SUBJECT TEST

Describe person's condition (odor, actions, etc.):

POST REPAIR MAINTENANCE, INSTALLED NEW CYLINDER

	THE PERSON V	VAS INFORMED OF	THE IMPLIED C	ONSENT TEST REQUEST	
BREATH-ALCOHOL ANALYSIS RECORD, REPORT OF THE PI	ERSON NAMED	N SECTION 1 ABOV	E, and LOG OF	TEST AND MAINTENANCE RECORD.	Section 2
EQUIPMENT I-8000		NUMBER 004896	BOARD	OF TESTS	
STANDARD Dry		0.080		PERMIT NUMBER 23020	
MANUFÁCTURER ILMO		LOT#	080A4	DATE CYLINDER INSTALLED 11/08/2023	DATE CYLINDER EXPIRES 06/05/2025
OBSERVATION START TEST DATE 11/08/2023				er) deprivation/observation I prior to and continuing	T-Walter -
RESULT TABLE  Test g/210L	Time			est administration.	OF THE STATE
Air Blank 0.00 Subject Test 0.00 Breath Volume 2.074 LITERS Air Blank 0.00 Wait Air Blank 0.00 Subject Test 0.00 Breath Volume 2.093 LITERS Air Blank 0.00 Cal Check 0.08 Air Blank 0.00	13:35		CE With the edure of the properator	Test Result BrAC g/210L 0.00	

#### What happens next?

Section 3

- The Arresting Officer is required to submit the Officer's Impaired Driving Affidavit to the State of Oklahoma.
- 2. Upon receipt and review of the Officer's Impaired Driving Affidavit and any associated test results, if applicable, the State of Oklahoma may revoke your Oklahoma driving privileges. If you hold a commercial driver license or were driving a commercial motor vehicle, the State of Oklahoma may disqualify your privilege to operate a commercial motor vehicle. You will continue to have driving privileges, if otherwise eligible, unless or until an Order of Revocation and/or Disqualification is issued by the State of Oklahoma.
- The State of Oklahoma will send you the Order of Revocation and/or Disqualification to the address you have on file with Service Oklahoma. To confirm or change your address on file, contact Service Oklahoma.
- The Revocation and/or Disqualification will be effective 55 days from the date of the Order, including 10 days allowed for mailing, 47 O.S. §§2-116, 753, 754.

  The Order of Revocation and/or Disqualification will contain important information about your rights to appeal the revocation and/or disqualification of your driving privileges. To appeal the applicable Order of Revocation and/or Disqualification your petition must be filed with the District Court of the County in which you were arrested within 40 days of the date of the order, including 10 days allowed for mailing. 47 O.S. §§2-116, 6-211.

		12 O.S. Section 426, "I state under p	enalty of perj				Section +
Date		Place (location when signed):			arresting officer:		
PI		STOPHER		BOARD OF T	ESTS		
O C	000000484	AGENCY ADDRESS INTER-AGENCY			ZIP CODE 00000	4054252460	
ОТН	ER WITNESSES						Section :
1	NAME		TITLE	ADDRESS		PHONE	
2	NAME		TITLE	ADDRESS		PHONE	

#### 8/17

# OFFICER'S IMPAIRED DRIVING AFFIDAVIT

												Section
11/08/2023		1300			55	(	CITY#	Ci	TATION#	The state of the s		
ARREST LOCATION BOT						CITY	0.10			OK	y LAHC	MA
BREATHTEST,	FRANK						OF BIRTH	1981	M	Y I	505	
3600 MARTIN	LUTHER	KING		OKLA	HOMA CI	TY			OK		ZIP COL	
H083463197					10/31/		STATE	D		RESTRI	CTIONS	ENDORSEMENTS
VEHICLE MAKE MOCK		11.51	MOCK				MOCK	-	1000	OK	N N	N N

On the above date, time, and location, the above named person was arrested, and I had reasonable grounds to believe the person was driving or in actual physical control of a motor vehicle upon a public road, highway, street, tumpike, public parking lot, or other public place, or upon a private road, street, alley or lane which provides access to one or more single or multi-family dwelling within this state while under the influence of alcohol and/or other intoxicating substances as prohibited by law. (Describe driving behavior or circumstances):

MOCK SUBJECT TEST

Describe person's condition (odor, actions, etc.):

POST REPAIR MAINTENANCE, INSTALLED NEW CYLINDER

		THE PERSON V	VAS INFORMED OF	THE IMPLIED CO	DISENT TE	ST REQUEST	
BREATH-ALCOHOL ANALYSIS	RECORD, REPORT OF THE P	ERSON NAMED	N SECTION 1 ABOV	/E, and LOG OF T	EST AND M	IAINTENANCE RECORD.	Section 2
I-8000			NUMBER 004896	SPECIALIST BOARD	OF T	rests	
STANDARD Dry			0.080		PERMIT I		
MANUFACTURER ILMO			14323	080A4	DATE CY	08/2023	DATE CYLINDER EXPIRES 06/05/2025
OBSERVATION START TEST	DATE L/08/2023				r) depriv	ation/observation	1.11.04.7
RESULTTABLE  Test  Air Blank Subject Test  Breath Volume Air Blank Wait Air Blank Subject Test	g/210L 0.00 0.00 2.074 LITERS 0.00 0.00	13:30 13:31 13:31 13:33 13:33 13:33	I administ	the time of to tered said I ace with the edure of the	est admi Breath- e applie	Alcohol Test in cable operating of Tests.	STATE STATE OF THE
Breath Volume Air Blank Cal Check Air Blank	2.093 LITERS 0.00 0.08 0.00	13:34 13:35 13:35	NAME PAPE, CHI BADGE # 000000	RISTOPHE	000	0023020	Test Result BrAC g/210L 0.00

#### What happens next?

Section 3

- 1. The Arresting Officer is required to submit the Officer's Impaired Driving Affidavit to the State of Oklahoma.
- Upon receipt and review of the Officer's Impaired Driving Affidavit and any associated test results, if applicable, the State of Oklahoma may revoke your Oklahoma driving privileges. If you hold
  a commercial driver license or were driving a commercial motor vehicle, the State of Oklahoma may disqualify your privilege to operate a commercial motor vehicle. You will continue to have
  driving privileges, if otherwise eligible, unless or until an Order of Revocation and/or Disqualification is issued by the State of Oklahoma.
- driving privileges, if otherwise eligible, unless or until an Order of Revocation and/or Disqualification is issued by the State of Oklahoma.

  3. The State of Oklahoma will send you the Order of Revocation and/or Disqualification to the address you have on file with Service Oklahoma. To confirm or change your address on file, contact Service Oklahoma.

. The Revocation and/or Disqualification will be effective 55 days from the date of the Order, including 10 days allowed for mailing, 47 O.S. §§2-116, 753, 754.

The Order of Revocation and/or Disqualification will contain important information about your rights to appeal the revocation and/or disqualification of your driving privileges. To appeal the applicable Order of Revocation and/or Disqualification your petition must be filed with the District Court of the County in which you were arrested within 40 days of the date of the order, including 10 days allowed for mailing. 47 O.S. §§2-116, 6-211.

In accordance with Title 12 O.S. Section 426, "I state under penalty of perjury under the laws of Oklahoma that the foregoing is true and correct."						
Date:	Place (location when signed):	Signature of arresting officer:				
PAPE, CH	RISTOPHER	BOARD OF TESTS				
000000048	AGENCY ADDRESS INTER-AGENCY	ZIP CODE				



# STATE OF OKLAHOMA

Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

BOARD OF TESTS FOR ALCOHOL AND DRUG INFLUENCE
Post Office Box 36307
Oklahoma City, Oklahoma 73136-2307
Phone: (405) 425-2460 Fax: (405) 425-2490
www.bot.ok.gov

Joshua Smith Director

### Certificate of Analysis - Reference Solution

ALPHA INSTRUMENT Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-002591 09/06/2023 10:21

Test	g/210L	Time
Air Blank	0.000	10:22
Cal Check	0.019	10:23
Air Blank	0.000	10:23
Cal Check	0.019	10:24
Air Blank	0.000	10:24
Cal Check	0.019	10:25
Air Blank	0.000	10:26
Cal Check	0.019	10:26
Air Blank	0.000	10:27
Cal Check	0.020	10:28
Air Blank	0.000	10:28
Cal Check	0.019	10:29
Air Blank	0.000	10:29
Cal Check	0.019	10:30
Air Blank	0.000	10:31
Cal Check	0.019	10:31
Air Blank	0.000	10:32
Cal Check	0.018	10:33
Air Blank	0.000	10:33
Cal Check	0.019	10:34
Air Blank	0.000	10:34
Cal Check Stats		
Average	0.0190	
Std Dev	0.0004	
Rel Std Dev(%)	2.4809	

NIST Traceable Breath-Alcohol Ed	quivalent (g/210L +/- 3%) Ratio of ethanol and w	ater free of test interfering substances.
Manufacturer: Guth	Manufactured Date: 2/28/2023	Concentration: 0.020
Lot No.: 23060	Expiration Date: 2/28/2025	Equipment Model: Guth 2100
prescribed by the State Director of Tests pursuant to 40 OAC 25-1-2.1. Approved breath-alcohol mea Conforming Products List of Calibrating Units on E	asurement equipment and reference solutions. Equ Breath Alcohol Testers (72 FR 34747), published by	as performed in accordance with the procedures brence solution meets the requirements for use pursuant ipment used appears on the current or supplemented the National Highway Traffic Safety Administration.
Signature Alleman Alleman	Date: 9/6/2023	3
Name Christopher Pape	Permit No.: 2302	0



CP RCVP 8/23/23

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 23060 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on March 1, 2023, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.0241% (w/vol) ethyl alcohol. The expiration date for this lot number is February 28, 2025 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.020 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN05122004 whose values are traceable to NIST.



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

### Certificate of Analysis - Reference Solution

ALPHA INSTRUMENT
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-002591
09/06/2023 11:15

Test	g/210L	Time
Air Blank	0.000	11:16
Cal Check	0.040	11:16
Air Blank	0.000	11:17
Cal Check	0.040	11:18
Air Blank	0.000	11:18
Cal Check	0.040	11:19
Air Blank	0.000	11:19
Cal Check	0.040	11:20
Air Blank	0.000	11:21
Cal Check	0.040	11:21
Air Blank	0.000	11:22
Cal Check	0.040	11:23
Air Blank	0.000	11:23
Cal Check	0.040	11:24
Air Blank	0.000	11:24
Cal Check	0.040	11:25
Air Blank	0.000	11:26
Cal Check	0.040	11:26
Air Blank	0.000	11:27
Cal Check	0.040	11:28
Air Blank	0.000	11:28
Cal Check Stats		
Average	0.0400	
Std Dev	0.0000	,
Rel Std Dev(%)	0.0000	

Manufacturer: Guth	Manufactured Date: 12/28/2022	Concentration: 0.040
Lot No.: 22460	Expiration Date: 12/28/2024	Equipment Model: Guth 2100
to 40 OAC 25-1-2.1. Approved breath-alcoho	of measurement equipment and reference solutions. Eq	erence solution meets the requirements for use pursuan
Conforming Products List of Calibrating Units Signature	s for Breath Alcohol Testers (72 FR 34747), published b	y the National Highway Traffic Safety Administration.



CP RCVD 8/31/23

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 22460 of
Alcohol Reference Solution for Simulator were analyzed by
gas chromatography on December 29, 2022, using a Perkin Elmer Gas
Chromatograph Autosystem XL S/N: 610N9030209, and found to contain
0.0486% (w/vol) ethyl alcohol. The expiration date for this lot
number is December 28, 2024 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.040 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN06032102 whose values are traceable to NIST.



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

## Certificate of Analysis - Reference Solution

ALPHA INSTRUMENT
Intoxilyzer - Alcohol Analyzer
Model 8000 SN 80-002591
09/06/2023 11:36

Test	g/210L	Time
Air Blank	0.000	11:37
Cal Check	0.099	11:37
Air Blank	0.000	11:38
Cal Check	0:100	11:39
Air Blank	0.000	11:39
Cal Check	0.100	11:40
Air Blank	0.000	11:40
Cal Check	0.100	11:41
Air Blank	0.000	11:42
Cal Check	0.101	11:42
Air Blank	0.000	11:43
Cal Check	0.100	11:44
Air Blank	0.000	11:44
Cal Check	0.100	11:45
Air Blank	0.000	11:45
Cal Check	0.100	11:46
Air Blank	0.000	11:47
Cal Check	0.100	11:47
Air Blank	0.000	11:48
Cal Check	0.100	11:49
Air Blank	0.000	11:49
Cal Check Stats		
Average	0.1000	
Std Dev	0.0004	
Rel Std Dev(%)	0.4714	

Manufacturer: Guth	Manufactured Date: 5/17/2023	Concentration: 0.100	
Lot No.: 23180	Expiration Date: 5/17/2025	Equipment Model: Guth 2100	
prescribed by the State Director of Tests to 40 OAC 25-1-2.1, Approved breath-ake	cohol measurement equipment and reference solutions.	reference solution meets the requirements for use pursuant Equipment used appears on the current or supplemented d by the National Highway Traffic Safety Administration.	
Signature	Date: 9/6/20	9/6/2023	
Name Christopher Pape		Permit No.: 23020	



CP RCVD 8/23/23

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 23180 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on May 18, 2023, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1220% (w/vol) ethyl alcohol. The expiration date for this lot number is May 17, 2025 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN11172002 whose values are traceable to NIST.



Dr. Jarrad Wagner, Ph.D., F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

### Certificate of Analysis - Reference Solution

ALPHA INSTRUMENT

Intoxilyzer - Alcohol Analyzer

Model 8000

SN 80-002591

1	0/	04	/20	23
---	----	----	-----	----

09:22

Test	g/210L	Time
Air Blank	0.000	09:22
Cal Check	0.198	09:23
Air Blank	0.000	09:24
Cal Check	0.198	09:24
Air Blank	0.000	09:25
Cal Check	0.199	09:26
Air Blank	0.000	09:26
Cal Check	0.198	09:27
Air Blank	0:000	09:27
Cal Check	0.198	09:28
Air Blank	0.000	09:29
Cal Check	0.198	09:29
Air Blank	0.000	09:30
Cal Check	0.198	09:31
Air Blank	0.000	09:31
Cal Check	0.198	09:32
Air Blank	0.000	09:33
Cal Check	0.199	09:33
Air Blank	0.000	09:34
Cal Check	0.199	09:35
Air Blank	0.000	09:35
Cal Check Stats		
Average	0.1983	
Std Dev	0.0004	
Rel Std Dev(%)	0.2435	

Manufacturer: Guth	Manufactured Date: 9/18/202		
Lot No.: 23340	Expiration Date: 9/18/2025	Equipment Model: Guth 2100	
prescribed by the State Director of Tests to 40 OAC 25-1-2.1. Approved breath-alc	pursuant to 40 OAC 30-1-3. The above NIST to ohol measurement equipment and reference so	8000 and was performed in accordance with the procedures aceable reference solution meets the requirements for use pursuan flutions. Equipment used appears on the current or supplemented published by the National Highway Traffic Safety Administration.	
Signature / Signature		Date: 10/4/2023	
Name Christopher Pape		o.: 23020	



RCVD 10/2/2023

Certified Alcohol Reference Solution for Simulator

P

Random Samples of Lot Number 23340 of
Alcohol Reference Solution for Simulator were analyzed by
gas chromatography on September 19, 2023, using a Perkin Elmer Gas
Chromatograph Autosystem XL S/N: 610N9030209, and found to contain
0.2399% (w/vol) ethyl alcohol. The expiration date for this lot
number is September 18, 2025 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.200 g/210L+/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN09142002 whose values are traceable to NIST.

Location Returned From Stillwater PD 2	<u>Unit Serial #</u> 80-004896	17/17 <u>Date:</u> 9/11/23
nventory Process  Unit checked in to BOT  Unlink Gas Canister fro	Old Cylinder #:	ventory
Check in Procedure  Connect unit to ethernet connection 172.019.004.  Download any stored records  Change unit location to Serial number  Update Unit Date and time  Previous time: $0/9/23$	Notes from the Field	ssues in the Field
Lab Evaluation  Breath Hose Hold Vacuum?:	Red (F7)  Receipt paper  Flow Sensor	Is regulator close to Tank Sensor Value?
Sent to CMI (If needed) Reason Sent  Date Sent Invoice Sent Date Returned		
Replace Break hose	Repairs Made:  Breath Lose	replaced
Inventory Process ☐ Bench Check Ted ☐ Linked Cylinder to Unit ☐ Mock Subject Te ☐ Changed Unit Status ☐ Certificate of Cal	st Completed	if unit sent to CMI der ate of Calibration