

## STATE OF OKLAHOMA BOARD OF TESTS FOR ALCOHOL AND DRUG INFLUENCE

Dr. Jarrad Wagner, Ph.D. F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

#### **ATTESTATION**

I, the undersigned <u>Joshua Smith</u> (Director/Records Custodian), attest under oath pursuant to 12 O.S. §2902 (11)(a) this record was made, at or near the time of the occurrence of the matters set forth by or from information transmitted by a person having knowledge of those matters; was kept in the course of the regularly conducted business activity of the Board; and was made pursuant to the regularly conducted activity. The duplicate record provided is listed below:

1) The attached "Service Overview" is a true and accurate copy of maintenance performed on the Intoxilyzer 8000, serial number 80-005386, in compliance with the agency's Maintenance Policy and Procedures.

This packet consists of 13 pages.

Signed \_\_\_\_\_

Name of Position: Director/Records Custodian

Date of Attestation: 05/24/2022

## Service Overview



3600 N. Martin Luther King Bldg #9 Oklahoma City, OK 73111 405-425-2460

BENCH CHECK DATE: 01-22-2019

**SERIAL:** <u>80-005386</u>

#### **DOCUMENTATION WITHIN SERVICE PACKET:**

- BOT Certificate of Calibration and Operation
- Bench Check
- Bench Check Technician's Report
- BOT Dry Gas Certificate of Analysis
- Mock Subject Test
- BOT Certificates of Analysis on Certified Solutions (4)



## STATE OF OKLAHOMA BOARD OF TESTS FOR ALCOHOL AND DRUG INFLUENCE

Administrative Office at:
Board of Tests for Alcohol and Drug Influence
Post Office Box 36307
Oklahoma City, Oklahoma 73136-2307

### **Certificate of Calibration and Operation**

This is to certify that the calibration of INTOXILYZER 8000, serial number 80-005386, manufactured by CMI Inc., a subsidiary of MPD, Inc. of Owensboro, Kentucky, was tested and found to be in compliance with National Highway Traffic Safety Administration Standard for Devices to Measure Breath Alcohol (F.R., No. 179 48705-48710 Sept. 17, 1993) and with International Recommendation OIML R 126 Edition 1998 (E). Calibration solutions are traceable to NIST standard material 1828.

In addition the operational procedures and functions of this Intoxilyzer have also been tested and found to be in compliance with the Oklahoma Board of Tests for Alcohol and Drug Influence - Rules and Regulations OAC Title 40: 25-1-2 and 30-1-3 and is hereby certified for use as an evidential breath testing device within this state.

CERTIFICATION TYPE	DATE
Manufacturer Calibration	10-13-11
Board of Tests Bench Check	10-19-11
Manufacturer Calibration	07-29-13
Board of Tests Bench Check	08-02-13
Board of Tests Bench Check	09-04-14
Board of Tests Bench Check	07-21-15
Board of Tests Bench Check	04-17-17
Board of Tests Bench Check	01-22-19
A	

OFFITIEIO ATION TO THE



Presentation of this form certifies that the Intoxilyzer, SN: 80-005386, functioned properly at the time of the breath test and hereby certifies the breath test result as valid.

Date:

Re: Request for Log of Test on Intoxilyzer 8000

Name:

#### Dear Sir:

The Log of Test information from a breath test conducted on a Intoxilyzer 8000 is printed in Section 2 of the Affidavit and is no longer retained as a separate form as with the Intoxilyzer 5000. The 8000 uses a nitrogen/ethanol dry gas for the external control test. The dry gas cylinder has the capacity to perform approximately 325 tests and does not require maintenance until its expiration date or can no longer provide the required pressure to perform a control test. The cylinder has an expiration date of two years from date filled at the manufacturer and is programmed in the 8000 upon installation, both dates are printed in Section 2 of the Affidavit. The Intoxilyzer is programmed to disable itself when the date of expiration is reached or the pressure monitor indicates 25 lbs psi.

Due to the reliability and stability of dry gas, maintenance on the 8000 is not required every 30 days or 25 tests, therefore maintenance is no longer performed by officers in the field. The Board of Tests for Alcohol and Drug Influence certifies each Intoxilyzer 8000 for field use by performing a bench check. The bench check is designed to check the operation and procedures programmed in the 8000 to ensure it is functioning properly for use in Oklahoma. A bench check is performed every time a new dry gas cylinder is installed or the instrument has returned from the manufacturer for repair. Once the instrument passes the bench check, all that is required to set it up for field use is to plug it in and turn it on and does not require a permit to set up. The Board of Tests has the ability to monitor each instrument once it is set up in the field.

Each Intoxilyzer has a database of valid permit numbers that it checks prior to allowing a breath test to be performed. The database is updated on a regular basis, if a permit number is invalid that officer will not be able to conduct a breath test on the 8000. If the 8000 prints an affidavit the breath test operator listed in Section 2 had a valid permit at the time of the test and the breath test is valid.

If you have any further questions please feel free to call me.

Sincerely,

J. Robert Blakeburn (OBA11059) State Director of Tests for Alcohol and Drug Influence

pohrt Blokelium.

# STATE OF OKLAHOMA - DEPARTMENT OF PUBLIC SAFETY CMI INTOXILYZER BENCH CHECK PROCEDURE REPORT

ENCH CHECK DATE 01/22/2019	108:37		NSTRUMENT LOCATION 80-00538	6	
rument type Intoxilyzer	INSTRUMENT			80-00	
ncoxilyzei	1 00	TEST RE	SULTS	00-00	3300
22-1-01-0					
Diagnost					
PASS					
Mouth Alcoho	ol Test				
PASS					
RFI Tes	t				
PASS					
About To	A .	Du	plicate Breath	Sampling Re	esults
Abort Te PASS	St	Walas.	7,44	T	Control
		Value	Test 1	Test 2	Test
Improper Sam	ole Test			Loon torous	tus cuasur
PASS		.02	0.018	0.017	0.019
Interferent Dete	ect Test	.04	0.037	0.038	0.039
PASS		1.0	0.000		
No Sample Giv	on Tost	.10	0.098	0.097	0.100
PASS	cii rest	.20	0.199	0.197	0.200
	Un 200		PΔ	SS	
Insufficient Sam	iple Test				
PASS					
Air Blank Contami	nation Test				
PASS					
Range Exceede	ed Test				
PASS					
		CERTIFYING T	ECHNICIAN		
I swear under penalty	of periury that in r	egards to the	above listed i	instrument I	have complied
the Intoxliyzer Bench					
ATURE		V 101 10 10 10 10 10 10 10 10 10 10 10 10		DATE	<u> </u>
/ XV)					
				0	11-22-2019
ATTU TOCHTA				PERMIT#	01.500=
IITH, JOSHUA				10000	016837

CMI Form OKBE1 09/2007

3600 N. Martin Luther King Bldg #9 Oklahoma City, OK 73111 405-425-2460

### **BENCH CHECK TECHNICIAN'S REPORT**

xilyzer Model		Serial #	
00		80-005386	1
S CANISTER		EVD Data	
	014	EXP Date 2-5-2021	
FERENCE	011	2 5 50 51	
	Simulato	r Model	
Guth 2100	Guth 2100	Guth 2100	Guth 2100
550504	<u>Seri</u>		
DR3591	DR3753	DR3594	DR3754
.02	Concer .04	.10	.20
.02	LO.		.20
18020	17410	18070	17340
7. 7. 7. 7. 7. 7. 7. 7. 7. 7. 7. 7. 7. 7	Manufactu	ured Date	
01-09-2018	12-06-2017	02-26-2018	10-09-2017
A Principal State of the	Expiration	on Date	
01-09-2020	12-06-2019	02-26-2020	10-09-2019
	Solution Com	mission Date	
01-17-19	01-17-19	01-17-19	01-17-10
	NCH CHECK ☑CYLINDER REPLA FING ☑ROUTINE N		
MMENTS			

Printed Name

g/210L

0.000

0.081

0.000

0.082

0.000

0.082

0.000

0.082

0.000

0.082

0.000

0.081

0.000

0.082

0.000

0.081

0.000

0.081

0.000

0.081

0.000

0.0815

0.0005

0.6466

SN 80-004923

09:51

Time

09:51

09:52

09:52

09:53

09:53

09:53

09:54

09:54

09:55

09:55

09:56

09:56

09:56

09:57

09:57

09:58

09:58

09:58

09:59

09:59

10:00

ALPHA INSTRUMENT 2

Model 8000

01/16/2019

Test

Air Blank

Cal Check

Air Blank

Std Deu

Cal Check Stats Average

Rel Std Deu(%)

Intoxilyzer - Alcohol Analyzer

Certificate of Analysis – Dry Gas					
01-16-2019	,080				
Date of Analysis	Labelled target value (g/210L)				
014	.0615				
Cylinder #	Average test result				



and Signature BRT Form

In accordance with BOT Policy and Procedure Statement BRT-2.1.0, Maintenance, the above referenced dry gas cylinder is suitable for use as an external reference with the Intoxilyzer 8000.

3 • Jacksonville, IL 62651-0790

21/-243-2103 • FAX: 21/-243-7634 • www.ilmoproducts.com

### Certificate of Analysis

Certificate ID:

11669

Part #:

BAC105L080T

Cylinder Size:

105L

Lot Number:

34418080A3

Expiration:

2/5/2021

Recyplia Polliplia

0.080 BAC (For the calibration of instruments used to determine breath alcohol concentration)

Contents:

105 Liters @ 1000 psig 70°F (21°C)

Analytical

Reported

Accuracy

Analytical

Component:

Concentration:

(U, k=2):

Method:

Ethanol Nitrogen 208 ppm Balance +/-0.002 BAC(G/210L) NDIR

[5.2 ppm]

Distributed by:

CMI Inc.

316 East Ninth Stree Owensboro, KY 423 Phone 866-835-069( www.alcoholtest.com Operator's Signature

\*NIST Traceable Reference Material Cylinder No. CC274523 / Job No. 09160306 Certified 362.2 µmol/mol Ethanol in Nitrogen

Store in dry area, away from sources of heat, ignition and direct sunlight. Do not allow storage area to exceed 52 °C (125 °F).

Specialty Gas Lab Tech

01-03-19

Date

PJLA
Calibration and Testing
Accreditation #61895

The calibration results within this certificate were obtained using equipment and standards capable of producing analytical results traceable to NIST, and apply only to the items contained on this certificate. ILMO Products Company makes no warranty or representation as to the suitability of the use of any information provided for any particular purpose. The information use is at the sole discretion and risk of the user. Liability shall be limited to established replacement cost of this material or service.

#### 8/13

# OFFICER'S AFFIDAVIT AND NOTICE OF REVOCATION/DISOUALIFICATION

ARREST DATE		L V OCA		COUNTY#	CIT			TATION#			S	Sectio
01/22/201	9	1000		55	9			MOCK				
RREST LOCATION					MOCK				COUNT	, LAHO	MA	
JBJECT NAME MOCK, SUB	JECT TEST				DATE OF	BIRTH 09/1:	0.61	SEX M	RACE	HEIGHT	WEI	GHT
DDRESS			CITY		<u> </u>	09/1	ST	ATE	Į VV	ZIP COD	L E	65
3600 MART	IN LUTHER	KING A	OKLA	HOMA CIT	ry	STATE	CLASS	OK	RESTRIC	731	.33 ENDORS	EMEN
r08345487	1			06/30/2		OK	A				M	
HICLE MAKE MOCK		MOCK	5		TAG <b>[M</b> ]	ÖCK			ate O <b>K</b>	CMV N	PLACARI <b>N</b>	)
hway, street, turnpike, put uence of alcohol and/or ot sscribe driving behavior or MOCK SUBJEC scribe person's condition (	T TEST. ROU	olic place, or upon a privat as prohibited by law.	e road, street, alley or lar	ne which provides acce	ess to one or	more single o	or multi-fai	mily dwelli	ng within th	his state wi	nile under	the
NONE		THE PERSON	WAS INFORMED OF TH	IE IMPLIED CONSEN	T TEST REC	)UEST						
	SIS RECORD, REPORT C				ND MAINTE	NANCE REC	ORD.				Se	ectio
QUIPMENT I – 8000				specialist BOARD OF	TES	TS						
andard Dry			TARGET VALUE	PER	RMIT NUMBE	R						
NUFACTURER			LOT#	DAT	E CYLINDEI	RINSTALLED	)			R EXPIRE		
LMO SERVATION START   TE	EST DATE	MOUTHPIECE	344180			/2019		02	2/05	/202	1	
	01/22/2019	9 G		e (or greater) de s observed prio								
rest	g/210L	Time	to th	e time of test ac	dministra	tion.			10° 1	THE S	TATE	
Air Blank Subject Test	0.00	10:36 10:37		red said Brea						N.		: \\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\
Breath Volum	e 1.382 LITER		1 -	e with the ap	-	-	ing	(			**	
Air Blank Wait	0.00	10:37 10:39	procec	lure of the B	oaru oi	rests.		,			翻:	20
Air Blank	0.00	10:39							10/	***		/
Subject Test Breath Volum	0.00 e 1.468 LITER	10:40	SIGNATURE OF C	PERATOR						1907		
Air Blank	0.00	10:40	NAME					┨		D	1.	
Cal Check Air Blank	0.08 0.00	10:40 10:41		SHUA						est Re AC g/2		
			BADGE # 000004		MIT # 0 0 0 0 1	L6837				_		
			AGENCY BOARD OF	трстс				1 L	(	0.0	0	
riod of one hundred state will be disquammercial motor vehing Administrative Hear son at 3600 N. MLK	eighty (180) days or alified if you refused icle. To appeal the cring within fifteen (4 Ave. or by mail to Let I personally hand-	more, if you refuse or failed the state's driver's license revo (15) days from the clegal Division, Departed a copy of	d or failed the state test(s) (operating a cation/disqualificatidate you received artment of Public Sathe above containi	e's test(s). In add a CMV or CDL ho ion, the Departme notice of revoca afety, P.O. Box 1 ing the Notice of I	lition, your older) OR ent of Pubertion or of 1415, OK	commerce a test resulic Safety lisqualific C, OK, 73	cial drivi ult of .04 must re cation. 136 erson n	ing privi 4 or mor eceive y Submit amed ir	lege in re while our writt your wr	operatin ten requi r <b>itten</b> red	ig a est for quest in	ectio
	e 12 O.S. Section 42									rrect."	Se	ction
e: 01/22/20	19 Place (locatio	n when signed):		Qi	ionature c	f arresting	ı officer					
1E	SHUA		*****	AGENCY				•				
GE	AGENCY ADDRESS	anic:		BOARD O		STS IP CODE		ONE				
000000435		JENCY				0000	U   "	4054	1252	460"		
IER WITNESSES:			TITLE	ADDRESS				PHONE			Se	ctio
								LIONE				
NAME			TITLE	ADDRESS				PHONE				

#### 9/13

# OFFICER'S AFFIDAVIT AND NOTICE OF REVOCATION/DISQUALIFICATION

	K	EVUCA	110	N/DISQUA	LIF	ICATIO	IN				Section
ARREST DATE 01/22/2019	***************************************	ARREST TIME		COUNTY 55	#	9 0		ATION #			
ARREST LOCATION W		1000		33	CITY	OCK	1.	1001	COUNT	, LAHO	 М Д
SUBJECT NAME MOCK, SUBJI	ECT TEST				D	ATE OF BIRTH		SEX	RACE	HEIGHT	WEIGHT
ADDRESS				CITY		08/09/19	STA		W	ZIP CODE	
3600 MARTIN	LUTHER	KING A		OKLAHOMA (		STATE (	CLASS	K	DESTRI	731	33 ENDORSEMENTS
T083454871				06/30		19 OK	A				M
VEHICLE MAKE MOCK		MOC]	K			TAG # MOCK			OK	N CMV	PLACARD <b>N</b>
On the above date, time, and loc highway, street, tumpike, public p filuence of alcohol and/or other (Describe driving behavior or circ MOCK SUBJECT	parking lot, or other publintoxicating substances umstances):	ic place, or upon a priva as prohibited by law.	ite road, stre	et, alley or lane which provide	s access to	one or more single or n	nulti-fam	ontrol of ily dwell	a motor ve ing within t	ehicle upon a this state wh	a public road, ille under the
Describe person's condition (odo	r, actions, etc.):			~~							
NONE											
		<b>TUP PP</b>									
BREATH-ALCOHOL ANALYSIS	RECORD, REPORT OF			RMED OF THE IMPLIED CO			PD.				Section 2
EQUIPMENT		SERIA	L NUMBER	SPECIALIST			·	****			30011011 2
I-8000 STANDARD		80	-005	386 BOARD		TESTS NUMBER					
Dry				.080	168	37					
MANUFACTURER TIMO			LOT	# 4418080A3		22/2019				/202	
DBSERVATION START TEST		MOUTHPIECE		5-minute (or greater			on .		2,03	7202	
1005   01	./22/2019	)   G		period was observed						77 00	
Test	g/210L	Time		to the time of te	st admi	nistration.			(3)	THE ST	ATE
			T 1	· · · · 1 · · 1 F	1	A1 1 175 /					1.1/2/
Air Blank Subject Test	0.00	10:36 10:37		ministered said E			1	- 1		ON	
Breath Volume	1.382 LITER		con	npliance with the	~ ~	-	ng	- 1	E		12+1 E
Air Blank	0.00	10:37	1	procedure of the	e Boar	d of Tests.		1			***
Wait Air Blank	0.00	10:39 10:39					İ		13/1	***	* 15
Subject Test	0.00	10:40					_			1907	
Breath Volume	1.468 LITER		SIGN	IATURE OF OPERATOR						1801	
Air Blank Cal Check	0.00 0.08	10:40 10:40	NAME					Г	7	Coat Das	14
Air Blank	0.00	10:41	SMIT		In			1		Test Res rAC g/2	
			BADGE #	0000435	PERMIT #	, 0016837				_	
			AGENCY		1				(	0.0	0
FFECTIVE THIRTY (30)	DAYS FROM DA	TE SERVED No		RD OF TESTS	t vour de	iver's license is	woko -	or do	nind for	^	Section 3
eriod of one hundred eig	hty (180) days or	more, if you refuse	ed or faile	d the state's test(s). In	addition	i, your commercia	vokea I drivin	or de a priv	ilea ior	а	Section 5
nis state will be disqualifi	ed if you refused o	or failed the state's	test(s) (c	perating a CMV or CD	L holder	OR a test result	of .04	or mo	re while	operating	ga
ommercial motor vehicle n Administrative Hearing	. To appeal the d	river's license revo	date vou	squalification, the Department	artment o	of Public Safety m	ust rec	eive y	our writ	ten reque	est for
erson at 3600 N. MLK Av	e. or by mail to L	egal Division, Dep	artment o	of Public Safety, P.O. B	ox 1141	5, OKC, OK, 7313	6	Jubilli	. your w	iiiteii req	uest III
ERVING OF NOTICE: 1							son na	med ir	n Section	n 1 above	€.
ate served: 01/22	2/2019	Signature of office	r:					_ Bad	dge #:		
accordance with Title 12										rrect."	Section 4
ate: 01/22/201						ture of arresting o	-				
AME SMITH, JOS	HUA			AGENCY BOARD	OF	TESTS					
	AGENCY ADDRESS	TENICV				ZIP CODE	PHON		4050	4.60"	
0000000433	INTER-AC	2 ロバイ I				00000	"4	:054	4252	460"	



Certificate of Analysis - Wet Bath 09-26-2018

Labelled target value (g/210L)

0.0204

**8286**0/13



590 NORTH 67th STREET O HARRISBURG, PA 1711 4511 9 TELEPHONE: 717-564-6470

**BOT Technician Name and Signature** 

BRT Form 6.0

In accordance with BRT - 2.1.0, Maintenance, the above referenced simulator solution is suitable for use as an external reference in maintenance of the intoxilyzer 8000.

Average test result (g/210L)

CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 18020 Alcohol Reference Solution for Simulator were analyzed by chromatography on January 11, 2018, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.0238% (w/vol) ethyl alcohol. The expiration date for this lot number is January 9, 2020 at 11:59 PM.

When used in a calibrated Simulator, operating at +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.020 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Signature

ALPHA INSTRUMENT

Model 8000

09/26/2018

Test

Air Blank

Cal Check

Air Blank

Cal Check

Air Blank

Cal Check

Air Blank

Cal Check

Air Blank Cal Check

Air Blank

Cal Check

Air Blank

Cal Check

Air Blank

Cal Check

air Blank

Tal Check

ir Blank

'al Check

lir Blank

al Check Stats Average

Rel Std Deu(%)

Intoxilyzer - Alcohol Analyzer

g/210L

.0.000 0.020

0.000

0.020

0.000

0.020

0.000

0.020

0.000

0.020

0.000

0.020

0.000

0.021

0.000

0.021

0.000

0.021

0.000

0.021

0,000

0.0204 0.0005

2.5312

SN 80-002591

15:36

Time

15:36

15:37

15:37

15:38

15:39

15:39

15:40

15:41

15:41

15:42

15:42

15:43

15:44

15:44

15:45

15:46

15:46

15:47

15:47

15:48

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN08031603 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.



TELEPHONE: 747 FC4 F 470 (esult (g/210L)



BRT Form 6.0

BOTTechnician Name and Signature

CERTIFICATE OF

In accordance with BRT – 2.1.0, Maintenance, the above referenced simulator solution is A suitable for use of an external reference in maintenance of the Intoxilyzer 8000.

Certified Alcohol Reference Solution for Simulator

RLPHR INSTRUMENT Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-802591 09/26/2018 15:54

Test g/210L Air Blank 15:54 Cal Check 0.040 15:55 Air Blank 0.000 15.56 Cal Check 0.039 15:56 Air Blank 0.000 15:57 Cal Check 0,039 15:58 Air Blank 15:58 Cal Check 0.040 15:59 Air Blank 0.000 15:59 Cal Check 0.040 16:00 Air Blank 0.000 16:01 Cal Check 0.040 16:01 Air Blank 0.000 16:02 Cal Check 0.039 16:03 Air Blank 0.000 16:03 Cal Check 0.040 16:04 Air Blank 0.000 16:04 Cal Check 0.040 16:05 Air Blank 0.000 16:06 Cal Check 0.039 Air Blank 0,000 Cal Check Stats Average 0.0396 Std Deu 0.0005 Rel Std Dev(2) 1,3040

Random Samples of Lot Number 17410 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on December 7, 2017, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.0481% (w/vol) ethyl alcohol. The expiration date for this lot number is December 6, 2019 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.040 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

16837 Signature


Operator

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN12181501 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.



590 NORTH 67th STREET AHARRISBURG, PA 17111-4511 TELEPHONE: 717-564-5470

In accordance with BRT - 2.1.0, Maintenance, the above referenced simulator solution is suitable for use as an external reference in maintenance of the intoxilyzer 8000.

CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 18070 Alcohol Reference Solution for Simulator were analyzed by chromatography on February 28, 2018, using a Perkin Elmer Gas gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1219% (w/vol) ethyl alcohol. The expiration date for this lot number is February 26, 2020 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Operato

ALPHA INSTRUMENT

Model 8000

09/27/2018

Test

Air Blank

Cal Check

Air Blank

Cal: Check

Air Blank

Std Dev

Cal Check Stats Average

Rel Std Dev(%)

Intoxilyzer - Alcohol Analyzer

g/210L

0,000

0.100

0.000

0.100

0.000

0.099

0.000

0.100

0.000

0.100

0,000

0.100

0.000

0.099

0.000

0.099

0.000

0.099

0.000

0.100

0.000

0.0996

0.0005

0.5184

SN 80-002591

08:28

Tine

08:28

08-29

08:29

08:30

08:31

08:31

08:32

08:33

08:33

08:34

08:34

08:35

08:36

08:36

08:37

18-38

08:38

08:39

08:39

08:40


08:41

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN04271602 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.



GUTH LABORATORIES, INC.

HARRISBURG, PA 17111-4511 • TELEPHONE: 717-564-5470

BRT Form 6.0

In accordance with BRT – 2.1.0, Maintenance, the above referenced simulator solution is suitable for use as an external reference in maintenance of the intoxityzer 8000.

#### CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 17340 of
Alcohol Reference Solution for Simulator were analyzed by
gas chromatography on October 10, 2017, using a Perkin Elmer Gas
Chromatograph Autosystem XL S/N: 610N9030209, and found to contain
0.2410% (w/vol) ethyl alcohol. The expiration date for this lot
number is October 9, 2019 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.200 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Operator of Tignature

ALPHA INSTRUMENT

Model 8000

09/27/2018

Air Blank

Cal Check

Rir Blank

Cal Check

air Blank

Cal Check

Rir Blank

:al Check

lir Blank

al Check

lir Blank

Std Deu

al Check Stats Average

Rel Std Dev(%)

Intoxilyzer - Alcohol Analyzer

0.000

0.196

0.000

0.196

0.000

0.197

0.000

0.197

0.000

0.197

0.000

0.197

0.000

0.197

0.000

0.197

0.000

0.197

0.000

0.197

0.000

0.1968

0.0004

0.2142

5N 80-002591

15:27

Time

15:27

15:28

15:29

15:29

15:30

15:31

15:31

15:32

15:33

15:33

15:34

15:35

15:35

15:36

15:36

15:37

15:38

15:38

15:39

15:40

15:40

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN08101505 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.