

## STATE OF OKLAHOMA BOARD OF TESTS FOR ALCOHOL AND DRUG INFLUENCE

Dr. Jarrad Wagner, Ph.D. F-ABFT Chairman

Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

#### ATTESTATION

I, the undersigned *Joshua Smith* (Director/Records Custodian), attest under oath pursuant to 12 O.S. §2902 (11)(a) this record was made, at or near the time of the occurrence of the matters set forth by or from information transmitted by a person having knowledge of those matters; was kept in the course of the regularly conducted business activity of the Board; and was made pursuant to the regularly conducted activity. The duplicate record provided is listed below:

1) The attached "Service Overview" is a true and accurate copy of maintenance performed on the Intoxilyzer 8000, serial number , in compliance with the agency's Maintenance Policy and Procedures.

This packet consists of pages.

Signed \_\_\_\_\_

Name of Position: Director/Records Custodian

Date of Attestation:

### Service Overview



3600 N. Martin Luther King Bldg #9 Oklahoma City, OK 73111 405-425-2460

BENCH CHECK DATE: 01-23-2020

**SERIAL:** 80-006148

### **DOCUMENTATION WITHIN SERVICE PACKET:**

- BOT Certificate of Calibration and Operation
- Bench Check
- Bench Check Technician's Report
- BOT Dry Gas Certificate of Analysis
- Mock Subject Test
- BOT Certificates of Analysis on Certified Solutions (4)



## STATE OF OKLAHOMA BOARD OF TESTS FOR ALCOHOL AND DRUG INFLUENCE

Administrative Office at:
Board of Tests for Alcohol and Drug Influence
Post Office Box 36307
Oklahoma City, Oklahoma 73136-2307

### Certificate of Calibration and Operation

This is to certify that the calibration of INTOXILYZER 8000, serial number 80-006148, manufactured by CMI Inc., a subsidiary of MPD, Inc. of Owensboro, Kentucky, was tested and found to be in compliance with National Highway Traffic Safety Administration Standard for Devices to Measure Breath Alcohol (F.R., No. 179 48705-48710 Sept. 17, 1993) and with International Recommendation OIML R 126 Edition 1998 (E). Calibration solutions are traceable to NIST standard material 1828.

In addition the operational procedures and functions of this Intoxilyzer have also been tested and found to be in compliance with the Oklahoma Board of Tests for Alcohol and Drug Influence - Rules and Regulations OAC Title 40: 25-1-2 and 30-1-3 and is hereby certified for use as an evidential breath testing device within this state.

CERTIFICATION TYPE	DATE
Manufacturer Calibration	01-24-14
Board of Tests Bench Check	02-04-14
Manufacturer Calibration	01-29-16
Board of Tests Bench Check	03-24-16
Board of Tests Bench Check	03-06-18
Board of Tests Bench Check	01-23-20



Presentation of this form certifies that the Intoxilyzer, SN: 80-006148, functioned properly at the time of the breath test and hereby certifies the breath test result as valid.

Date:

Re: Request for Log of Test on Intoxilyzer 8000

Name:

#### Dear Sir:

The Log of Test information from a breath test conducted on a Intoxilyzer 8000 is printed in Section 2 of the Affidavit and is no longer retained as a separate form as with the Intoxilyzer 5000. The 8000 uses a nitrogen/ethanol dry gas for the external control test. The dry gas cylinder has the capacity to perform approximately 325 tests and does not require maintenance until its expiration date or can no longer provide the required pressure to perform a control test. The cylinder has an expiration date of two years from date filled at the manufacturer and is programmed in the 8000 upon installation, both dates are printed in Section 2 of the Affidavit. The Intoxilyzer is programmed to disable itself when the date of expiration is reached or the pressure monitor indicates 25 lbs psi.

Due to the reliability and stability of dry gas, maintenance on the 8000 is not required every 30 days or 25 tests, therefore maintenance is no longer performed by officers in the field. The Board of Tests for Alcohol and Drug Influence certifies each Intoxilyzer 8000 for field use by performing a bench check. The bench check is designed to check the operation and procedures programmed in the 8000 to ensure it is functioning properly for use in Oklahoma. A bench check is performed every time a new dry gas cylinder is installed or the instrument has returned from the manufacturer for repair. Once the instrument passes the bench check, all that is required to set it up for field use is to plug it in and turn it on and does not require a permit to set up. The Board of Tests has the ability to monitor each instrument once it is set up in the field.

Each Intoxilyzer has a database of valid permit numbers that it checks prior to allowing a breath test to be performed. The database is updated on a regular basis, if a permit number is invalid that officer will not be able to conduct a breath test on the 8000. If the 8000 prints an affidavit the breath test operator listed in Section 2 had a valid permit at the time of the test and the breath test is valid.

If you have any further questions please feel free to call me.

Sincerely,

J. Robert Blakeburn (OBA11059) State Director of Tests for Alcohol

O. pohr &Blokelune

and Drug Influence

# STATE OF OKLAHOMA - BOARD OF TESTS FOR ALCOHOL AND DRUG INFLUENCE CMI INTOXILYZER BENCH CHECK PROCEDURE REPORT

CH CHECK DATE 1/23/2020	13:18:0		INSTRUMENT LOCATION 80-006148			
RUMENT TYPE	INSTRUMENT MO		INSTRUMENT SERIAL NUMBER 80-006148			
ntoxilyzer	I-8000	TEST RE	SIII TS	80-00	6148	
		TESTINE	30113			
Diagnostic PASS	es					
Mouth Alcoho	l Test					
RFI Test						
PASS						
Abort Tes		Du	plicate Breath	Sampling Re	esults	
PASS	it.	Value	Test 1	Test 2	Control Test	
Improper Samp	le Test		20102	1577555	3 33 33	
PASS	634	.02	0.017	0.018	0.019	
Interferent Dete	ct Test	.04	0.038	0.037	0.038	
PASS		.10	0.097	0.096	0.097	
No Sample Give	en Test					
PASS		.20	0.200	0.199	0.201	
Insufficient Samp	ple Test		PA	SS		
PASS	717 0 TO					
Air Blank Contamin	nation Test					
PASS						
Range Exceede	d Test					
PASS						
	CI	ERTIFYING T	ECHNICIAN			
swear under penalty						
he Intoxliyzer Bench	Checking Procedure	approved	by the State o	f Oklahoma.		
ATURE				DATE		
				01	-23-202	
-y						
				PERMIT #		



3600 N. Martin Luther King Bldg #9 Oklahoma City, OK 73111 405-425-2460

### **BENCH CHECK TECHNICIAN'S REPORT**

000 AS CANISTER DT #		80-006148	}
OT#			
0219080A5		EXP Date	
	# 020	01-05-7	2022
EFERENCE			
Guth 2100	Simulato Guth 2100		Cuth 2100
Gulfi 2 100		Guth 2100	Guth 2100
DR3591	DR3753	DR3594	DR3754
	Concer	The state of the s	
.02	.04	.10	.20
	The state of the s	T#	
19050	19080	19301	19200
	Manufacti		
02-04-2019	03-04-2019	10-30-2019	08-06-2019
	Expiration		
02-04-2021	03-04-2021	10-30-2021	08-06-2021
	Solution Com	mission Date	
1-22-2020	4-22-2020	01-22-2020	01-22-2020
	CYLINDER REPLATING ROUTINE		
$\sim$			

q/210L

0.000

0.081

0.000

0.081

0.000

0.081

0.000

0.081

0.000

0.081

0.000

0.081

0.000

0.081

0.000

0.080

0.000

0.080

0.000

0.081

0.000

0.0808

0.0004

0.5218

5N 80-004923

14:01

Time

14:01

14:01

14:02

14:02

14:03

14:03

14:04

14:04

14:04

14:05

14:05

14:06

14:06

14:06

14:07

14:07

14:08

14:08

14:09

14:09

14:10

ALPHA INSTRUMENT 2

Model 8000

01/12/2020

Test

Air Blank

Cal Check

Air Blank

Std Deu

Cal Check Stats Average

Rel Std Deu(%)

Intoxilyzer - Alcohol Analyzer

Certificate	of Analysis – Dry Gas	
01-12-2020	.080	
Date of Analysis	Labelled target value (g/210L)	
H-020	1881	
Cylinder#	Average test result	
1884 Snith		
BOT Technician Name and Signature		specialty gases

In accordance with BOT Policy and Procedure Statement BATE 220 Management P.O. Box 790 • Jacksonville, IL 62651-0790 referenced dry gas cylinder is suitable for use as an external reference with the Intoxilyzer 8000.

217-245-2183 • Fax: 217-243-7634 • www.ilmoproducts.com

### Certificate of Analysis

Certificate ID:

12515

Part #:

BAC105L080T

Cylinder Size:

105L

Lot Number:

30219080A5

Expiration:

1/5/2022

0.080 BAC (For the calibration of instruments used to determine breath alcohol concentration)

105 Liters @ 1000 psig 70°F (21°C)

Analytical

Reported

Accuracy

Analytical

Component:

Nitrogen

Concentration:

(U, k=2):

Method:

Ethanol +/-0.002 BAC(G/210L) NDIR 208 ppm

[5.2 ppm] Balance

Distributed by:

CMI Inc.

316 East Ninth Stree Owensboro, KY 423 Phone 866-835-0690 www.alcoholtest.con

\*Traceable to:

Certified Reference Material - 262.4 µmol/mol

Ethanol in Nitrogen - Serial No. GN0015026 Lot No. 050319E11

Store in dry area, away from sources of heat, ignition and direct sunlight. Do not allow storage area to exceed 52 °C (125 °F).

Specially Gas Lab Tech



The calibration results within this certificate were obtained using equipment and standards capable of producing analytical results traceable to NIST, and apply only to the items contained on this certificate. ILMO Products Company makes no warranty or representation as to the suitability of the use of any information provided for any particular purpose. The information use is at the sole discretion and risk of the user. Liability shall be limited to established replacement cost of this material or service.

#### 8/13

# OFFICER'S AFFIDAVIT AND NOTICE OF REVOCATION/DISOUALIFICATION

۷/	ч	3
וכ	-	

		E VOCA		DQUAL	/II. I.C	AII	)				Section
ARREST DATE 01/23/2020		ARREST TIME 1350		COUNTY # 55	- 1	Y # 8 <b>4</b>		TION#			***************************************
ARREST LOCATION NONE					CITY	·····	1 -4	<b></b>	COUNT		
SUBJECT NAME					DATE O	F BIRTH		SEX	RACE	HEIGHT	WEIGHT
NORTH, GAYI	LON LEVE	RNE	CITY		08,	/09/19	61 STAT	M	W	510	165
3600 MARTIN	I LUTHER	KING A		HOMA CI	-		0:			7313	3 3
DRIVER LICENSE # T083454871				06/30/		STATE OK	CLASS A		RESTRI		NDORSEMENTS <b>M</b>
VEHICLE MAKE		MODEL	-	100/30/	TA	3#		STA		CMV P	LACARD
MOCK On the above date, time, and local	ation, the above named	MOCK		nds to believe the per		IOCK	hysical cor		Motor ve		N nublic road
highway, street, turnpike, public pinfluence of alcohol and/or other is (Describe driving behavior or circic MOCK SUBJECT  Describe person's condition (odor NONE	intoxicating substances umstances): TEST. ROU'	as prohibited by law.		·				y dwellir	ng within t	his state whil	le under the
110112		THE PERSON	WAS INFORMED OF TI	HE IMPLIED CONSE	NT TEST RE	OUEST					
BREATH-ALCOHOL ANALYSIS	RECORD, REPORT O			· · · · · · · · · · · · · · · · · · ·			RD.			***************************************	Section 2
I - 8000			NUMBER - 006148	BOARD O	F TES	TTC					
STANDARD		00	TARGET VALUE	PE	RMIT NUMB						
Dry			0.080 LOT#		16837	RINSTALLED		DATE	CYI INDE	R EXPIRES	
ILMO			302190	080A5 (	01/23	/2020				/2022	
OBSERVATION START TEST 1430 01	DATE -/23/2020	MOUTHPIECE		e (or greater) d							
RESULT TABLE	,,			is observed pri		_	;			THE ST	AT
Test	g/210L	Time	to tr	ne time of test a	administr	ation.			( St. )	· : : : : : : : : : : : : : : : : : : :	
Air Blank	0.00	14:57	I administe	red said Bre	ath-Alc	ohol Tes	t in				
Subject Test Breath Volume	0.00	14:58	complianc	e with the ap	pplicab	le operati	ng		0 **		
Air Blank	1.460 LITER 0.00	.s 14:59	procee	dure of the H	Board o	f Tests.		- /	A ***		****   E
Wait Air Blank	0.00	15:00 15:01						`	15/1/		(A)
Subject Test	0.00	15:01	SIGNATURE OF	ODERATOR						1907	
Breath Volume Air Blank	1.355 LITER 0.00	.S 15:03	SIGNATURE OF	OPERATOR							
Cal Check	0.07	15:03	NAME TO A	71117					]	Test Res	ult
Air Blank	0.00	15:04	SMITH, JOS BADGE#	PE	RMIT#					rAC g/2	
			00000004	435 (	00000	16837			(	0.0	0
			BOARD OF	TESTS						•••	
EFFECTIVE THIRTY (30) DAY days or more, if you refused or (operating a CMV or CDL holde Safety must receive your writter The approved hearing requee Division, Department of Public a fine of not more than Five Huperson's license, permit, or other approved than the control of t	failed the state's test er) OR a test result of n request for an admi st form is available a Safety, P.O. Box 114 ndred Dollars (\$500.0 er evidence of driving	s). In addition, your co .04 or more while opers instrative hearing on the thtps://www.ok.gov/15, OKC, OK 73136. It .10), or by both such fine privilege is in the custo	mmercial driving privile ating a commercial moi e Department approv dps/documents/Reqh shall be a misdemeand and imprisonment, for dy of a law enforcement	ege in this state will I tor vehicle. To appered form, within fifted dearing.pdf. Submitor punishable by important applies any person to applient officer or the Depuishable in the Depuishable will be subject to the Depuishable in the Depuishable will be subject to the subject	be disqualified the driver the dr	ed if you refuse 's license revor s from the date ng request form or not less than val or a replace	ed or failed cation/dis- you recein in person seven (7 ement lice	I the sta qualifica ved not n at 360 n days r	ate's test( ation, the ice of rev 0 N. MLA nor more sperate a	(s) Departmen vocation or o K Ave. or by than six (6) motor vehice	disqualification. mail to Legal months, or by cle while the
SERVING OF NOTICE: Date served: 01/23									Sectio ge #:	n 1 above	ı.
In accordance with Title 1	2 O.S. Section 42	6, "I state under <b>pe</b>	nalty of perjury u	inder the laws of	Oklahom	a that the fo	regoing	is true	and co	orrect."	Section 4
Date: 01/23/202	O Place (location	n when signed):			Signature	of arresting	officer:_				
NAME SMITH, JOSH	UA			BOARD (	OF TE	STS					
BADGE 0000000435	AGENCY ADDRESS INTER-A	GENCY				ZIP CODE 0000	PHON 4 (		2524	160	
OTHER WITNESSES:											Section 5
NAME 1			TITLE	ADDRESS			P	HONE			
NAME			TITLE	ADDRESS			P	HONE			
2											

### 9/13

# OFFICER'S AFFIDAVIT AND NOTICE OF REVOCATION/DISOUALIFICATION

	RE	VOCA	rion/disc	QUALI	FIC	CATION	N				Section 1
ARREST DATE 01/23/2020		ST TIME 550		COUNTY#		TY# 34	CITAT	ONE			
ARREST LOCATION NONE	1 = 5			CIT			110		COUNT		
SUBJECT NAME NORTH, GAYI	ON LEVERNE				DATE	) BIRTH / 09/196		EX M	RACE <b>W</b>	HEIGHT 510	WEIGHT 165
ADDRESS 3600 MARTIN			OKT <sub>1</sub> AHC	MA CITY		, ,	STATE	= '		ZIP CODE	E
DRIVER LICENSE # T083454871			EX	PIRATION DATE 06/30/20		STATE CLA	ASS		RESTRI		ENDORSEMENTS M
VEHICLE MAKE MOCK		MOCK		0,00,20	TA	G# MOCK		STA		CMV N	PLACARD N
On the above date, time, and lock highway, street, turnpike, public p influence of alcohol and/or other i (Describe driving behavior or circu MOCK SUBJECT	arking lot, or other public place toxicating substances as prounds and the substances are prounds and the substances. TEST. ROUTIN	e, or upon a private hibited by law.	e road, street, alley or lane wh	nich provides access	to one o	or more single or mul	lti-family	trol of a	motor ve	hicle upon his state wh	a public road, nile under the
NONE		THE PERSON V	VAS INFORMED OF THE IN	APLIED CONSENT	TEST RI	EQUEST					
BREATH-ALCOHOL ANALYSIS EQUIPMENT	RECORD, REPORT OF THE			LOG OF TEST ANI	TAIAM C	ENANCE RECORD.	•				Section 2
I-8000 STANDARD		80-	-006148 BO	DARD OF	TE:	STS					
Dry MANUFACTURER			0.080	16	837			DATE	OVI INIDE	D EVDIDE	
ILMO			30219080	0A5 01	1/23	1/2020		01	/05	7202	2
OBSERVATION START TEST 1430 01 RESULT TABLE Test	9/210L	MOUTHPIECE G Time		r greater) depr bserved prior me of test adr	to and	continuing	1		OF .	THE S	PARE OF
Air Blank Subject Test Breath Volume Air Blank Wait Air Blank	0.00 0.00 1.460 LITERS 0.00	14:57 14:58 14:59 15:00 15:01	I administered compliance w procedur		licab	le operating	- 1	TA PA			N THE STATE OF THE
Subject Test Breath Volume	0.00 1.355 LITERS	15:03	SIGNATURE OF OPER	RATOR			-			1907	
Air Blank Cal Check Air Blank	0.00 0.07 0.00	15:03 15:04 15:04	NAME SMITH, JOSHU BADGE# 00000043 AGENCY BOARD OF T	5 0 0		16837			Bı	Test RerAC g/2	210L
EFFECTIVE THIRTY (30) DAY days or more, if you refused or (operating a CMV or CDL holde Safety must receive your writtet The approved hearing reques Division, Department of Public of a fine of not more than Five Hu person's license, permit, or othe SERVING OF NOTICE:	failed the state's test(s). In r) OR a test result of .04 or n request for an administrat it form is available at http. Safety, P.O. Box 11415, Ohdred Dollars (\$500.00), or evidence of driving privile	addition, your cormore while operative hearing on the s://www.ok.gov/cC, OK 73136. It is by both such fine toge is in the custod	ation is given that your driv immercial driving privilege in ting a commercial motor ve e Department approved fr ips/documents/ReqHeari shall be a misdemeanor pu and imprisonment, for any by of a law enforcement off	er's license is revon this state will be of chicle. To appeal to prm, within fifteen mg.pdf. Submit you nishable by impriss person to apply fo icer or the Departn	disqualifiche drive (15) day our heari conment r a rene nent.	ed if you refused o r's license revocati is from the date you ng request form in for not less than se wal or a replaceme	or failed ion/disq u receiv person even (7) ent licen	the stat qualificat yed notic at 3600 days notic se to op	te's test( tion, the ce of rev O N. MLH or more perate a	s) Departme cocation or Ave. or b than six (6 motor veh	disqualification. by mail to Legal b) months, or by sicle while the
Date served: 01/23	10000	ature of officer				ion to the perso				II I ADOV	e.
In accordance with Title 1										rrect."	Section 4
Date: 01/23/202				Sig		of arresting off	-	11,			
SMITH, JOSH				ENCY BOARD OF	TE	ESTS					
0000000435	AGENCY ADDRESS INTER-AGEI	ICY				ZIP CODE 00000	PHONE <b>4</b> 0		2524	160	

EUTH LABORATORIES, INC.

250 NORTH 67th STREET • HARRISBURG, PA 17111-4511 • TELEPHONE: 717-564-5470

BOT Technician Name and Signature

Lot Number

BRT Form 6.0

In accordance with BRT – 2.2.0, Maintenance, the above referenced simulator solution is suitable for use as an external reference in maintenance of the Intoxilyzer 8000.

### CERTIFICATE OF ANALYSIS

RLPHA INSTRUMENT 2 Intoxilyzer - Alcohol Analyzer Model 8000 . SN 80-804923 12/10/2019 14:52

5M 801-0049233 Certified Alcohol Reference Solution for Simulator

Test g/210L Time Air Blank 14:53 Cal Check 0.019 14:53 Air Blank 0.000 14:54 Cal Check 0.020 14:55 Air Blank 0,000 14:55 Cal Check 0.019 14:56 Air Blank 0.000 14:56 Cal Check 0.020 14:57 Air Blank 0.000 14:58 Cal Check 0.019 14:58 Air Blank 0.000 14:59 Cal Check 0.020 15:00 Air Blank 0.000 15:00 Cal Check 0.019 15:01 Air Blank 0.000 15:01 Cal Check 0.020 15:02 Air Blank 0.000 15:02 Cal Check 0.020 15:03 Air Blank 0.000 15:04 Cal Check 0.020 15:04 Air Blank 0.000 15:05 Cal Check Stats Average 0.0196 Std Deu 0.0005

Random Samples of Lot Number 19050 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on February 5, 2019, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.0236% (w/vol) ethyl alcohol. The expiration date for this lot number is February 4, 2021 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.020 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

1937 Perator's Signature

2,6345

Rel Std Deu(%)

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN08031603 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.

Certificate of	Analysis - Wet Bath
12,10-19	,940
Date of Analysis	Labelled tar (R value (g/210L)
19080	1, 040
Lot Number	Aldra echyrsolic (219)

ABORATORIES, INC.

590 NORTH 67th STREET 9 HARRISBURG, PA 17111-4511 9 TELEPHONE: 717-564-5470

BOT Technician Name and Signature

BRT Form 6.0

In accordance with BRT -2.2.0, *Maintenance*, the above referenced simulator solution is suitable for use as an external reference in maintenance of the Intoxilyzer 8000.

#### CERTIFICATE OF ANALYSIS

ALPHA INSTRUMENT 2 Intoxilyzer - Alcohol Analyzer · Model 8000 SN 80-004923 12/10/2019 15:16

Certified Alcohol Reference Solution for Simulator

Test g/210L Time Air Blank 0.000 15:17 Cal Check 0,040 15:17 Air Blank 0.000 15:18 Cal Check 0.040 15:19 0.000 Air Blank 15:19 Cal Check 0.040 15:20 Air Blank 0.000 15:20 Cal Check 0.040 15:21 Air Blank 0.000 15:22 Cal Check 0.040 15:22 Air Blank 0,000 15:23 Cal Check 0.040 15:24 Air Blank 0,000 15:24 Cal Check 15:25 0.040 Air Blank 0.000 15:25 Cal Check 0.039 15:26 Air Blank 0.000 15:27 Cal Check 0.039 15:27 Air Blank 0,000 15:28 Cal Check 0.039 15:29 Air Blank 0.000 15:29 Cal Check Stats 0.0397 Average Std Dev 0.0004 Rel Std Deu(%) 1,2167

Random Samples of Lot Number 19080 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on March 5, 2019, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.0481% (w/vol) ethyl alcohol. The expiration date for this lot number is March 4, 2021 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.040 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

16837 Operator's Signature

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN12181501 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.

Certificate of	Analysis – Wet Bath	
12-10-19	,100	And Charles
Date of Analysis	Labelled target value (g/210L)	1
19301	JE 72 . 100	
Lat Niverland	1 A. S. F. T. T. B. F. C. W. H. (a./2101)	

● HARRISBURG, PA 17111-4511 ● TELEPHONE: 717-564-5470

In accordance with BRT - 2.2.0, Maintenance, the above referenced simulator solution is suitable for use as an external reference in maintenance of the Intoxilyzer 8000.

### CERTIFICATE OF ANALYSIS

ALPHA INSTRUMENT 2

JOSH

Intoxilyzer - Alcohol Analyzer

BOT Technician Name and Signature

12/10/2019

SN 80-004923 Model 8000

Test	g/210L	Time
Air Blank	0.000	15:37
Cal Check	0,100	15:37
Air Blank	0.000	15:38
Cal Check	0.101	15:39
Air Blank	0.000	15:39
Cal Check	0.100	15:40
Air Blank	0.000	15:40
Cal Check	0.100	15:41
Air Blank	0.000	15:42
Cal Check	0,100	15: 42
Air Blank	0.000	15:43
Cal Check	0,099	15:44
Air Blank	0,000	15:44
Cal Check	0.100	15:45
Air BlanK	0.000	15:45
Cal Check	0.099 .	15:46
Air Blank	0.000	15:47
Cal Check	0.100	15:47
Air Blank	0.000	15:48
Cal Check	0.100	15:49
Air Blank	0.000	15:49
Cal Check Stats		
Average	0.0999	
Std Dev`	0.0005	
Rel Std Deu(%)	0.5682	

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 19301 Alcohol Reference Solution for Simulator were analyzed by chromatography on November 1, 2019, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1210% (w/vol) ethyl alcohol. The expiration date for this lot number is October 30, 2021 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

operator's Signature

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN04271602 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.

Certificate	of Analysis – Wet Bath	1000
12-11-19	,200	( SOF
Date of Analysis	Labelled tarket value (g/210L)	
19260	1290	1967
Lot Number	Average Atea resolit (/2/0)	ABORATORIE
Task Smith	590 NORTH 67th STRE	ET O HARRISBURG, PA 17111.

In accordance with BRT – 2.2.0, *Maintenance*, the above referenced simulator solution is suitable for use as an external reference in maintenance of the Intoxilyzer 8000.

### CERTIFICATE OF ANALYSIS

ALPHA INSTRUMENT 2 Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-004923 12/11/2019 08:03

BOT Technician Name and Signature

5N 89-004923 Certified Alcohol Reference Solution for Simulator

Test (	g/210L	. Time
Air Blank	0.000	08:04
Cal Check	0,198	08:04
Air Blank	0.000	08: 05
Cal Check	0.198	08:06
Air Blank	0.000	08: 06
Cal Check	0.197	08: 07
Air Blank	0,000	08; 08
Cal Check	0.197	08:08
Air Blank	0.000	08:09
Cal Check	0.197	08:10
Air Blank	0.000	08:10
Cal Check	0.197	08:11
Air Blank	0.000	08:11
Cal Check	0.197	08:12
Air Blank	0.000	08:13
Cal Check	0.197	08:13
Air Blank	0,000	08:14
Cal Check	0.196	09:15
Air Blank	0.000	08:15
Cal Check	0.198	08:16
Air Blank	0,000	08;17
Cal Check Stats		
Average	0.1972	
Std Deu	0.0006	

Rel Std Deu(%) 0.3207

Random Samples of Lot Number 19200 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on August 7, 2019, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.2418% (w/vol) ethyl alcohol. The expiration date for this lot number is August 6, 2021 at 11:59 PM.

4511 @ TELEPHONE: 717-564-5470

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.200 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

16837 erator's Signature

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN08101505 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.