

STATE OF OKLAHOMA BOARD OF TESTS FOR ALCOHOL AND DRUG INFLUENCE

Dr. Jarrad Wagner, Ph.D. F-ABFT Chairman Post Office Box 36307 Oklahoma City, Oklahoma 73136-2307 Phone: (405) 425-2460 Fax: (405) 425-2490 www.bot.ok.gov

Joshua Smith Director

ATTESTATION

I, the undersigned *Joshua Smith* (Director/Records Custodian), attest under oath pursuant to 12 O.S. §2902 (11)(a) this record was made, at or near the time of the occurrence of the matters set forth by or from information transmitted by a person having knowledge of those matters; was kept in the course of the regularly conducted business activity of the Board; and was made pursuant to the regularly conducted activity. The duplicate record provided is listed below:

1) The attached "Service Overview" is a true and accurate copy of maintenance performed on the Intoxilyzer 8000, serial number 80-006163, in compliance with the agency's Maintenance Policy and Procedures.

This packet consists of 13 pages.

Signed _____

Name of Position: Director/Records Custodian

Date of Attestation: 04/20/2020



Board of Toots for Alcohol and Dray Influence Benetising as startic second as means and weater atoms of two and weater streng



3600 N. Martin Luther King Bldg #9 Oklahoma City, OK 73111 405-425-2460

BENCH CHECK DATE: 04-20-2020

SERIAL: 80-006163

DOCUMENTATION WITHIN SERVICE PACKET:

- BOT Certificate of Calibration and Operation
- Bench Check
- Bench Check Technician's Report
- BOT Dry Gas Certificate of Analysis
- Mock Subject Test
- BOT Certificates of Analysis on Certified Solutions (4)



STATE OF OKLAHOMA BOARD OF TESTS FOR ALCOHOL AND DRUG INFLUENCE

Administrative Office at: Board of Tests for Alcohol and Drug Influence Post Office Box 36(307 Oklahoma City, Oklahoma 7'3136-2307

Certificate of Calibration and Operation

This is to certify that the calibration of **INTOXILYZER 8000, serial number** <u>80-006163</u>, manufactured by CMI Inc., a subsidiary of MPD, Inc. of Owensboro, Kentucky, was tested and found to be in compliance with National Highway Traffic Safety Administration Standard for Devices to Measure Breath Alcohol (F.R., No. 179 48705-48710 Sept. 17, 1993) and with International Recommendation OIML R 126 Edition 1998 (E). Calibration solutions are traceable to NIST standard material 1828.

In addition the operational procedures and functions of this Intoxilyzer have also been tested and found to be in compliance with the Oklahoma Board of Tests for Alcohol and Drug Influence - Rules and Regulations OAC Title 40: 25-1-2 and 30-1-3 and is hereby certified for use as an evidential breath testing device within this state.

CERTIFICATION TYPE	DATE
Manufacturer Calibration	01-23-14
Manufacturer Calibration	03-10-14
Board of Tests Bench Check	06-12-14
Manufacturer Calibration	06-28-16
Board of Tests Bench Check	12-12-16
Board of Tests Bench Check	10-18-17
Board of Tests Bench Check	10-29-18
Board of Tests Bench Check	04-20-20



Presentation of this form certifies that the Intoxilyzer, SN: <u>80-006163</u>, functioned properly at the time of the breath test and hereby certifies the breath test result as valid.

Date:

Re: Request for Log of Test on Intoxilyzer 8000

Name:

Dear Sir:

The Log of Test information from a breath test conducted on a Intoxilyzer 8000 is printed in Section 2 of the Affidavit and is no longer retained as a separate form as with the Intoxilyzer 5000. The 8000 uses a nitrogen/ethanol dry gas for the external control test. The dry gas cylinder has the capacity to perform approximately 325 tests and does not require maintenance until its expiration date or can no longer provide the required pressure to perform a control test. The cylinder has an expiration date of two years from date filled at the manufacturer and is programmed in the 8000 upon installation, both dates are printed in Section 2 of the Affidavit. The Intoxilyzer is programmed to disable itself when the date of expiration is reached or the pressure monitor indicates 25 lbs psi.

Due to the reliability and stability of dry gas, maintenance on the 8000 is not required every 30 days or 25 tests, therefore maintenance is no longer performed by officers in the field. The Board of Tests for Alcohol and Drug Influence certifies each Intoxilyzer 8000 for field use by performing a bench check. The bench check is designed to check the operation and procedures programmed in the 8000 to ensure it is functioning properly for use in Oklahoma. A bench check is performed every time a new dry gas cylinder is installed or the instrument has returned from the manufacturer for repair. Once the instrument passes the bench check, all that is required to set it up for field use is to plug it in and turn it on and does not require a permit to set up. The Board of Tests has the ability to monitor each instrument once it is set up in the field.

Each Intoxilyzer has a database of valid permit numbers that it checks prior to allowing a breath test to be performed. The database is updated on a regular basis, if a permit number is invalid that officer will not be able to conduct a breath test on the 8000. If the 8000 prints an affidavit the breath test operator listed in Section 2 had a valid permit at the time of the test and the breath test is valid.

If you have any further questions please feel free to call me.

Sincerely,

O. pohrt Blackelun

J. Robert Blakeburn (OBA11059) State Director of Tests for Alcohol and Drug Influence

STATE OF OKLAHOMA - BOARD OF TESTS FOR ALCOHOL AND DRUG INFLUENCE CMI INTOXILYZER BENCH CHECK PROCEDURE REPORT

ENCH CHECK DATE TIME PROCEDURE BEG/ 04/20/2020 12:59:0			80-00616	3					
STRUMENT TYPE Intoxilyzer		IMENT MODEL		INSTRUMENT SERIAL NUMBER					
	-	TEST RE	SULTS	50 00	0100				
		100 A							
Diagnostic PASS	cs								
Mouth Alcoho PASS	ol Test								
RFI Test PASS	t		9		de l'an				
Abort Tes	st	Duj	Duplicate Breath Sampling Results						
PASS		Value	Test 1	Test 2	Control Test				
Improper Samp PASS	ole Test	. 02	0.017	0.017	0.018				
Interferent Dete	ect Test	.04	0.036	0.036	0.037				
PASS		.10	0.093	0.093	0.095				
No Sample Give PASS	en Test	.20	0.195	0.193	0.197				
Insufficient Sam PASS	ple Test		PA	SS					
Air Blank Contamir PASS	nation Test								
Range Exceede PASS	ed Test								
INDO									
		CERTIFYING T	ECHNICIAN						
I swear under penalty	of perjury that			instrument. I	have complied				
the Intoxliyzer Bench									
	4			DATE					
Cox				٥Y	20/2020				
ME				PERMIT #					

		Denvel of Tools for Alashal	
		Board of Tests for Alcohol Enhancing public safety by ensuring the ac of blood and breath	seuracy and scientific reliability
BENCH CH	IECK TEC	HNICIAN'S	
Date: 04 / 20 / 20	Start Time:	<u>600</u> End Time	1450
INSTRUMENT		Serial #	
8000		80-00610	03
GAS CANISTER		EXP Date	4
05620080A2	\$024	04-0	5-2022
REFERENCE	Simulato	or Model	
Guth 2100	Guth 2100	Guth 2100	Guth 2100
		ial <u>#</u>	
DR3591	DR3753	DR3594	DR3754
00	Exception products	ntration	00
.02	.04	.10 T#	.20
19050	19080	19301	19200
,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,	25 5 5 5	ured Date	10200
02-04-2019	03-04-2019	10-30-2019	08-06-2019
		on Date	
02-04-2021	03-04-2021	10-30-2021	08-06-2021
	Solution Com		
01110 01			011 +2 20
04-02-20	04-02-20	04-02-20	04-02-20

REASON FOR BENCH CHECK

□ POST REPAIR INCLUDER REPLACEMENT □ TROUBLESHOOTING INCLUDER MAINTENANCE COMMENTS

Technician Signature / Permit # Josh Printed Name

6/13

Certificate	of Analysis – Dry Gas	
04-19-20	.888	TAP
Date of Analysis	Labelled target value (g/210L)	Name and Address of the owner owner owner owner owner owner
#024	,989	ILMA
Cylinder#	Average test result	
Josh Smith	R.	
BOT Technician Name and Signati	Ire	specialty gases

referenced dry gas cylinder is suitable for use as an external 7e beau ganter in the in a soor 790 . Jacksonville, IL 62651-0790

217-245-2183 • Fax: 217-243-7634 • www.ilmoproducts.com

Certificate of Analysis

ALPHA INSTRUMENT 2 Intoxilyzer - Alcohol Analyzer Model 8000 5N 80-004923 04/19/2020 09:45

7/13

Test	g/210L	Time
Air Blank	0.000	09:45
Cal Check	0.080	09:46
Air Blank	0.000	09:46
Cal Check	0.079	09:47
Air Blank	0.000	09:47
Cal Check	0.080	09:47
Air Blank	0.000	09:48
Cal Check	0.080	09:48
Air Blank	0.000	09:49
Cal Check	0.080	09:49
Air Blank	0.000	09:50
Cal Check	0.080	09:50
Air Blank	0.000	09:50
Cal Check	0.080	09:51
Air Blank	0.000	09:51
Cal Check	0.080	09:52
Air Blank	0.000	09:52
Cal Check	0.080	09:52
Air Blank	0.000	09:53
Cal Check	0.080	09:53
Air Blank	0.000	09:54
al Check Stats		
Average	0.0799	
Std Dev	0.0003	
Rel Std Deu(%)	0.3957	

0.080 BAC (For the calibration of instruments used to determine breath alcohol concentration)

105 Liters @ 1000 psig 70°F (21°C) Contents: ALL MORE TO JA

Balance

12811

105L

BAC105L080T

05620080A2

4/5/2022

Component: Ethanol

Certificate ID:

Cylinder Size:

Lot Number:

Expiration:

Part #:

Nitrogen

Analytical Reported **Concentration:** 208 ppm

Analytical Accuracy (U, k=2): Method:

+/-0.002 BAC(G/210L) NDIR [5.2 ppm]

Distributed by:

CMI Inc. 316 East Ninth St Owensboro, KY 4 Phone 866-835-0(www.alcoholtest.

Operator's Signature

*Traceable to: Certified Reference Material - 262.4 µmol/mol Ethanol in Nitrogen - Serial No. GN0015026 Lot No. 050319E11

Store in dry area, away from sources of heat, ignition and direct sunlight. Do not allow storage area to exceed 52 °C (125 °F).

Specialty Gas Lab Tech

03-11-2020

Issuance Date



The calibration results within this certificate were obtained using equipment and standards capable of producing analytical results traceable to NIST, and apply only to the items rine contained on this certificate. ILMO Products Company makes no warranty or representation as to the suitability of the use of any information provided for any particular purpose. The information use is at the sole discretion and risk of the user. Liability shall be limited to established replacement cost of this material or service.

OFFICER'S AFFIDAVIT AND NOTICE OF REVOCATION/DISQUALIFICATION

8/13

		C								Section
ARREST DATE	ARREST TIME	COUNTY #	CITY	'#	CIT	ATION #				
04/21/2020	0800	55	9	9	1	JONE				
ARREST LOCATION		CI	ΓY				COUNT	Y		
N		1	JONE				NO	\mathbf{NE}		
SUBJECT NAME			DATE OF			SEX	RACE	HEIGHT		/EIGHT
NORTH, GAYLON I	JEVERNE		08/	09/1	.961	M	W	51	0	165
ADDRESS		CITY			ST.	ATE		ZIP COI	DE	
3600 MARTIN LUT	'HER KING A	OKLAHOMA CIT	Y			ЭK		73:	133	
DRIVER LICENSE #		EXPIRATION DATE		STATE	CLASS		RESTR	ICTIONS	ENDO	RSEMENT
T083454871		06/30/20	019	OK	A				M	
VEHICLE MAKE	MODEL		TAG	#		ST	ATE	CMV	PLACA	\RD
N	N		N			C	ΟK	N	N	
On the above date, time, and location, the ab highway, street, turnpike, public parking lot, o influence of alcohol and/or other intoxicating s	r other public place, or upon a private ro substances as prohibited by law.									
(Describe driving behavior or circumstances):										

MOCK SUBJECT TEST

Describe person's condition (odor, actions, etc.):

ROUTINE MAINTENANCE. NEW GAS CANSITER INSTALL

THE PERSON WAS INFORMED OF THE IMPLIED CONSENT TEST REQUEST

BREATH-ALCOHOL ANAL	YSIS RECORD, REPORT OF THI	E PERSON NAMED	IN SECTION 1 ABOVE	and LOG OF TE	ST AND MAINTE	NANCE RECORD.			Section 2
EQUIPMENT				PECIALIST					
I-8000		80-		BOARD					
STANDARD			TARGET VALUE			ER			
Dry			0.080		16837				
MANUFACTURER			056200	80A2	$\frac{1}{04/20}$	/2020		TE CYLINDER EXPIRES $04/05/2022$	
OBSERVATION START	EST DATE	MOUTHPIECE				n/observation		.,,	
	04/21/2020	G							
RESULT TABLE				s observed p		-		THE STA	
Test	g/210L	Time	to th	e time of tes	t administra	ation.		0 ⁴ + * / * *	
			Ladminista	nod and D	maath Ala	ahal Test i		· · · · · · · · · · · · · · · · · · ·	
Air Blank	0.00	08:41	I administe					AAAA	
Subject Test	0.00	08:41	complianc	e with the	applicabl	e operating	g		7. E
Breath Volum			-	dure of the		· ·		2 ++ 1	(+**)后
Air Blank	0.00	08:42	proced	iure or me	Dualu U	I Tests.		P[+ 1/ + 1/ + 1/ + 1/ + 1/ + 1/ + 1/ + 1	**/9
Wait Air Blank	0.00	08:43 08:44						VOL	47/51
Subject Test	0.00	08:44							
Breath Volum		00.44	SIGNATURE OF (DPERATOR				1907	
Air Blank	0.00	08:44							
Cal Check	0.08	08:45	NAME				[Test Resu	lt
Air Blank	0.00	08:45	SMITH, JOS					BrAC g/21	
			BADGE #			1 (0) 7		DIAC g/21	
			00000004 AGENCY	±35	00000	10837		0.00	n
			BOARD OF	TESTS			L		
days or more, if you refuse (operating a CMV or CDL H as a result of this arrest by You may be eligible to pa revocation as a result of request IDAP within thirty	DAYS FROM DATE SERVED. Id or failed the state's test(s). In holder) OR a test result of .04 o filing a petition in the District C articipate in the Impaired Driv this arrest from appearing or y (30) days of the date listed	n addition, your con or more while opera court of the County in yer Accountability on your driving reco below. Mail the re	nmercial driving privile ting a commercial mol n which you were arre Program (IDAP) adm ord. Participation in equest for IDAP to th	ege in this state w for vehicle. You r ested. Your petiti inistered by the IDAP may reduct e Department of	ill be disqualifie may appeal any on must be filed Department c te the amount Public Safety	ed if you refused or v Departmental active d within thirty (30) of f Public Safety (E of fees you will b at P.O. Box 1141	failed the ion agains days of th DPS). Co e require 5, Oklaho	e state's test(s) st your driver license e date listed below. mpletion of IDAP may pr d to pay to the State. Yo oma City, OK 73136.	ou must
	E: I personally hand delive			ing the Notice	of Revocatio	on and the Noti	ce of De	epartment Impaired D	riving
	m to the person named in								
Date served: 04/	21/2020 Sign	nature of officer					E	Badge #:	
In accordance with Tit	21/2020 Sign le 12 O.S. Section 426, "I	l state under pe	nalty of perjury u	inder the laws	of Oklahom	a that the foreg	oing is t	rue and correct."	Section 4
Date: 04/21/20	020 Place (location w	hen signed):			Signature	of arresting offic	cer:		
NAME				AGENCY					
SMITH, JOS				BOARD	OF TE	STS			
BADGE 0000000435	AGENCY ADDRESS	NCY					PHONE 405	4252460	
OTHER WITNESSES	•								Section 5
NAME			TITLE	ADDRESS	····		PHO	NE	
1				1					
NAME			TITLE	ADDRESS			PHO	NE	
2									

OFFICER'S AFFIDAVIT AND NOTICE OF REVOCATION/DISQUALIFICATION

9/13

			-							Section
ARREST DATE	ARREST TIME		COUNTY #	CIT	Y#	C	TATION #			
04/21/2020	0800		55	9	9		NONE	;		
ARREST LOCATION			CIT	Y				COUNT	Y	
N			1	JONE				NO	NE	
SUBJECT NAME					FBIRTH		SEX	RACE	HEIGHT	WEIGHT
NORTH, GAYLON LEVE	RNE			08	/09/:	1961	M	W	510) 165
ADDRESS		CITY				S	TATE		ZIP COD	E
3600 MARTIN LUTHER	KING A	OKLAHC	MA CIT	Y			OK		731	133
DRIVER LICENSE #			PIRATION DATE		STATE	CLASS	3	RESTRI	CTIONS	ENDORSEMENT
T083454871		0	6/30/20	019	OK	A				M
VEHICLE MAKE	MODEL			TAC	G #		ST	ATE	CMV	PLACARD
N	N			N	1		0	ΟK	N	N

On the above date, time, and location, the above named person was arrested, and I had reasonable grounds to believe the person was driving or in actual physical control of a motor vehicle upon a public road, highway, street, turnpike, public parking lot, or other public place, or upon a private road, street, alley or lane which provides access to one or more single or multi-family dwelling within this state while under the influence of alcohol and/or other intoxicating substances as prohibited by law. (Describe driving behavior or circumstances):

MOCK SUBJECT TEST

Describe person's condition (odor, actions, etc.):

ROUTINE MAINTENANCE. NEW GAS CANSITER INSTALL

THE PERSON	WAS INFORMED (OF THE IMPLIED	CONSENT TEST	REQUEST

BREATH-ALCOHOL ANAL	YSIS RECORD, REPORT OF THE	PERSON NAMED	IN SECTION 1 ABOVE	E, and LOG OF TE	ST AND MAINTENANCE RECORD.	Section 2
EQUIPMENT		SERIAL	NUMBER	SPECIALIST		
I-8000		80-	-006163	BOARD	OF TESTS	
STANDARD			TARGET VALUE		PERMIT NUMBER	
Dry			0.080		16837	
MANUFACTURER			LOT #		DATE CYLINDER INSTALLED	DATE CYLINDER EXPIRES
ILMO			056200	080A2	04/20/2020	04/05/2022
	EST DATE	MOUTHPIECE	A 15-minut	e (or greater) deprivation/observation	
	04/21/2020	G			prior to and continuing	
RESULT TABLE						THE STAT
Test	g/210L	Time	to t	he time of te	st administration.	
Air Blank	0.00	08:41	I administe	ered said E	Breath-Alcohol Test in	
Subject Test	0.00	08:41	complianc	e with the	applicable operating	
Breath Volum			-		~~ ~ ~	
Air Blank	0.00	08:42	proce	dure of the	e Board of Tests.	H ++/
Wait		08:43				
Air Blank	0.00	08:44				
Subject Test	0.00	08:44	SIGNATURE OF	OPERATOR	······	1907
Breath Volum						
Air Blank	0.00	08:44	NAME			
Cal Check	0.08	08:45	SMITH, JO	SHIIA		Test Result
Air Blank	0.00	08:45	BADGE #		PERMIT #	BrAC g/210L
			0000000	435	0000016837	
			AGENCY			0.00
			BOARD OF	TESTS		
EFFECTIVE THIRTY (30)	DAYS FROM DATE SERVED.	Notice of Revoca	tion is given that you	ur driver's license	is revoked or denied for a period of o	ne hundred eighty (180) Section 3
days or more, if you refuse	ed or failed the state's test(s). Ir	addition, your cor	nmercial driving privil	ege in this state	will be disqualified if you refused or fai	led the state's test(s)
(operating a CMV or CDL	holder) OR a test result of .04 or	more while opera	ting a commercial mo	otor vehicle. You	may appeal any Departmental action	against your driver license
You may be eligible to n	articinate in the Impaired Driv	er Accountability	Program (IDAP) adu	ested. Your peti ministered by th	tion must be filed within thirty (30) day	s of the date listed below. b). Completion of IDAP may prevent a
revocation as a result of	this arrest from appearing on	your driving reco	ord. Participation in	IDAP may redu	ce the amount of fees you will be re	equired to pay to the State. You must
					of Public Safety at P.O. Box 11415, (
SERVING OF NOTIC	E: I personally hand deliv	ered a copy of	the above contair	ning the Notice	e of Revocation and the Notice	of Department Impaired Driving
Accountability Progra	m to the person named in	Section 1 abov	e.	J		
	21/2020 Sigr					Badge #:
In accordance with Til	le 12 O.S. Section 426, "I	state under pe	nalty of perjury	under the law	s of Oklahoma that the foregoir	g is true and correct." Section 4
Date: 04/21/2	020 Place (location wh	nen signed):			_ Signature of arresting officer	

NAME	AGENCY
SMITH, JOSHUA	BOARD OF TESTS
BADGE AGENCY ADDRESS	ZIP CODE PHONE
000000435 INTER-AGENCY	00000 4054252460

Certificate	of Analysis – Wet Bath	
12-10-19 late of Analysis	Labelled targenaiue (g/210L)	
19.050 st Number	A. 829	
JOSH SMITH	Contraction of the second seco	ABORATORIES, INC.

CERTIFICATE OF ANALYSIS

ALPHA INSTRUMENT 2 Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-004923 Certified Alcohol Reference Solution for Simulator 12/10/2019 14:52

Ţ	est	g/210L	Time	
A	ir Blank	0.000	 14:53	
С	al Check	0.019	14:53	
Ĥ	ir Blank	0,000	14:54	
C	al Check	0.020	14:55	
A	ir Blank	0.000	14:55	
C	al Check	0.019	14:56	
A	ir Blank	0.000	14:56	
C	al Check	0.020	14:57	
ß	ir Blank	0.000	14:58	
Ca	al Check	0.019	14:58	
R	ir Blank	0.000	14:59	
Ca	al Check	0.020	15:00	
81	r Blank	0.000	15:00	
Ca	I Check	0.019	15:01	
AI	r Blank	0.000	15:01	
Ca	I Check	0.020	15:02	
-AI	r Blank	0.000	15:02	
Ca	1 Check	0.020	15:03	
Ri	r Blank	0.000	15:04	
Ca	1 Check	0.020	15:04	
Ai	r Blank	0.000	15:05	
Ca	I Check Stats			
A	verage	0.0196		
S	td Deu	0.0005		
R	el Std Deu(%)	2.6345		

Random Samples of Lot Number Alcohol Reference Solution for Simulator were analyzed by chromatography on February 5, 2019, using a Perkin Elmer Gas gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.0236% (w/vol) ethyl alcohol. / The expiration date for this lot number is February 4, 2021 at 11:59 PM.

19050

of

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol 5:04 analysis instrument reading of 0.020 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability: Testing was conducted using Cerilliant Reference Standard lot number FN08031603 whose values are traceable to NIST. All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.

rator's Signature

Certificate of A	nalysis – Wet Bath	
12.10-19	,840	
	abellet tar Ralue (g/210L)	
19080	1040	1107
ot Number	Are recommended of the	Aboratories, Inc.
JOSH SMITH 34	590 NORTH 67th ST	REET . HARRISBURG, PA 17111- 4511 . TELEPHONE: 717-564-5470
OT Technician Name and Signature		BRT Form 6.0
In percentance with PPT 220 A	Maintenance, the above referenced sir	
	erence in maintenance of the Intoxilyze	
		The second se

RLPHA INSTRUMENT 2 Intoxilyzer - Alcohol Analyzer Model 8000 SN 80-004923

12/10/2019

15:16 Certified Alcohol Reference Solution for Simulator

Test	g/210L	Tine	
Air Blank	0.000	 15:17	
Cal Check	0.040	15:17	
Air Blank	0.000	15:18	
Cal Check	0.040	15:19	
Rir Blank	0,000	15:19	
Cal Check	0.040	15:20	
Air Blank	0.000	15:20	
Cal Check	0.040	15:21	
Air Blank	0.000	15:22	
Cal Check	0.040	15:22	
Air Blank	0.000	15:23	
Cal Check	0.040	15:24	
Air Blank	0.000	15:24	
Cal Check	0.040	15:25	
Air Blank	0.000	15:25	
Cal Check	0.039	15:26	
Air Blank	0.000	15:27	
Cal Check	0.039	15:27	
Air Blank	0,000	15:28	
Cal Check	0.039	15:29)
Air Blank	0.000	15:29	1
Cal Check Stats			1
Average	0.0397		1
Std Dev	0.0004		
Rel Std Deu(%)	1.2167		1

tor's Signatur

Random Samples of Lot Number 19080 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on March 5, 2019, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.0481% (w/vol) ethyl alcohol. The expiration date for this lot number is March 4, 2021 at 11:59 PM.

11/13

When used in a calibrated Simulator, operating at $34^{\circ}C$ +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.040 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN12181501 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.

Certificate of A	Analysis – Wet Bath		
12-10-19	Labelled target_value (g/210L)		
1939	100 . 100		
JOSH SMITH	SEO NORTH GTILL	ABORATORIES,	INC. • TELEPHONE: 717-564-5470
T Technician Name and Signature		BRT Form 6.0	
TTechnician Name and Signature accordance with BRT – 2.2.0,	Main ance, the above referenced sin	BRT Form 6.0 nulator solution is	

CERTIFICATE OF ANALYSIS

ALPHA INSTRUMENT 2 Intoxiluzer - Alcohol	Analuzer		
Model 8000		80-004923	÷
12/10/2019		15:36	

Certified Alcohol Reference Solution for Simulator

Test	9/210L	Time	
Air Blank	0.000	15:37	
Cal Check	0.100	15:37	
Air Blank	0.000	15:38	
Cal Check	0.101	15:39	
Air Blank	0.000	15:39	
Cal Check	0.100	15:40	
Air Blank	0.000	15:40	
Cal Check	0.100	15:41	
Air Blank	0.000	15:42	
Cal Check	0.100	15:42	
Rir Blank	0.000	15:43	
Cal Check	0.099	15:44	
Air Blank	0.000	15:44	
Cal Check	0.100	15:45	
AIC Blank	0.000	15:45	
Cal Check	0.099	15:46	
Air Blank	0.000	15:47	
Cal Check	0.100	15:47	
Air Blank	0.000	15:48	
Cal Check	0,100	15:49	
Air Blank	0.000	15:49	
Cal Check Stats			
Average	0.0999		
Std Deu	0.0005		
Rel Std Deu(%)	0.5682		

Random Samples of Lot Number 19301 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on November 1, 2019, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1210% (w/vol) ethyl alcohol. The expiration date for this lot number is October 30, 2021 at 11:59 PM.

When used in a calibrated Simulator, operating at $34^{\circ}C$ +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

s Signature

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN04271602 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.

12/13

Certificate	of Analysis – Wet Bath	
2-11-19	,200	
ate of Analysis	Labertod taren value (g/210L)	
19200	1200	- And -
ot Number		LABORATORIES, INC.
Jost SMITH	590 NORTH 67th	
OT Technician Name and Signat	ure	BRT Form 6.0
a pagagdapage with PPT 2	.2.0, Maintenance, the above referenced	simulator solution is

CERTIFICATE OF ANALYSIS

ALPHA INSTRUMENT 2			1
Intoxilyzer - Alcohol	Analyzer		1 -
Model 8000	SN	80-004923	C
12/11/2019		08:03	

Certified Alcohol Reference Solution for Simulator

Test c	g/210L	- Time
Air Blank	0.000	08:04
Cal Check	0.198	08:04
Air Blank	0.000	08:05
Cal Check	0.198	08:06
Air Blank	0.000	08;06
Cal Check	0.197	08:07
Air Blank	0,000	08:08
Cal Check	0.197	08:08
Air Blank	0.000	08:09
Cal Check	0.197	08:10
Air Blank	0.000	08:10
Cal Check	0.197	08:11
Air Blank	0.000	08:11
Cal Check	0.197	08:12
Air Blank	0.000	08:13
Cal Check	0.197	08:13
Air Blank	0.000	08:14
Cal Check	0.196	08:15
Air Blank	0.000	08:15
Cal Check	0.198	08:16
Air Blank	0.000	08:17
Cal Check Stats		
Average	0.1972	
Std Deu	0.0006	
Rel Std Deu(%)	0.3207	

Random Samples of Lot Number 19200 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on August 7, 2019, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.2418% (w/vol) ethyl alcohol. The expiration date for this lot number is August 6, 2021 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.200 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN08101505 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.